# limma

April 19, 2010

01.Introduction

Introduction to the LIMMA Package

## **Description**

LIMMA is a library for the analysis of gene expression microarray data, especially the use of linear models for analysing designed experiments and the assessment of differential expression. LIMMA provides the ability to analyse comparisons between many RNA targets simultaneously in arbitrary complicated designed experiments. Empirical Bayesian methods are used to provide stable results even when the number of arrays is small. The normalization and data analysis functions are for two-colour spotted microarrays. The linear model and differential expression functions apply to all microarrays including Affymetrix and other multi-array oligonucleotide experiments.

### **Details**

There are three types of documentation available:

- (1) The LIMMA User's Guide can be reached through the "User Guides and Package Vignettes" links at the top of the L
- (2) An overview of limma functions grouped by purpose is contained in the numbered chapters at the top of the LIMMA
- (3) The LIMMA contents page gives an alphabetical index of detailed help topics.

The function changeLog displays the record of changes to the package.

## Author(s)

Gordon Smyth

## References

Smyth, G. K., Yang, Y.-H., Speed, T. P. (2003). Statistical issues in microarray data analysis. *Methods in Molecular Biology* 224, 111-136.

Smyth, G. K. (2004). Linear models and empirical Bayes methods for assessing differential expression in microarray experiments. *Statistical Applications in Genetics and Molecular Biology*, Volume 3, Article 3. http://www.bepress.com/sagmb/vol3/iss1/art3

Smyth, G. K. (2005). Limma: linear models for microarray data. In: *Bioinformatics and Computational Biology Solutions using R and Bioconductor*. R. Gentleman, V. Carey, S. Dudoit, R. Irizarry, W. Huber (eds), Springer, New York, 2005.

2 02.Classes

02.Classes

Classes Defined by this Package

## **Description**

This package defines the following data classes.

- RGList A class used to store raw intensities as they are read in from an image analysis output file, usually by read.maimages.
- MAList Intensities converted to M-values and A-values, i.e., to with-spot and whole-spot contrasts on the log-scale. Usually created from an RGList using MA.RG or normalizeWithinArrays. Objects of this class contain one row for each spot. There may be more than one spot and therefore more than one row for each probe.
- EListRaw A class to store raw intensities for one-channel microarray data. May or may not be background corrected. Usually created by read.maimages.
- EList A class to store normalized log2 expression values for one-channel microarray data. Usually created by normalizeBetweenArrays.
- MArrayLM Store the result of fitting gene-wise linear models to the normalized intensities or log-ratios. Usually created by lmFit. Objects of this class normally contain only one row for each unique probe.
- TestResults Store the results of testing a set of contrasts equal to zero for each probe. Usually created by decideTests. Objects of this class normally contain one row for each unique probe.

All these data classes obey many analogies with matrices. In the case of RGList, MAList, EListRaw and EList, rows correspond to spots or probes and columns to arrays. In the case of MarrayLM, rows correspond to unique probes and the columns to parameters or contrasts. The functions summary, dim, length, ncol, nrow, dimnames, rownames, colnames have methods for these classes. Objects of any of these classes may be subsetted. Multiple data objects may be combined by rows (to add extra probes) or by columns (to add extra arrays).

Furthermore all of these classes may be coerced to actually be of class matrix using as .matrix, although this entails loss of information. Fitted model objects of class MArrayLM can be coerced to class data.frame using as.data.frame.

The first three classes belong to the virtual class LargeDataObject. A show method is defined for LargeDataOjects which uses the utility function printHead.

### Author(s)

Gordon Smyth

03.ReadingData 3

03.ReadingData

Reading Microarray Data from Files

### **Description**

This help page gives an overview of LIMMA functions used to read data from files.

## **Reading Target Information**

The function readTargets is designed to help with organizing information about which RNA sample is hybridized to each channel on each array and which files store information for each array.

### **Reading Intensity Data**

The first step in a microarray data analysis is to read into R the intensity data for each array provided by an image analysis program. This is done using the function read.maimages.

read.maimages optionally constructs quality weights for each spot using quality functions listed in QualityWeights.

If the data is two-color, then read.maimages produces an RGList object. If the data is one-color (single channel) then an EListRaw object is produced. In either case, read.maimages stores only the information required from each image analysis output file. read.maimages uses utility functions removeExt, read.imagene and read.columns. There are also a series of utility functions which read the header information from image output files including readGPRHeader, readImaGeneHeader and readGenericHeader.

read.ilmn reads probe or gene summary profile files from Illumina BeadChips, and produces an ElistRaw object.

The function as.MAList can be used to convert a marrayNorm object to an MAList object if the data was read and normalized using the marray and marrayNorm packages.

#### Reading the Gene List

Most image analysis software programs provide gene IDs as part of the intensity output files, for example GenePix, Imagene and the Stanford Microarray Database do this. In other cases the probe ID and annotation information may be in a separate file. The most common format for the probe annotation file is the GenePix Array List (GAL) file format. The function readGAL reads information from a GAL file and produces a data frame with standard column names.

The function <code>getLayout</code> extracts from the GAL-file data frame the print layout information for a spotted array. The functions <code>gridr</code>, <code>gridc</code>, <code>spotr</code> and <code>spotc</code> use the extracted layout to compute grid positions and spot positions within each grid for each spot. The function <code>printorder</code> calculates the printorder, plate number and plate row and column position for each spot given information about the printing process. The utility function <code>getSpacing</code> converts character strings specifying spacings of duplicate spots to numeric values.

The Australian Genome Research Facility in Australia often produces GAL files with composite probe IDs or names consisting of multiple strings separated by a delimiter. These can be separated into name and annotation information using strsplit2.

If each probe is printed more than once of the arrays in a regular pattern, then uniquegenelist will remove duplicate names from the gal-file or gene list.

4 04.Background

### **Identifying Control Spots**

The functions readSpotTypes and controlStatus assist with separating control spots from ordinary genes in the analysis and data exploration.

### **Manipulating Data Objects**

cbind, rbind, merge allow different RGList or MAList objects to be combined. cbind combines data from different arrays assuming the layout of the arrays to be the same. merge can combine data even when the order of the probes on the arrays has changed. merge uses utility function makeUnique.

### Author(s)

Gordon Smyth

04.Background

**Background Correction** 

## **Description**

This page deals with background correction methods for two-color microarray data.

Usually one doesn't need to explicitly ask for background correction of the intensities because this is done by default by normalizeWithinArrays, which subtracts the background from the foreground intensities before applying the normalization method. This default background correction method can be over-ridden by using backgroundCorrect which offers a number of alternative background correct methods to simple subtraction. The function backgroundCorrect is used to correct the RGList before applying normalizeWithinArrays.

The movingmin method of backgroundCorrect uses utility functions ma3x3.matrix and ma3x3.spottedarray.

The normexp method of backgroundCorrect uses utility functions normexp.fit and normexp.signal.

neqc performs normexp background correction and quantile normalization using control probes. This method uses utility functions normexp.fit.control to conduct the background correction.

kooperberg is a Bayesian background correction tool designed specifically for GenePix data. kooperberg is not currently used as the default method for GenePix data because it is computationally intensive. It requires several additional columns from the GenePix data files which can be read in using read.maimages and specifying the other.columns argument.

## Author(s)

Gordon Smyth

05.Normalization 5

05.Normalization Normalization of Microarray Data

### **Description**

This page gives an overview of the LIMMA functions available to normalize data from spotted two-colour microarrays. Smyth and Speed (2003) give an overview of the normalization techniques implemented in the functions.

Usually data from spotted microarrays will be normalized using normalizeWithinArrays. A minority of data will also be normalized using normalizeBetweenArrays if diagnostic plots suggest a difference in scale between the arrays.

In rare circumstances, data might be normalized using normalizeForPrintorder before using normalizeWithinArrays.

All the normalization routines take account of spot quality weights which might be set in the data objects. The weights can be temporarily modified using modifyWeights to, for example, remove ratio control spots from the normalization process.

If one is planning analysis of single-channel information from the microarrays rather than analysis of differential expression based on log-ratios, then the data should be normalized using a single channel-normalization technique. Single channel normalization uses further options of the normalizeBetweenArrays function. For more details see the LIMMA User's Guide which includes a section on single-channel normalization.

 ${\tt normalizeBetweenArrays} \ {\tt uses} \ {\tt utility} \ {\tt functions} \ {\tt normalizeMedianAbsValues} \ {\tt and} \ {\tt normalizeQuantile} \ {\tt none} \ {\tt of} \ {\tt which} \ {\tt need} \ {\tt to} \ {\tt be} \ {\tt called} \ {\tt directly} \ {\tt by} \ {\tt users}.$ 

normalizeWithinArrays uses utility functions MA.RG, loessFit and normalizeRobustSpline.

removeBatchEffect can be used to remove a batch effect, associated with hybridization time or some other technical variable, prior to unsupervised analysis.

## Author(s)

Gordon Smyth

### References

Smyth, G. K., and Speed, T. P. (2003). Normalization of cDNA microarray data. In: *METHODS: Selecting Candidate Genes from DNA Array Screens: Application to Neuroscience*, D. Carter (ed.). Methods Volume 31, Issue 4, December 2003, pages 265-273. http://www.statsci.org/smyth/pubs/normalize.pdf

06.LinearModels Linear Models for Microarrays

6 06.LinearModels

## Description

This page gives an overview of the LIMMA functions available to fit linear models and to interpret the results. This page covers models for two color arrays in terms of log-ratios or for single-channel arrays in terms of log-intensities. If you wish to fit models to the individual channel log-intensities from two colour arrays, see 07.SingleChannel.

The core of this package is the fitting of gene-wise linear models to microarray data. The basic idea is to estimate log-ratios between two or more target RNA samples simultaneously. See the LIMMA User's Guide for several case studies.

### **Fitting Models**

The main function for model fitting is lmFit. This is recommended interface for most users. lmFit produces a fitted model object of class MArrayLM containing coefficients, standard errors and residual standard errors for each gene. lmFit calls one of the following three functions to do the actual computations:

lm.series Straightforward least squares fitting of a linear model for each gene.

mrlm An alternative to lm.series using robust regression as implemented by the rlm function in the MASS package.

gls.series Generalized least squares taking into account correlations between duplicate spots (i.e., replicate spots on the same array) or related arrays. The function duplicateCorrelation is used to estimate the inter-duplicate or inter-block correlation before using gls.series.

All the functions which fit linear models use link{getEAW} to extract data from microarray data objects, and unwrapdups which provides an unified method for handling duplicate spots.

### **Forming the Design Matrix**

1mFit has two main arguments, the expression data and the design matrix. The design matrix is essentially an indicator matrix which specifies which target RNA samples were applied to each channel on each array. There is considerable freedom in choosing the design matrix - there is always more than one choice which is correct provided it is interpreted correctly.

Design matrices for Affymetrix or single-color arrays can be created using the function model.matrix which is part of the R base package. The function modelMatrix is provided to assist with creation of an appropriate design matrix for two-color microarray experiments. For direct two-color designs, without a common reference, the design matrix often needs to be created by hand.

### **Making Comparisons of Interest**

Once a linear model has been fit using an appropriate design matrix, the command makeContrasts may be used to form a contrast matrix to make comparisons of interest. The fit and the contrast matrix are used by contrasts.fit to compute fold changes and t-statistics for the contrasts of interest. This is a way to compute all possible pairwise comparisons between treatments for example in an experiment which compares many treatments to a common reference.

### **Assessing Differential Expression**

After fitting a linear model, the standard errors are moderated using a simple empirical Bayes model using eBayes or treat. ebayes is an older version of eBayes. A moderated t-statistic and a log-odds of differential expression is computed for each contrast for each gene. treat tests whether log-fold-changes are greater than a threshold rather than merely different to zero.

07.SingleChannel 7

eBayes and eBayes use internal functions squeezeVar, fitFDist, tmixture.matrix and tmixture.vector.

The function <code>zscoreT</code> is sometimes used for computing z-score equivalents for t-statistics so as to place t-statistics with different degrees of freedom on the same scale. <code>zscoreGamma</code> is used the same way with standard deviations instead of t-statistics. These functions are for research purposes rather than for routine use.

## **Summarizing Model Fits**

After the above steps the results may be displayed or further processed using:

toptable **or** topTable Presents a list of the genes most likely to be differentially expressed for a given contrast.

topTableF Presents a list of the genes most likely to be differentially expressed for a given set of contrasts.

volcanoplot Volcano plot of fold change versus the B-statistic for any fitted coefficient.

plotlines Plots fitted coefficients or log-intensity values for time-course data.

write.fit Writes an MarrayLM object to a file. Note that if fit is an MArrayLM object, either write.fit or write.table can be used to write the results to a delimited text file.

For multiple testing functions which operate on linear model fits, see 08. Tests.

#### **Model Selection**

selectModel provides a means to choose between alternative linear models using AIC or BIC information criteria.

### Author(s)

Gordon Smyth

#### References

Smyth, G. K. (2004). Linear models and empirical Bayes methods for assessing differential expression in microarray experiments. *Statistical Applications in Genetics and Molecular Biology*, **3**, No. 1, Article 3. http://www.bepress.com/sagmb/vol3/iss1/art3

Smyth, G. K., Michaud, J., and Scott, H. (2005). The use of within-array replicate spots for assessing differential expression in microarray experiments. Bioinformatics 21(9), 2067-2075.

07. SingleChannel Individual Channel Analysis of Two-Color Microarrays

### **Description**

This page gives an overview of the LIMMA functions fit linear models to two-color microarray data in terms of the log-intensities rather than log-ratios.

The function intraspotCorrelation estimates the intra-spot correlation between the two channels. The regression function lmscFit takes the correlation as an argument and fits linear models to the two-color data in terms of the individual log-intensities. The output of lmscFit is

8 08.Tests

an MArrayLM object just the same as from lmFit, so inference proceeds in the same way as for log-ratios once the linear model is fitted. See 06.LinearModels.

The function targetsA2C converts two-color format target data frames to single channel format, i.e, converts from array-per-line to channel-per-line, to facilitate the formulation of the design matrix.

### Author(s)

Gordon Smyth

08.Tests

Hypothesis Testing for Linear Models

### **Description**

LIMMA provides a number of functions for multiple testing across both contrasts and genes. The starting point is an MArrayLM object, called fit say, resulting from fitting a linear model and running eBayes and, optionally, contrasts.fit. See O6.LinearModels or O7.SingleChannel for details.

### Multiple testing across genes and contrasts

The key function is decideTests. This function writes an object of class TestResults, which is basically a matrix of -1, 0 or 1 elements, of the same dimension as fit\$coefficients, indicating whether each coefficient is significantly different from zero. A number of different multiple testing strategies are provided. The function calls other functions classifyTestsF, classifyTestsP and classifyTestsT which implement particular strategies. The function FStat provides an alternative interface to classifyTestsF to extract only the overall moderated F-statistic.

selectModel chooses between linear models for each probe using AIC or BIC criteria. This is an alternative to hypothesis testing and can choose between non-nested models.

A number of other functions are provided to display the results of decideTests. The functions heatDiagram (or the older version heatdiagram displays the results in a heat-map style display. This allows visual comparison of the results across many different conditions in the linear model.

The functions vennCounts and vennDiagram provide Venn diagram style summaries of the results.

Summary and show method exists for objects of class TestResults.

The results from decideTests can also be included when the results of a linear model fit are written to a file using write.fit.

## **Gene Set Tests**

Competitive gene set testing for an individual gene set is provided by wilcoxGST or geneSetTest, which permute genes. The gene set can be displayed using barcodeplot or barcodeplot2, the latter for directional gene sets with both up and down genes.

Self-contained gene set testing for an individual is provided by roast, which uses rotation technology, analogous to permuting arrays.

Gene set enrichment analysis for a battery of gene sets is provided by romer.

09.Diagnostics 9

The functions alias2Symbol and alias2SymbolTable are provided to help match gene sets with microarray probes by way of official gene symbols.

#### **Global Tests**

The function genas can test for associations between two contrasts in a linear model.

Given a set of p-values, the function convest can be used to estimate the proportion of true null hypotheses.

When evaluating test procedures with simulated or known results, the utility function auROC can be used to compute the area under the Receiver Operating Curve for the test results for a given probe.

### Author(s)

Gordon Smyth

09.Diagnostics

Diagnostics and Quality Assessment

## **Description**

This page gives an overview of the LIMMA functions available for microarray quality assessment and diagnostic plots.

This package provides an anova method which is designed for assessing the quality of an array series or of a normalization method. It is not designed to assess differential expression of individual genes. anova uses utility functions bwss and bwss.matrix.

The function arrayWeights estimates the empirical reliability of each array following a linear model fit.

Diagnostic plots can be produced by

- imageplot Produces a spatial picture of any spot-specific measure from an array image. If the log-ratios are plotted, then this produces an in-silico representation of the well known false-color TIFF image of an array. imageplot3by2 will write imageplots to files, six plots to a page.
- plotFB Plots foreground versus background log-intensies for a two-color array.
- plotMA MA-plots. One of the most useful plots of a two-color array. plotMA3by2 will write MA-plots to files, six plots to a page. mdplot can also be useful for comparing two one-channel microarrays.
- plotPrintTipLoess Produces a grid of MA-plots, one for each print-tip group on an array, together with the corresponding lowess curve. Intended to help visualize print-tip loess normalization.
- plotPrintorder For an array, produces a scatter plot of log-ratios or log-intensities by print order.
- plotDensities Individual channel densities for one or more arrays. An essential plot to accompany between array normalization, especially quantile normalization.
- plotMDS Multidimensional scaling plot for a set of arrays. Useful for visualizing the relationship between the set of samples.

10 alias2Symbol

plotSA Sigma vs A plot. After a linear model is fitted, this checks constancy of the variance with respect to intensity level.

plotPrintTipLoess uses utility functions gridr and gridc. plotDensities uses utility function RG.MA.

### Author(s)

Gordon Smyth

10.Other

Other Functions

## Description

This page describes some functions not covered in the previous numbered pages, so far only blockDiag and poolVar which are not used in the package yet but are part of the development of methods to handle technical and biological replicates.

### Author(s)

Gordon Smyth

alias2Symbol

Convert Gene Alias to Official Gene Symbols

## **Description**

Maps gene alias names to official gene symbols.

### Usage

```
alias2Symbol(alias, species = "Hs", expand.symbols = FALSE)
alias2SymbolTable(alias, species = "Hs")
```

## Arguments

alias character vector of gene aliases

species character string specifying the species. Possible values are "Dm", "Hs", "Mm" or "Rn".

expand.symbols

logical, should those elements of alias which are already official symbols be expanded if they are aliases for other symbols.

anova.MAList-method 11

#### **Details**

Aliases are mapped via NCBI Entrez Gene identity numbers using Bioconductor organism packages. Species are "Dm" for fly, "Hs" for human, "Mm" for mouse and "Rn" for rat. The user needs to have the appropriate Bioconductor organism package installed.

alias2Symbol maps a set of aliases to a set of symbols, without necessarily preserving order. The output vector may be longer or shorter than the original vector, because some aliases might not be found and some aliases may map to more than one symbol. alias2SymbolTable maps each alias to a gene symbol and returns a table with one row for each alias. If an alias maps to more than one symbol, then the first one found will be returned.

#### Value

Character vector of gene symbols.

alias2SymbolTable returns a vector of the same length and order as alias, including NA values where no gene symbol was found. alias2Symbol returns an unordered vector which may be longer or shorter than alias.

### Author(s)

Gordon Smyth and Yifang Hu

#### See Also

This function is often used to assist gene set testing, see 08. Tests.

## **Examples**

```
if(require("org.Hs.eg.db")) alias2Symbol(c("PUMA", "NOXA", "BIM"))
```

```
anova.MAList-method
```

ANOVA Table - method

### **Description**

Analysis of variance method for objects of class MAList. Produces an ANOVA table useful for quality assessment by decomposing between and within gene sums of squares for a series of replicate arrays. This method produces a single ANOVA Table rather than one for each gene and is not used to identify differentially expressed genes.

## Usage

```
anova(object, design=NULL, ndups=2,...)
```

12 arrayWeightsQuick

## **Arguments**

object object of class MAList. Missing values in the M-values are not allowed.

design numeric vector or single-column matrix containing the design matrix for linear model. The length of the vector or the number of rows of the matrix should agree with the number of columns of M.

ndups number of duplicate spots. Each gene is printed ndups times in adjacent spots on each array.

... other arguments are not used

### **Details**

This function aids in quality assessment of microarray data and in the comparison of normalization methodologies. It applies only to replicated two-color experiments in which all the arrays are hybridized with the same RNA targets, possibly with dye-swaps, so the design matrix should have only one column. The function has not been heavily used and is somewhat experimental.

#### Value

An object of class anova containing rows for between genes, between arrays, gene x array interaction, and between duplicate with array sums of squares. Variance components are estimated for each source of variation.

#### Note

This function does not give valid results in the presence of missing M-values.

## Author(s)

Gordon Smyth

### See Also

```
MAList-class, bwss.matrix, anova.
```

An overview of quality assessment and diagnostic functions in LIMMA is given by 09.Diagnostics.

arrayWeightsQuick Array Quality Weights

### **Description**

Estimates relative quality weights for each array in a multi-array experiment with replication.

## Usage

```
arrayWeightsQuick(y, fit)
```

#### **Arguments**

y the data object used to estimate fit. Can be of any class which can be coerced to matrix, including matrix, MAList, marrayNorm or ExpressionSet.

fit MArrayLM fitted model object

array Weights 13

#### **Details**

Estimates the relative reliability of each array by measuring how well the expression values for that array follow the linear model.

This is a quick and dirty version of arrayWeights.

### Value

Numeric vector of weights of length ncol (fit).

### Author(s)

Gordon Smyth

#### References

Ritchie, M. E., Diyagama, D., Neilson, van Laar, R., J., Dobrovic, A., Holloway, A., and Smyth, G. K. (2006). Empirical array quality weights for microarray data. BMC Bioinformatics. (Accepted 11 April 2006)

### See Also

See arrayWeights. An overview of LIMMA functions for reading data is given in 03.ReadingData.

### **Examples**

```
## Not run:
fit <- lmFit(y, design)
arrayWeightsQuick(y, fit)
## End(Not run)</pre>
```

arrayWeights

Array Quality Weights

### **Description**

Estimates relative quality weights for each array in a multi-array experiment.

## Usage

```
arrayWeights(object, design = NULL, weights = NULL, method = "genebygene", maxit arrayWeightsSimple(object, design = NULL, maxiter = 100, tol = 1e-6, maxratio =
```

## **Arguments**

object of class numeric, matrix, MAList, marrayNorm, ExpressionSet or PLMset containing log-ratios or log-values of expression for a series of mi-

croarrays.

design the design matrix of the microarray experiment, with rows corresponding to

arrays and columns to coefficients to be estimated. Defaults to the unit vector

meaning that the arrays are treated as replicates.

14 arrayWeights

weights	optional numeric matrix containing prior weights for each spot.
method	character string specifying the estimating algorithm to be used. Choices are "genebygene" and "reml".
maxiter	maximum number of iterations allowed.
tol	convergence tolerance.
maxratio	maximum ratio between largest and smallest weights before iteration stops
trace	logical variable. If true then output diagnostic information at each iteration of the '"reml" algorithm, or at every 1000th iteration of the '"genebygene" algorithm.

#### Details

The relative reliability of each array is estimated by measuring how well the expression values for that array follow the linear model.

The method is described in Ritchie et al (2006). A heteroscedastic model is fitted to the expression values for each gene by calling the function lm.wfit. The dispersion model is fitted to the squared residuals from the mean fit, and is set up to have array specific coefficients, which are updated in either full REML scoring iterations, or using an efficient gene-by-gene update algorithm. The final estimates of these array variances are converted to weights.

The data object object is interpreted as for lmFit. In particular, the arguments design and weights will be extracted from the data object if available and do not normally need to be set explicitly in the call; if any of these are set in the call then they will over-ride the slots or components in the data object.

arrayWeightsSimple is a fast version of arrayWeights with method="reml", no prior weights and no missing values.

#### Value

A vector of array weights.

#### Author(s)

Matthew Ritchie and Gordon Smyth

### References

Ritchie, M. E., Diyagama, D., Neilson, van Laar, R., J., Dobrovic, A., Holloway, A., and Smyth, G. K. (2006). Empirical array quality weights in the analysis of microarray data. BMC Bioinformatics 7, 261. http://www.biomedcentral.com/1471-2105/7/261/abstract

## See Also

An overview of linear model functions in limma is given by 06.LinearModels.

## **Examples**

```
if(require("sma")) {
# Subset of data from ApoAI case study in Limma User's Guide
# Can be loaded from sma package if it is installed
data(MouseArray)
# Avoid non-positive intensities
RG <- backgroundCorrect(mouse.data, method="half")</pre>
```

as.data.frame

```
MA <- normalizeWithinArrays(RG, mouse.setup)
MA <- normalizeBetweenArrays(MA, method="Aq")
targets <- data.frame(Cy3=I(rep("Pool",6)),Cy5=I(c("WT","WT","WT","KO","KO","KO")))
design <- modelMatrix(targets, ref="Pool")
arrayw <- arrayWeightsSimple(MA, design)
fit <- lmFit(MA, design, weights=arrayw)
fit2 <- contrasts.fit(fit, contrasts=c(-1,1))
fit2 <- eBayes(fit2)
# Use of array weights increases the significance of the top genes
topTable(fit2)
}
```

as.data.frame

Turn a Microarray Linear Model Object into a Dataframe

## **Description**

Turn a MArrayLM object into a data.frame.

## Usage

```
## S3 method for class 'MArrayLM':
as.data.frame(x, row.names = NULL, optional = FALSE, ...)
```

### **Arguments**

```
x an object of class MArrayLM

row.names NULL or a character vector giving the row names for the data frame. Missing values are not allowed.

optional logical. If TRUE, setting row names and converting column names (to syntactic names) is optional.

... additional arguments to be passed to or from methods.
```

## Details

This method combines all the components of x which have a row for each probe on the array into a data.frame.

### Value

A data.frame.

### Author(s)

Gordon Smyth

## See Also

```
as.data.frame in the base package.
```

02. Classes gives an overview of data classes used in LIMMA. 06. Linear Models gives an overview of linear model functions in LIMMA.

16 as.matrix

as.MAList

Convert marrayNorm Object to an MAList Object

## Description

Convert marrayNorm Object to an MAList Object

## Usage

```
as.MAList(object)
```

## **Arguments**

```
object an marrayNorm object
```

#### Value

Object of class MAList

### Author(s)

Gordon Smyth

#### See Also

02. Classes gives an overview of all the classes defined by this package.

as.matrix

Turn a Microarray Data Object into a Matrix

## **Description**

Turn a microarray data object into a numeric matrix by extracting the expression values.

## Usage

```
## S3 method for class 'MAList':
as.matrix(x,...)
```

## Arguments

```
    an object of class RGList, MAList, EList, MArrayLM, marrayNorm, PLMset, ExpressionSet, LumiBatch or vsn.
    additional arguments, not used for these methods.
```

#### **Details**

These methods extract the matrix of log-ratios, for MAList or marrayNorm objects, or the matrix of expression values for other expression objects such as EList or ExressionSet. For MArrayLM objects, the matrix of fitted coefficients is extracted.

These methods involve loss of information, so the original data object is not recoverable.

asMatrixWeights 17

#### Value

A numeric matrix.

### Author(s)

Gordon Smyth

#### See Also

```
as.matrix in the base package or exprs in the Biobase package.
```

02. Classes gives an overview of data classes used in LIMMA.

asMatrixWeights

asMatrixWeights

## **Description**

Convert probe-weights or array-weights to a matrix of weights.

## Usage

```
asMatrixWeights (weights, dim)
```

## **Arguments**

weights numeric matrix of weights, rows corresponding to probes and columns to arrays.

Or vector of probe weights. Or vector of array weights.

dim numeric dimension vector of length 2, i.e., the number of probes and the number

of arrays.

### **Details**

This function converts a vector or probe-weights or a vector of array-weights to a matrix of the correct size. Probe-weights are repeated across rows while array-weights are repeated down the columns. If weights has length equal to the number of probes, it is assumed to contain probe-weights. If it has length equal to the number of arrays, it is assumed to contain array-weights. If the number of probes is equal to the number of arrays, then weights is assumed to contain array-weights if it is a row-vector of the correct size, i.e., if it is a matrix with one row.

This function is used internally by the linear model fitting functions in limma.

### Value

Numeric matrix of dimension dim.

### Author(s)

Gordon Smyth

18 auROC

### See Also

```
modifyWeights.
```

An overview of functions in LIMMA used for fitting linear models is given in 06.LinearModels.

### **Examples**

```
asMatrixWeights(1:3,c(4,3))
asMatrixWeights(1:4,c(4,3))
```

auROC

Area Under Receiver Operating Curve

## Description

Compute exact area under the ROC for empirical data.

## Usage

```
auROC(truth, stat=NULL)
```

### **Arguments**

truth logical vector, or numeric vector of 0 and 1, indicating whether each case is a

true positive.

stat numeric vector containing test statistics used to rank cases, from largest to small-

est. If NULL, then truth is assumed to be already sorted in decreasing test

statistic order.

### **Details**

This function computes the exact area under an empirical ROC curve. Cases are ranked by stat from largest to smallest. The number of true and false discoveries are determined by how well the true states represented by truth match up with the observed statistics given by stat.

#### Value

Numeric vector giving area under the curve, 1 being perfect and 0 being the minimum, or  $\mathtt{NULL}$  if truth has zero length.

### Author(s)

Gordon Smyth

### See Also

See 08.Tests for other functions for testing and processing p-values.

See also AUC in the ROC package.

avedups 19

### **Examples**

```
auROC(c(1,1,0,0,0))
truth <- rbinom(30,size=1,prob=0.2)
stat <- rchisq(30,df=2)
auROC(truth,stat)</pre>
```

avedups

Average Over Duplicate Spots

## **Description**

Condense a microarray data object so that values for within-array replicate spots are replaced with their average.

## Usage

```
## Default S3 method:
avedups(x, ndups=2, spacing=1, weights=NULL)
## S3 method for class 'MAList':
avedups(x, ndups=x$printer$ndups, spacing=x$printer$spacing, weights=x$weights)
```

## **Arguments**

x a matrix-like object, usually a matrix or an MAList object.

ndups number of within-array replicates for each probe.

spacing number of spots to step from a probe to its duplicate.

weights numeric matrix of spot weights.

## Details

A new data object is computed in which each probe is represented by the (weighted) average of its duplicate spots. For a MAList object, the components  ${\tt M}$  and  ${\tt A}$  are both averaged in this way.

If x is of mode "character", then the duplicate values are assumed to be equal and the first is taken as the average.

#### Value

A data object of the same class as x with 1/ndups as many rows.

## Author(s)

Gordon Smyth

## See Also

dimnames in the base package.

02. Classes gives an overview of data classes used in LIMMA.

20 avereps

avereps

Average Over Irregular Replicate Spots

## **Description**

Condense a microarray data object so that values for within-array replicate spots are replaced with their average.

## Usage

```
## Default S3 method:
avereps(x, ID=rownames(x))
## S3 method for class 'MAList':
avereps(x, ID=NULL)
```

## **Arguments**

x a matrix-like object, usually a matrix or an MAList object.

ID probe identifier.

### **Details**

A new data object is computed in which each probe is represented by the average of its replicate spots. For an MAList object, the components M and A are both averaged in this way, as weights and any matrices found in objectSother.

For an MAList object, ID defaults to MA\$genes\$ID is that exists, otherwise to rownames (MA\$M).

If x is of mode "character", then the replicate values are assumed to be equal and the first is taken as the average.

### Value

A data object of the same class as x with a row for each unique value of ID.

## Author(s)

Gordon Smyth

## See Also

dimnames in the base package.

02. Classes gives an overview of data classes used in LIMMA.

backgroundCorrect 21

backgroundCorrect Correct Intensities for Background

#### **Description**

Background correct microarray expression intensities.

#### Usage

backgroundCorrect(RG, method=NULL, offset=0, printer=RG\$printer, normexp.method=

### **Arguments**

RG a numeric matrix, EListRaw or RGList object.

method character string specifying correction method. Possible values are "none",

"subtract", "half", "minimum", "movingmin", "edwards" or "normexp".

If RG is a matrix or an EListRaw object, possible values are restricted to "none" or "normexp". Defaults to "subtract" for two-color data ob-

jects or "normexp" for single-channel objects.

offset numeric value to add to intensities

printer a list containing printer layout information, see PrintLayout-class. Ig-

nored if RG is a matrix.

normexp.method

character string specifying parameter estimation strategy used by normexp, ignored for other methods. Possible values are "saddle", "mle", "rma" or

"rma75".

verbose logical. If TRUE, progress messages are sent to standard output

### **Details**

This function implements the background correction methods reviewed or developed in Ritchie et al (2007) and Silver at al (2009). Ritchie et al (2007) recommend method="normexp" whenever RG contains local background estimates. Silver et al (2009) shows that either normexp.method="mle" or normexp.method="saddle" are excellent options for normexp. If RG contains morphological background estimates instead (available from SPOT or GenePix image analysis software), then method="subtract" performs well.

If method="none" then no correction is done, i.e., the background intensities are treated as zero. If method="subtract" then the background intensities are subtracted from the foreground intensities. This is the traditional background correction method, but is not necessarily recommended. If method="movingmin" then the background estimates are replaced with the minimums of the backgrounds of the spot and its eight neighbors, i.e., the background is replaced by a moving minimum of 3x3 grids of spots.

The remaining methods are all designed to produce positive corrected intensities. If method="half" then any intensity which is less than 0.5 after background subtraction is reset to be equal to 0.5. If method="minimum" then any intensity which is zero or negative after background subtraction is set equal to half the minimum of the positive corrected intensities for that array. If method="edwards" a log-linear interpolation method is used to adjust lower intensities as in Edwards (2003). If method="normexp" a convolution of normal and exponential distributions is fitted to the foreground intensities using the background intensities as a covariate, and the expected signal given the

22 backgroundCorrect

observed foreground becomes the corrected intensity. This results in a smooth monotonic transformation of the background subtracted intensities such that all the corrected intensities are positive.

The normexp method is available in a number of variants depending on how the model parameters are estimated, and these are selected by normexp.method. Here "saddle" gives the saddle-point approximation to maximum likelihood from Ritchie et al (2007) and improved by Silver et al (2009), "mle" gives exact maximum likelihood from Silver at al (2009), "rma" gives the background correction algorithm from the RMA-algorithm for Affymetrix microarray data as implemented in the affy package, and "rma75" gives the RMA-75 method from McGee and Chen (2006). In practice "mle" performs well and is nearly as fast as "saddle", but "saddle" is the default for backward compatibility. See normexp.fit for more details.

The offset can be used to add a constant to the intensities before log-transforming, so that the log-ratios are shrunk towards zero at the lower intensities. This may eliminate or reverse the usual 'fanning' of log-ratios at low intensities associated with local background subtraction.

Background correction (background subtraction) is also performed by the normalizeWithinArrays method for RGList objects, so it is not necessary to call backgroundCorrect directly unless one wants to use a method other than simple subtraction. Calling backgroundCorrect before normalizeWithinArrays will over-ride the default background correction.

#### Value

An RGList object in which components R and G are background corrected and components Rb and Gb are removed.

### Author(s)

Gordon Smyth

### References

Edwards, D. E. (2003). Non-linear normalization and background correction in one-channel cDNA microarray studies *Bioinformatics* 19, 825-833.

McGee, M., and Chen, Z. (2006). Parameter estimation for the exponential-normal convolution model for background correction of Affymetrix GeneChip data. *Stat Appl Genet Mol Biol*, Volume 5, Article 24.

Ritchie, M. E., Silver, J., Oshlack, A., Silver, J., Holmes, M., Diyagama, D., Holloway, A., and Smyth, G. K. (2007). A comparison of background correction methods for two-colour microarrays. *Bioinformatics* 23, 2700-2707. http://bioinformatics.oxfordjournals.org/cgi/content/abstract/btm412

Silver, J., Ritchie, M. E., and Smyth, G. K. (2009). Microarray background correction: maximum likelihood estimation for the normal-exponential convolution model. *Biostatistics* 10, 352-363. http://biostatistics.oxfordjournals.org/cgi/content/abstract/kxn042

#### See Also

An overview of background correction functions is given in 04.Background.

#### **Examples**

```
RG <- new("RGList", list(R=c(1,2,3,4),G=c(1,2,3,4),Rb=c(2,2,2,2),Gb=c(2,2,2,2))) backgroundCorrect(RG) backgroundCorrect(RG, method="half") backgroundCorrect(RG, method="minimum") backgroundCorrect(RG, offset=5)
```

blockDiag 23

blockDiag

Block Diagonal Matrix

## **Description**

Form a block diagonal matrix from the given blocks.

## Usage

```
blockDiag(...)
```

## Arguments

... numeric matrices

## Value

A block diagonal matrix with dimensions equal to the sum of the input dimensions

### Author(s)

Gordon Smyth

### See Also

10.Other

## **Examples**

```
a <- matrix(1,3,2)
b <- matrix(2,2,2)
blockDiag(a,b)</pre>
```

bwss.matrix

Between and within sums of squares for matrix

## Description

Sums of squares between and within the columns of a matrix. Allows for missing values. This function is called by the anova method for MAList objects.

## Usage

```
bwss.matrix(x)
```

## **Arguments**

Х

a numeric matrix.

24 bwss

### **Details**

This is equivalent to a one-way analysis of variance where the columns of the matrix are the groups. If x is a matrix then bwss.matrix(x) is the same as bwss(x,col(x)) except for speed of execution.

## Value

## A list with components

bss sums of squares between the column means.

wss sums of squares within the column means.

bdf degrees of freedom corresponding to bss.

wdf degrees of freedom corresponding to wss.

## Author(s)

Gordon Smyth

## See Also

```
bwss, anova. MAList
```

bwss

Between and within sums of squares

### **Description**

Sums of squares between and within groups. Allows for missing values.

## Usage

```
bwss(x,group)
```

## Arguments

x a numeric vector giving the responses.

group a vector or factor giving the grouping variable.

## **Details**

This is equivalent to one-way analysis of variance.

## Value

## A list with components

bss	sums of squares between the group means.
WSS	sums of squares within the groups.
bdf	degrees of freedom corresponding to bss.
wdf	degrees of freedom corresponding to wss.

cbind 25

#### Author(s)

Gordon Smyth

#### See Also

```
bwss.matrix
```

cbind

Combine RGList, MAList, EList or EListRaw Objects

## **Description**

Combine a set of RGList, MAList, EList or EListRaw objects.

## Usage

```
## S3 method for class 'RGList':
cbind(..., deparse.level=1)
## S3 method for class 'RGList':
rbind(..., deparse.level=1)
```

## **Arguments**

```
... RGList, MAList, EList or EListRaw objects. deparse.level not currently used, see cbind in the base package
```

## **Details**

cbind combines data objects assuming the same probes in the same order but different arrays. rbind combines data objects assuming equivalent arrays, i.e., the same RNA targets, but different probes.

For cbind, the matrices of expression data from the individual objects are cbinded. The data frames of target information, if they exist, are rbinded. The combined data object will preserve any additional components or attributes found in the first object to be combined. For rbind, the matrices of expression data are rbinded while the target information, in any, is unchanged.

## Value

An RGList, MAList, EList or EListRaw object holding data from all the arrays and all genes from the individual objects.

### Author(s)

Gordon Smyth

### See Also

```
cbind in the base package.
```

03.ReadingData gives an overview of data input and manipulation functions in LIMMA.

26 changeLog

## **Examples**

```
M <- A <- matrix(11:14,4,2)
rownames(M) <- rownames(A) <- c("a","b","c","d")
colnames(M) <- colnames(A) <- c("A1","A2")
MA1 <- new("MAList",list(M=M,A=A))

M <- A <- matrix(21:24,4,2)
rownames(M) <- rownames(A) <- c("a","b","c","d")
colnames(M) <- colnames(A) <- c("B1","B2")
MA2 <- new("MAList",list(M=M,A=A))</pre>
cbind(MA1,MA2)
```

changeLog

Limma Change Log

## Description

Write as text the most recent changes from the limma package changelog.

## Usage

```
changeLog(n=20)
```

## **Arguments**

n

integer, number of lines to write of changelog.

## Value

No value is produced, but a number of lines of text are written to standard output.

## Author(s)

Gordon Smyth

#### See Also

01.Introduction

designI2M 27

designI2M

Convert Individual Channel Design Matrix to M-A Format

## **Description**

Convert a design matrix in terms of individual channels to ones in terms of M-values or A-values for two-color microarray data.

## Usage

```
designI2M(design)
designI2A(design)
```

### **Arguments**

design

numeric model matrix with one row for each channel observation, i.e., twice as many rows as arrays

## **Details**

If design is a model matrix suitable for modelling individual log-intensities for two color microarray data, then designI2M computes the corresponding model matrix for modelling M-values (log-ratios) and designI2A computes the model matrix for modelling A-values (average log-intensities).

Note that the matrices designI2M(design) or designI2A(design) may be singular if not all of the coefficients are estimable from the M or A-values. In that case there will be columns containing entirely zeros.

### Value

numeric model matrix with half as many rows as design

### Author(s)

Gordon Smyth

### See Also

```
model.matrix in the stats package.
```

An overview of individual channel linear model functions in limma is given by 07. Single Channel.

## **Examples**

```
X \leftarrow cbind(1,c(1,1,1,1,0,0,0),c(0,0,0,0,1,1,1,1))
designI2M(X)
designI2A(X)
```

28 classifyTests

classifyTests	Multiple Testing Genewise Across Contrasts	
---------------	--	--

### **Description**

For each gene, classify a series of related t-statistics as up, down or not significant.

## Usage

```
classifyTestsF(object, cor.matrix=NULL, df=Inf, p.value=0.01, fstat.only=FALSE)
classifyTestsT(object, t1=4, t2=3)
classifyTestsP(object, df=Inf, p.value=0.05, method="holm")
FStat(object, cor.matrix=NULL)
```

## **Arguments**

object	numeric matrix of t-statistics or an ${\tt MArrayLM}$ object from which the t-statistics may be extracted.
cor.matrix	covariance matrix of each row of t-statistics. Defaults to the identity matrix.
df	numeric vector giving the degrees of freedom for the t-statistics. May have length 1 or length equal to the number of rows of tstat.
p.value	numeric value between 0 and 1 giving the desired size of the test
fstat.only	logical, if TRUE then return the overall F-statistic as for FStat instead of classifying the test results
t1	first critical value for absolute t-statistics
t2	second critical value for absolute t-statistics
method	character string specifying p-value adjustment method. See p.adjust for possible values.

### **Details**

Note that these functions do not adjust for multiple testing across genes. The adjustment for multiple testing is across the contrasts rather than the more usual control across genes. The functions described here are called by decideTests. Most users should use decideTests rather than using these functions directly.

These functions implement multiple testing procedures for determining whether each statistic in a matrix of t-statistics should be considered significantly different from zero. Rows of tstat correspond to genes and columns to coefficients or contrasts.

FStat computes the gene-wise F-statistics for testing all the contrasts equal to zero. It is equivalent to classifyTestsF with fstat.only=TRUE.

classifyTestsF uses a nested F-test approach giving particular attention to correctly classifying genes which have two or more significant t-statistics, i.e., are differential expressed under two or more conditions. For each row of tstat, the overall F-statistic is constructed from the t-statistics as for FStat. At least one constrast will be classified as significant if and only if the overall F-statistic is significant. If the overall F-statistic is significant, then the function makes a best choice as to which t-statistics contributed to this result. The methodology is based on the principle that any t-statistic should be called significant if the F-test is still significant for that row when all the larger t-statistics are set to the same absolute size as the t-statistic in question.

contrasts.fit 29

classifyTestsT and classifyTestsP implement simpler classification schemes based on threshold or critical values for the individual t-statistics in the case of classifyTestsT or p-values obtained from the t-statistics in the case of classifyTestsP. For classifyTestsT, classifies any t-statistic with absolute greater than t2 as significant provided that at least one t-statistic for that gene is at least t1 in absolute value. classifyTestsP applied p-value adjustment from p.adjust to the p-values for each gene.

If tstat is an MArrayLM object, then all arguments except for p.value are extracted from it.

cor.matrix is the same as the correlation matrix of the coefficients from which the t-statistics are calculated. If cor.matrix is not specified, then it is calculated from design and contrasts if at least design is specified or else defaults to the identity matrix. In terms of design and contrasts, cor.matrix is obtained by standardizing the matrix t (contrasts) %\*% solve(t(design) %\*% design) %\*% contrasts to a correlation matrix.

#### Value

An object of class TestResults. This is essentially a numeric matrix with elements -1, 0 or 1 depending on whether each t-statistic is classified as significantly negative, not significant or significantly positive respectively.

FStat produces a numeric vector of F-statistics with attributes df1 and df2 giving the corresponding degrees of freedom.

#### Author(s)

Gordon Smyth

### See Also

An overview of multiple testing functions is given in 08. Tests.

### **Examples**

```
tstat <- matrix(c(0,5,0,0,2.5,0,-2,-2,2,1,1,1), 4, 3, byrow=TRUE) classifyTestsF(tstat)

# See also the examples for contrasts.fit and vennDiagram
```

contrasts.fit

Compute Contrasts from Linear Model Fit

### **Description**

Given a linear model fit to microarray data, compute estimated coefficients and standard errors for a given set of contrasts.

## Usage

```
contrasts.fit(fit, contrasts=NULL, coefficients=NULL)
```

30 contrasts.fit

#### **Arguments**

an MArrayLM object or a list object produced by the function lm.series or equivalent. Must contain components coefficients and stdev.unscaled.

contrasts

numeric matrix with row corresponding to coefficients in fit and columns containing contrasts. May be a vector if there is only one contrast.

coefficients vector indicating which coefficients are to be kept in the revised fit object. An alternative way to specify the contrasts.

### **Details**

This function accepts input from any of the functions lmFit, lm.series, mrlm, gls.series or lmscFit. The function re-orientates the fitted model object from the coefficients of the original design matrix to any set of contrasts of the original coefficients. The coefficients, unscaled standard deviations and correlation matrix are re-calculated in terms of the contrasts.

The idea of this function is to fit a full-rank model using lmFit or equivalent, then use contrasts.fit to obtain coefficients and standard errors for any number of contrasts of the coefficients of the original model. Unlike the design matrix input to lmFit, which normally has one column for each treatment in the experiment, the matrix contrasts may have any number of columns and these are not required to be linearly independent. Methods of assessing differential expression, such as eBayes or classifyTestsF, can then be applied to fitted model object.

The coefficients argument provides a simpler way to specify the contrasts matrix when the desired contrasts are just a subset of the original coefficients.

Warning. For efficiency reasons, this function does not re-factorize the design matrix for each probe. A consequence is that, if the design matrix is non-orthogonal and the original fit included quality weights or missing values, then the unscaled standard deviations produced by this function are approximate rather than exact. The approximation is usually acceptable. The results are always exact if the original fit was a oneway model.

### Value

An list object of the same class as fit, usually MArrayLM. This is a list with components

coefficients numeric matrix containing the estimated coefficients for each contrast for each probe.

stdev.unscaled

numeric matrix conformal with coef containing the unscaled standard deviations for the coefficient estimators.

cov.coefficients:

numeric matrix giving the unscaled covariance matrix of the estimable coefficients

... any other components input in fit

#### Author(s)

Gordon Smyth

### See Also

An overview of linear model functions in limma is given by 06.LinearModels.

controlStatus 31

#### **Examples**

```
# Simulate gene expression data: 6 microarrays and 100 genes
# with one gene differentially expressed in first 3 arrays
M <- matrix(rnorm(100*6,sd=0.3),100,6)
M[1,1:3] <- M[1,1:3] + 2
# Design matrix corresponds to oneway layout, columns are orthogonal
design <- cbind(First3Arrays=c(1,1,1,0,0,0),Last3Arrays=c(0,0,0,1,1,1))
fit <- lmFit(M,design=design)
# Would like to consider original two estimates plus difference between first 3 and last
contrast.matrix <- cbind(First3=c(1,0),Last3=c(0,1),"Last3-First3"=c(-1,1))
fit2 <- contrasts.fit(fit,contrast.matrix)
fit2 <- eBayes(fit2)
# Large values of eb$t indicate differential expression
results <- classifyTestsF(fit2)
vennCounts(results)</pre>
```

controlStatus

Set Status of each Spot from List of Spot Types

### **Description**

Determine the type (or status) of each spot in the gene list.

### Usage

```
controlStatus(types, genes, spottypecol="SpotType", regexpcol, verbose=TRUE)
```

#### **Arguments**

types	dataframe containing spot type specifiers, usually input using readSpotTypes
genes	dataframe containing the microarray gene list, or an ${\tt RGList}, {\tt MAList}$ or ${\tt MArrayList}$ containing genes
spottypecol	integer or name specifying column of types containing spot type names
regexpcol	vector of integers or column names specifying columns of types containing regular expressions. Defaults to any column names in common between types and genes.
verbose	logical, if TRUE then progess on pattern matching is reported to the standard output channel

### **Details**

This function constructs a vector of status codes by searching for patterns in the gene list. The data frame genes contains gene IDs and should have as many rows as there are spots on the microarrays. Such a data frame is often read using readGAL. The data frame types has as many rows as you want to distinguish types of spots in the gene list. This data frame should contain a column or columns, the regexpcol columns, which have the same names as columns in genes and which contain patterns to match in the gene list. Another column, the spottypecol, contains the names of the spot types. Any other columns are assumed to contain plotting parameters, such as colors or symbols, to be associated with the spot types.

The patterns in the regexpcol columns are simplified regular expressions. For example, AA\* means any string starting with AA, \*AA means any code ending with AA, AA means exactly these

32 convest

two letters, \*AA\* means any string containing AA, AA. means AA followed by exactly one other character and AA\. means exactly AA followed by a period and no other characters. Any other regular expressions are allowed but the codes  $^$  for beginning of string and \$ for end of string should not be included.

Note that the patterns are matched sequentially from first to last, so more general patterns should be included first. For example, it is often a good idea to include a default spot-type as the first line in types with pattern \* for all regexpcol columns and default plotting parameters.

### Value

Character vector specifying the type (or status) of each spot on the array. Attributes contain plotting parameters associated with each spot type.

### Author(s)

Gordon Smyth

#### See Also

An overview of LIMMA functions for reading data is given in 03.ReadingData.

## **Examples**

```
genes <- data.frame(ID=c("Control", "Control", "Control", "Control", "AA1", "AA2", "AA3", "AA4")
Name=c("Ratio 1", "Ratio 2", "House keeping 1", "House keeping 2", "Gene 1", "Gene 2", "Gene 3'
types <- data.frame(SpotType=c("Gene", "Ratio", "Housekeeping"), ID=c("*", "Control", "Control
status <- controlStatus(types, genes)
```

convest

Estimate Proportion of True Null Hypotheses

## Description

Returns an estimate of the proportion of true null hypotheses using a convex decreasing density estimate on a vector of p-values.

## Usage

```
convest(p, niter = 100, doplot = FALSE, doreport = FALSE)
```

## Arguments

р	numeric vector of p-values, calculated using any test of your choice. Missing values are not allowed
niter	number of iterations to be used in fitting the convex, decreasing density for the p-values. Default is 100.
doplot	logical, should updated plots of fitted convex decreasing p-value density be produced at each iteration? Default is FALSE.
doreport	logical, should the estimated proportion be printed at each iteration? Default is FALSE.

decideTests 33

#### **Details**

The proportion of true null hypotheses is often denoted pi0.

#### Value

Numeric value in the interval [0,1] representing he estimated proportion of true null hypotheses.

### Author(s)

Egil Ferkingstad and Mette Langaas

#### References

Langaas, M., Ferkingstad, E., and Lindqvist, B. (2005). Estimating the proportion of true null hypotheses, with application to DNA microarray data. *Journal of the Royal Statistical Society Series* B, 67, 555-572. Preprint at http://www.math.ntnu.no/~mettela/

### See Also

See 08.Tests for other functions for producing or interpretting p-values.

### **Examples**

```
# First simulate data, use no.genes genes and no.ind individuals,
# with given value of pi0. Draw from normal distribution with mean=0
# (true null) and mean=mean.diff (false null).

no.genes <- 1000
no.ind <- 20
pi0 <- 0.9
mean.diff <- 1
n1 <- round(pi0*no.ind*no.genes)
n2 <- round((1-pi0)*no.ind*no.genes)
x <- matrix(c(rnorm(n1,mean=0),rnorm(n2,mean=mean.diff)),ncol=no.ind,byrow=TRUE)
# calculate p-values using your favorite method, e.g.
pvals <- ebayes(lm.series(x))$p.value
# run the convex decreasing density estimator to estimate pi0
convest(pvals,niter=100,doplot=interactive())</pre>
```

decideTests

Multiple Testing Across Genes and Contrasts

## Description

Classify a series of related t-statistics as up, down or not significant. A number of different multiple testing schemes are offered which adjust for multiple testing down the genes as well as across contrasts for each gene.

## Usage

```
decideTests(object,method="separate",adjust.method="BH",p.value=0.05,lfc=0)
```

34 decideTests

## **Arguments**

object MArrayLM object output from eBayes from which the t-statistics may be ex-

tracted.

method character string specify how probes and contrasts are to be combined in the mul-

tiple testing strategy. Choices are "separate", "global", "hierarchical",

"nestedF" or any partial string.

adjust.method

character string specifying p-value adjustment method. Possible values are "none",

"BH", "fdr" (equivalent to "BH"), "BY" and "holm". See p.adjust for

details.

p.value numeric value between 0 and 1 giving the desired size of the test

lfc minimum log2-fold-change required

#### **Details**

These functions implement multiple testing procedures for determining whether each statistic in a matrix of t-statistics should be considered significantly different from zero. Rows of tstat correspond to genes and columns to coefficients or contrasts.

The setting method="separate" is equivalent to using topTable separately for each coefficient in the linear model fit, and will give the same lists of probes if adjust.method is the same. method="global" will treat the entire matrix of t-statistics as a single vector of unrelated tests. method="hierarchical" adjusts down genes and then across contrasts. method="nestedF" adjusts down genes and then uses classifyTestsF to classify contrasts as significant or not for the selected genes. Please see the limma User's Guide for a discussion of the statistical properties of these methods.

#### Value

An object of class TestResults. This is essentially a numeric matrix with elements -1, 0 or 1 depending on whether each t-statistic is classified as significantly negative, not significant or significantly positive respectively.

If lfc>0 then contrasts are judged significant only when the log2-fold change is at least this large in absolute value. For example, one might choose lfc=log2 (1.5) to restrict to 50% changes or lfc=1 for 2-fold changes. In this case, contrasts must satisfy both the p-value and the fold-change cutoff to be judged significant.

### Author(s)

Gordon Smyth

### See Also

An overview of multiple testing functions is given in 08.Tests.

dimnames 35

dimnames	Retrieve the Dimension Names of an RGList, MAList, EList, EListRaw or MArrayLM Object

## Description

Retrieve the dimension names of a microarray data object.

## Usage

```
## S3 method for class 'RGList':
dimnames(x)
## S3 replacement method for class 'RGList':
dimnames(x) <- value</pre>
```

## Arguments

X	an object of class RGList, MAList, EList, EListRaw or (not for assign-
	ment) MArrayLM
value	a possible value for dimnames (x): see dimnames

### **Details**

The dimension names of a microarray object are the same as those of the most important matrix component of that object.

A consequence is that rownames and colnames will work as expected.

## Value

Either NULL or a list of length 2. If a list, its components are either 'NULL' or a character vector the length of the appropriate dimension of  $\mathbf{x}$ .

## Author(s)

Gordon Smyth

## See Also

dimnames in the base package.

02. Classes gives an overview of data classes used in LIMMA.

36 dim

dim

Retrieve the Dimensions of an RGList, MAList or MArrayLM Object

## **Description**

Retrieve the number of rows (genes) and columns (arrays) for an RGList, MAList or MArrayLM object.

### Usage

```
## S3 method for class 'RGList':
dim(x)
## S3 method for class 'RGList':
length(x)
```

## **Arguments**

Х

an object of class RGList, MAList or MArrayLM

### **Details**

Microarray data objects share many analogies with ordinary matrices in which the rows correspond to spots or genes and the columns to arrays. These methods allow one to extract the size of microarray data objects in the same way that one would do for ordinary matrices.

A consequence is that row and column commands nrow(x), ncol(x) and so on also work.

## Value

Numeric vector of length 2. The first element is the number of rows (genes) and the second is the number of columns (arrays).

## Author(s)

Gordon Smyth

## See Also

dim in the base package.

02. Classes gives an overview of data classes used in LIMMA.

## **Examples**

```
M <- A <- matrix(11:14,4,2)
rownames(M) <- rownames(A) <- c("a","b","c","d")
colnames(M) <- colnames(A) <- c("A1","A2")
MA <- new("MAList",list(M=M,A=A))
dim(M)
ncol(M)
nrow(M)
length(M)</pre>
```

dupcor 37

dupcor	Correlation Between Duplicates	

## **Description**

Estimate the correlation between duplicate spots (regularly spaced replicate spots on the same array) or between technical replicates from a series of arrays.

## Usage

#### Arguments

object	a numeric matrix of expression values, or any data object from which as .matrix will extract a suitable matrix such as an MAList, marrayNorm or ExpressionSet object. If object is an MAList object then the arguments design, ndups, spacing and weights will be extracted from it if available and do not have to be specified as arguments. Specifying these arguments explicitly will over-rule any components found in the data object.
design	the design matrix of the microarray experiment, with rows corresponding to arrays and columns to comparisons to be estimated. The number of rows must match the number of columns of object. Defaults to the unit vector meaning that the arrays are treated as replicates.
ndups	a positive integer giving the number of times each gene is printed on an array. nrow (object) must be divisible by ndups. Will be ignored if block is specified.
spacing	the spacing between the rows of object corresponding to duplicate spots, spacing=1 for consecutive spots
block	vector or factor specifying a blocking variable
trim	the fraction of observations to be trimmed from each end of tanh (all.correlations) when computing the trimmed mean.
weights	an optional numeric matrix of the same dimension as object containing weights for each spot. If smaller than object then it will be filled out the same size.

### **Details**

When block=NULL, this function estimates the correlation between duplicate spots (regularly spaced within-array replicate spots). If block is not null, this function estimates the correlation between repeated observations on the blocking variable. Typically the blocks are biological replicates and the repeated observations are technical replicates. In either case, the correlation is estimated by fitting a mixed linear model by REML individually for each gene. The function also returns a consensus correlation, which is a robust average of the individual correlations, which can be used as input for functions lmFit or gls.series.

At this time it is not possible to estimate correlations between duplicate spots and between technical replicates simultaneously. If block is not null, then the function will set ndups=1.

For this function to return statistically useful results, there must be at least two more arrays than the number of coefficients to be estimated, i.e., two more than the column rank of design.

38 dupcor

The function may take long time to execute as it fits a mixed linear model for each gene for an iterative algorithm. It is not uncommon for the function to return a small number of warning messages that correlation estimates cannot be computed for some individual genes. This is not a serious concern providing that there are only a few such warnings and the total number of genes is large. The consensus estimator computed by this function will not be materially affected by a small number of genes.

#### Value

### A list with components

```
consensus.correlation
```

the average estimated inter-duplicate correlation. The average is the trimmed mean of the individual correlations on the atanh-transformed scale.

cor

same as consensus.correlation, for compatibility with earlier versions of the software

atanh.correlations

numeric vector of length nrow (object) /ndups giving the individual genewise atanh-transformed correlations.

#### Author(s)

Gordon Smyth

# References

Smyth, G. K., Michaud, J., and Scott, H. (2005). The use of within-array replicate spots for assessing differential expression in microarray experiments. *Bioinformatics* 21(9), 2067-2075. http://www.statsci.org/smyth/pubs/dupcor.pdf

# See Also

These functions use mixedModel2Fit from the statmod package.

An overview of linear model functions in limma is given by 06.LinearModels.

# **Examples**

```
## Also see lmFit examples

## Not run:
corfit <- duplicateCorrelation(MA, ndups=2, design)
all.correlations <- tanh(corfit$atanh.correlations)
boxplot(all.correlations)
fit <- lmFit(MA, design, ndups=2, correlation=corfit$consensus)

## End(Not run)</pre>
```

ebayes 39

ebayes

Empirical Bayes Statistics for Differential Expression

## **Description**

Given a series of related parameter estimates and standard errors, compute moderated t-statistics, moderated F-statistic, and log-odds of differential expression by empirical Bayes shrinkage of the standard errors towards a common value.

### Usage

```
ebayes(fit, proportion=0.01, stdev.coef.lim=c(0.1,4))
eBayes(fit, proportion=0.01, stdev.coef.lim=c(0.1,4))
treat(fit, lfc=0)
```

## Arguments

fit	an MArrayLM fitted model object produced by lmFit or contrasts.fit, or an unclassed list produced by lm.series, gls.series or mrlm containing components coefficients, stdev.unscaled, sigma and df.residual
proportion	numeric value between 0 and 1, assumed proportion of genes which are differentially expressed
stdev.coef.l	im
	numeric vector of length 2, assumed lower and upper limits for the standard deviation of log2-fold-changes for differentially expressed genes
lfc	the minimum log2-fold-change which is considered material

### **Details**

These functions is used to rank genes in order of evidence for differential expression. They use an empirical Bayes method to shrink the probe-wise sample variances towards a common value and to augmenting the degrees of freedom for the individual variances (Smyth, 2004). The functions accept as input argument fit a fitted model object from the functions lmFit, lm.series, mrlm or gls.series. The fitted model object may have been processed by contrasts.fit before being passed to eBayes to convert the coefficients of the design matrix into an arbitrary number of contrasts which are to be tested equal to zero. The columns of fit define a set of contrasts which are to be tested equal to zero.

The empirical Bayes moderated t-statistics test each individual contrast equal to zero. For each probe (row), the moderated F-statistic tests whether all the contrasts are zero. The F-statistic is an overall test computed from the set of t-statistics for that probe. This is exactly analogous the relationship between t-tests and F-statistics in conventional anova, except that the residual mean squares and residual degrees of freedom have been moderated between probes.

The estimates s2.prior and df.prior are computed by fitFDist. s2.post is the weighted average of s2.prior and sigma^2 with weights proportional to df.prior and df.residual respectively. The lods is sometimes known as the B-statistic. The F-statistics F are computed by classifyTestsF with fstat.only=TRUE.

eBayes doesn't compute ordinary (unmoderated) t-statistics by default, but these can be easily extracted from the linear model output, see the example below.

40 ebayes

ebayes is the earlier and leaner function. eBayes is intended to have a more object-orientated flavor as it produces objects containing all the necessary components for downstream analysis.

treat computes empirical Bayes moderated-t p-values relative to a minimum required fold-change threshold. Use topTreat to summarize output from treat. Instead of testing for genes which have log-fold-changes different from zero, it tests whether the log2-fold-change is greater than lfc in absolute value (McCarthy and Smyth, 2009). treat is concerned with p-values rather than posterior odds, so it does not compute the B-statistic lods. The idea of thresholding doesn't apply to F-statistics in a straightforward way, so moderated F-statistics are also not computed.

#### Value

eBayes produces an object of class MArrayLM with the following components, see MArrayLM-class. ebayes produces an ordinary list without F or F.p.value. treat produces an MArrayLM object, but without lods, var.prior, F or F.p.value.

t	numeric vector or matrix of moderated t-statistics
p.value	numeric vector of p-values corresponding to the t-statistics
s2.prior	estimated prior value for sigma^2
df.prior	degrees of freedom associated with s2.prior
s2.post	vector giving the posterior values for sigma^2
lods	numeric vector or matrix giving the log-odds of differential expression
var.prior	estimated prior value for the variance of the log2-fold-change for differentially expressed gene
F	numeric vector of moderated F-statistics for testing all contrasts defined by the columns of $\texttt{fit}$ simultaneously equal to zero
F.p.value	numeric vector giving p-values corresponding to F

## Author(s)

Gordon Smyth

# References

McCarthy, D. J., and Smyth, G. K. (2009). Testing significance relative to a fold-change threshold is a TREAT. *Bioinformatics*. http://bioinformatics.oxfordjournals.org/cgi/content/abstract/btp053

Loennstedt, I., and Speed, T. P. (2002). Replicated microarray data. Statistica Sinica 12, 31-46.

Smyth, G. K. (2004). Linear models and empirical Bayes methods for assessing differential expression in microarray experiments. *Statistical Applications in Genetics and Molecular Biology*, Volume 3, Article 3. http://www.bepress.com/sagmb/vol3/iss1/art3

#### See Also

```
squeezeVar, fitFDist, tmixture.matrix.
```

An overview of linear model functions in limma is given by 06.LinearModels.

EList-class 41

#### **Examples**

```
See also lmFit examples
  Simulate gene expression data,
  6 microarrays and 100 genes with one gene differentially expressed
set.seed(2004); invisible(runif(100))
M \leftarrow matrix(rnorm(100*6, sd=0.3), 100, 6)
M[1,] \leftarrow M[1,] + 1
fit <- lmFit(M)
# Ordinary t-statistic
par(mfrow=c(1,2))
ordinary.t <- fit$coef / fit$stdev.unscaled / fit$sigma</pre>
qqt(ordinary.t,df=fit$df.residual,main="Ordinary t")
abline(0,1)
# Moderated t-statistic
eb <- eBayes(fit)
qqt(eb$t,df=eb$df.prior+eb$df.residual,main="Moderated t")
abline(0,1)
# Points off the line may be differentially expressed
par(mfrow=c(1,1))
```

EList-class

Expression List - class

### **Description**

Simple list-based classes for storing expression values (E-values) for a set of one-channel microarrays. EListRaw holds expression values on the raw scale. EList holds expression values on the log scale, usually after background correction and normalization. EListRaw objects are normally created by read.maimages. In the future, EList objects are likely to be created normalizeBetweenArrays.

# **Slots/List Components**

EList objects can be created by new("EList", E) where E is a list. These classes contains no slots (other than .Data), but objects should contain a list component E as follows:

E: numeric matrix containing the E-values (raw or log-2 expression ratios). Rows correspond to spots and columns to an

Optional components include:

weights: numeric matrix of same dimensions as E containing relative spot quality weights. Elements should be non-nother: list containing other matrices, all of the same dimensions as E.

data.frame containing probe information. Should have one row for each probe. May have any number of collargets: data.frame containing information on the target RNA samples. Rows correspond to arrays. May have any number of collargets:

Valid EList or EListRaw objects may contain other optional components, but all probe or array information should be contained in the above components.

42 exprs.MA

#### Methods

These classes inherit directly from class list so any operation appropriate for lists will work on objects of this class. In addition, EList objects can be subsetted and combined. EList objects will return dimensions and hence functions such as dim, nrow and ncol are defined. ELists also inherit a show method from the virtual class LargeDataObject, which means that ELists will print in a compact way.

### Author(s)

Gordon Smyth

### See Also

02. Classes gives an overview of all the classes defined by this package.

ExpressionSet is a more formal class in the Biobase package.

exprs.MA

Extract Log-Expression Matrix from MAList

## **Description**

Extract the matrix of log-expression values from an MAList object.

## Usage

```
exprs.MA(MA)
```

## **Arguments**

MA

an MAList object.

#### **Details**

Converts M and A-values to log-expression values. The output matrix will have two columns for each array, in the order green, red for each array.

This contrasts with as .matrix.MAList which extracts the M-values only, or RG.MA which converts to expression values in RGList form.

## Value

A numeric matrix with twice the columns of the input.

## Author(s)

Gordon Smyth

## See Also

02. Classes gives an overview of data classes used in LIMMA.

fitFDist 43

E; LDD; -L	Moment Estimation	f Caalad E Distribution
fitFDist	Moment Estimation of	of Scaled F-Distribution

# Description

Moment estimation of the parameters of a scaled F-distribution given one of the degrees of freedom. This function is called internally by <code>ebayes</code> and is not usually called directly by a user.

# Usage

```
fitFDist(x,df1)
```

# Arguments

X	numeric vector or array of positive values representing a sample from an F-
	distribution.

df1 the first degrees of freedom of the F-distribution. May be an integer or a vector of the same length as x.

# Details

The function estimates scale and df2 under the assumption that x is distributed as scale times an F-distributed random variable on df1 and df2 degrees of freedom.

## Value

A list containing the components

scale scale factor for F-distribution

df2 the second degrees of freedom of the F-distribution

# Author(s)

Gordon Smyth

## See Also

```
ebayes, trigammaInverse
```

44 genas

fitted.MArrayLM

Fitted Values Method for MArrayLM Fits

## **Description**

Obtains fitted values from a fitted microarray linear model object.

# Usage

```
## S3 method for class 'MArrayLM':
fitted(object, design = object$design, ...)
```

# **Arguments**

```
object a fitted object of class inheriting from "MArrayLM".

design numeric design matrix.

further arguments passed to or from other methods.
```

### Value

A numeric matrix of fitted values.

## Author(s)

Gordon Smyth

## See Also

fitted

genas

Genuine Association of Gene Expression Profiles

### **Description**

Calculates biological correlation between two gene expression profiles.

## Usage

```
genas(fit, coef=c(1,2))
```

# Arguments

fit	an MArravLM fitted model object produced by lmFit or contrasts.fit
LIL	an mallavin nucu mouci obicci bioducca ov lini il di conclasts. Il c

and followed by eBayes

numeric vector of length 2 to indicate which contrasts/columns in the fit object

are to be used

genas 45

#### **Details**

The biological correlation between the true log fold changes of pairs of comparisons is computed. This method is to be applied when multiple groups (such as treatment groups, mutants or knockouts) are compared back to the same control group.

This method is an extension of the empirical Bayes method of limma. It aims to separate the technical correlation, which comes from comparing multiple treatment/mutant/knock-out groups to the same control group, from biological correlation, which is the true correlation of the gene expression profiles between two treatment/mutant/knock-out groups.

The fit object should include only differentially expressed genes. One approach is to calculate the true proportion of differentially expressed genes using convest on the F p-value produced by lmFit. Any reasonable set of genes displaying some degree of differential expression should be adequate.

### Value

## Author(s)

Belinda Phipson and Gordon Smyth

## See Also

```
lmFit, eBayes, contrasts.fit
```

## **Examples**

```
library(limma)
# Simulate gene expression data,
# 6 microarrays with 100 genes on each array
set.seed(2004)
y<-matrix(rnorm(600),ncol=6)

# two experimental groups and one control group with two replicates each
group<-factor(c("A","A","B","B","control","control"))
design<-model.matrix(~0+group)
colnames(design)<-c("A","B","control")

# fit a linear model
fit<-lmFit(y,design)
contrasts<-makeContrasts(A-control,B-control,levels=design)
fit2<-contrasts.fit(fit,contrasts)</pre>
```

46 geneSetTest

```
fit2<-eBayes(fit2)
# calculate biological correlation between the gene expression profiles of (A vs control)
genas(fit2)</pre>
```

geneSetTest Mean-rank Gene Set Test

# Description

Test whether a set of genes is highly ranked relative to other genes in terms of a given statistic. Genes are assumed to be independent.

# Usage

```
geneSetTest(selected, statistics, alternative="mixed", type="auto", ranks.only=TRUE,
wilcoxGST(selected, statistics, alternative="mixed")
barcodeplot(selected, statistics, labels=c("Up", "Down"),...)
barcodeplot2(selected, statistics, selected2=NULL, labels=c("Up", "Down"),...)
```

## **Arguments**

selected	index vector for the gene set. This can be a vector of indices, or a logical vector of the same length as statistics or, in general, any vector such that statistic[selected] gives the statistic values for the gene set to be tested.
selected2	index vector for a second gene set. Usually used to specify down-regulated genes when selected is used for up-regulated genes.vector specifying the elements of statistic in the test group.
statistics	numeric vector giving the values of the test statistic for every gene or probe in the reference set, usually every probe on the microarray.
alternative	character string specifying the alternative hypothesis, must be one of "mixed" (default), "either", "up" or "down". two.sided, "greater" and "less" are also permitted as synonyms for "either", "up" and "down" respectively.
type	character string specifying whether the statistics are t-like (" $t$ "), F-like " $f$ " or whether the function should make an educated guess ("auto")
ranks.only	logical, if TRUE only the ranks of the statistics are used.
nsim	$number \ of \ random \ samples \ to \ take \ in \ computing \ the \ p-value. \ Not \ used \ if \ \verb"ranks.only=TRUE".$
labels	character vector of length 2 of labels associated with large and small statistics respectively.
	other arguments are passed to plot.

geneSetTest 47

#### **Details**

wilcoxGST is a synonym for geneSetTest with ranks.only=TRUE and type="t". This test procedure was developed by Michaud et al (2008), who called it *mean-rank gene-set enrichment*.

These functions compute a p-value to test the hypothesis that the selected set of genes tends to be more highly ranked in terms of some test statistic compared to randomly selected genes. The statistic might be any statistic of interest, for example t-statistics or F-statistics for differential expression. They are competitive tests in the sense that genes in the test set are compared to other genes (Goeman and Buhlmann, 2007). By contrast, a self-contained gene set test such as roast tests for differential expression for the test set only without regard to other genes on the array. Like all gene set tests, these functions can be used to detect differential expression for a group of genes, even when the effects are too small or there is too little data to detect the genes individually. The also provides a means to compare the results between different experiments.

Because it is based on permuting genes, geneSetTest assumes that the different genes (or probes) are independent. (Strictly speaking, it assumes that the genes in the set are no more correlated on average than randomly selected genes.) This assumption may be reasonable if the gene set is relatively small and if there is relatively little genotypic variation in the data, for example if the data is obtained from genetically identical inbred mice. The independence assumption may be misleading if the gene set is large or if the data contains a lot of genotypic variation, for example for human cancer samples. These assumptions, when valid, permit a much quicker and more powerful significance test to be conducted.

The statistics are usually a set of probe-wise statistics arising for some comparison from a microarray experiment. They may be t-statistics, meaning that the genewise null hypotheses would be rejected for large positive or negative values, or they may be F-statistics, meaning that only large values are significant. Any set of signed statistics, such as log-ratios, M-values or moderated t-statistics, are treated as t-like. Any set of unsigned statistics, such as F-statistics, posterior probabilities or chi-square tests are treated as F-like. If type="auto" then the statistics will be taken to be t-like if they take both positive and negative values and will be taken to be F-like if they are all of the same sign.

There are four possible alternatives to test for. alternative=="up" means the genes in the set tend to be up-regulated, with positive t-statistics. alternative=="down" means the genes in the set tend to be down-regulated, with negative t-statistics. alternative=="either" means the set is either up or down-regulated as a whole. alternative=="mixed" test whether the genes in the set tend to be differentially expressed, without regard for direction. In this case, the test will be significant if the set contains mostly large test statistics, even if some are positive and some are negative.

The latter three alternatives are appropriate if you have a prior expection that all the genes in the set will react in the same direction. The "mixed" alternative is appropriate if you know only that the genes are involved in the relevant pathways, possibly in different directions. The "mixed" is the only meaningful alternative with F-like statistics.

The test statistic used for the gene-set-test is the mean of the statistics in the set. If ranks.only is TRUE the only the ranks of the statistics are used. In this case the p-value is obtained from a Wilcoxon test. If ranks.only is FALSE, then the p-value is obtained by simulation using nsim random selected sets of genes.

barcodeplot is a graphical representation of the Wilcox gene set test using ranks. barcodeplot2 displays two sets. It is particularly used to displaying directional sets, when there are separate sets of genes expected to go up and down respectively.

48 getEAWP

#### Value

geneSetTest and wilcoxGST return a numeric value giving the estimated p-value. barcodeplot and barcodeplot2 return no value but produce a plot as a side effect.

## Author(s)

Gordon Smyth and Di Wu

#### References

Goeman, JJ, and Buhlmann P (2007). Analyzing gene expression data in terms of gene sets: methodological issues. *Bioinformatics* 23, 980-987.

Michaud, J, Simpson, KM, Escher, R, Buchet-Poyau, K, Beissbarth, T, Carmichael, C, Ritchie, ME, Schutz, F, Cannon, P, Liu, M, Shen, X, Ito, Y, Raskind, WH, Horwitz, MS, Osato, M, Turner, DR, Speed, TP, Kavallaris, M, Smyth, GK, and Scott, HS (2008). Integrative analysis of RUNX1 downstream pathways and target genes. *BMC Genomics* 9, 363. http://www.biomedcentral.com/1471-2164/9/363

#### See Also

```
roast, romer, wilcox.test
```

An overview of tests in limma is given in 08. Tests.

### **Examples**

```
stat <- rnorm(100)
sel <- 1:10
wilcoxGST(sel,stat)
barcodeplot(sel,stat)
sel2 <- 11:20
barcodeplot2(sel,stat,sel2)</pre>
```

getEAWP

Extract Basic Data from Microarray Data Objects

### **Description**

Given a microarray data object of any known class, get the expression values, weights, probe annotation and A-values, which are needed for linear modelling. This function is called by the linear modelling functions in LIMMA.

# Usage

```
getEAWP(object)
```

# **Arguments**

object

a microarray data object. An object of class list, MAList, EList, marrayNorm, PLMset, vsn, or any class inheriting from ExpressionSet, or any object that can be coerced to a numeric matrix.

getLayout 49

#### **Details**

In the case of two-color objects, the Amean is computed from the matrix of A-values. For singlechannel objects, Amean is computed from the matrix of expression vales. PLMset, vsn and ExpressionSet are assumed to be single-channel for this purpose.

If object is a matrix, it is assumed to contain log-intensities if the values are all positive and log-ratios otherwise. Amean is computed in the former case but not the latter.

#### Value

#### A list with components

numeric matrix of log-ratios or log-intensities exprs

numeric matrix of weights weights data.frame of probe-annotation probes

Amean numeric vector of average log-expression for each probe

#### Author(s)

Gordon Smyth

### See Also

02. Classes gives an overview of data classes used in LIMMA.

Extract the Print Layout of an Array from the GAL File

# getLayout

## **Description**

From the Block, Row and Column information in a genelist, determine the number of grid rows and columns on the array and the number of spot rows and columns within each grid.

## Usage

```
getLayout(gal, guessdups=FALSE)
getLayout2(galfile)
getDupSpacing(ID)
```

# **Arguments**

gal	data.frame containing t	the GAL, i.e	e., giving the	e position and	gene identifier of
-----	-------------------------	--------------	----------------	----------------	--------------------

each spot

name or path of GAL file galfile

logical, if TRUE then try to determine number and spacing of duplicate spots, quessdups

i.e., within-array replicates

ID vector or factor of gene IDs 50 getSpacing

#### **Details**

A GenePix Array List (GAL) file is a list of genes and associated information produced by an Axon microarray scanner. The function <code>getLayout</code> determines the print layout from a data frame created from a GAL file or gene list. The data.frame must contain columns <code>Block</code>, <code>Column</code> and <code>Row</code>. (The number of tip columns is assumed to be either one or four.)

On some arrays, each probe may be duplicated a number of times (ndups) at regular intervals (spacing) in the GAL file. getDupSpacing determines valid values for ndups and spacing from a vector of IDs. If guessdups=TRUE, then getLayout calls getDupSpacing.

The function getLayout2 attempts to determine the print layout from the header information of an actual GAL file.

### Value

A printlayout object, which is a list with the following components. The last two components are present only if guessdups=TRUE.

ngrid.r	integer, number of grid rows on the arrays
ngrid.c	integer, number of grid columns on the arrays
nspot.r	integer, number of rows of spots in each grid
nspot.c	integer, number of columns of spots in each grid
ndups	integer, number of times each probe is printed on the array
spacing	integer, spacing between multiple printings of each probe

### Author(s)

Gordon Smyth and James Wettenhall

## See Also

An overview of LIMMA functions for reading data is given in 03.ReadingData.

# **Examples**

```
# gal <- readGAL()
# layout <- getLayout(gal)</pre>
```

getSpacing

Get Numerical Spacing

## **Description**

Convert character to numerical spacing measure for within-array replicate spots.

## Usage

```
getSpacing(spacing, layout)
```

gls.series 51

## **Arguments**

spacing	character string or integer. Acceptable character strings are "columns", "rows",
	"subarrays" or "topbottom". Integer values are simply passed through.
layout	list containing printer layout information

#### **Details**

"rows" means that duplicate spots are printed side-by-side by rows. These will be recorded in consecutive rows in the data object.

"columns" means that duplicate spots are printed side-by-sidy by columns. These will be separated in the data object by layout\$nspot.r rows.

"subarrays" means that a number of sub-arrays, with identical probes in the same arrangement, are printed on each array. The spacing therefore will be the size of a sub-array.

"topbottom" is the same as "subarrays" when there are two sub-arrays.

#### Value

Integer giving spacing between replicate spots in the gene list.

### Author(s)

Gordon Smyth

### See Also

An overview of LIMMA functions for reading data is given in 03.ReadingData.

# **Examples**

```
getSpacing("columns", list(ngrid.r=2, ngrid.c=2, nspot.r=20, nspot.c=19))
getSpacing("rows", list(ngrid.r=2, ngrid.c=2, nspot.r=20, nspot.c=19))
getSpacing("topbottom", list(ngrid.r=2, ngrid.c=2, nspot.r=20, nspot.c=19))
```

gls.series

Fit Linear Model to Microarray Data by Generalized Least Squares

## **Description**

Fit a linear model genewise to expression data from a series of microarrays. The fit is by generalized least squares allowing for correlation between duplicate spots or related arrays. This is a utility function for lmFit.

### Usage

```
gls.series(M, design=NULL, ndups=2, spacing=1, block=NULL, correlation=NULL, weights=N
```

52 gls.series

## **Arguments**

numeric matrix containing log-ratio or log-expression values for a series of mi-Μ croarrays, rows correspond to genes and columns to arrays. numeric design matrix defining the linear model, with rows corresponding to design arrays and columns to comparisons to be estimated. The number of rows must match the number of columns of M. Defaults to the unit vector meaning that the arrays are treated as replicates. ndups positive integer giving the number of times each gene is printed on an array. nrow (M) must be divisible by ndups. the spacing between the rows of M corresponding to duplicate spots, spacing=1 spacing for consecutive spots block vector or factor specifying a blocking variable on the arrays. Same length as ncol(M). correlation numeric value specifying the inter-duplicate or inter-block correlation. an optional numeric matrix of the same dimension as M containing weights for weights each spot. If it is of different dimension to M, it will be filled out to the same size. other optional arguments to be passed to dupcor. series.

#### **Details**

This is a utility function used by the higher level function lmFit. Most users should not use this function directly but should use lmFit instead.

This function is for fitting gene-wise linear models when some of the expression values are correlated. The correlated groups may arise from replicate spots on the same array (duplicate spots) or from a biological or technical replicate grouping of the arrays. This function is normally called by lmFit and is not normally called directly by users.

Note that the correlation is assumed to be constant across genes. If correlation=NULL then a call is made to duplicateCorrelation to estimated the correlation.

#### Value

### A list with components

coefficients numeric matrix containing the estimated coefficients for each linear model. Same number of rows as M, same number of columns as design.

stdev.unscaled

numeric matrix conformal with coef containing the unscaled standard deviations for the coefficient estimators. The standard errors are given by stdev.unscaled

\* sigma.

sigma numeric vector containing the residual standard deviation for each gene.

df.residual numeric vector giving the degrees of freedom corresponding to sigma

correlation inter-duplicate or inter-block correlation

QR decomposition of the generalized linear squares problem, i.e., the decom-

position of design standardized by the Choleski-root of the correlation matrix

defined by correlation

### Author(s)

qr

#### Gordon Smyth

gridr 53

### See Also

```
duplicateCorrelation.
```

An overview of linear model functions in limma is given by 06.LinearModels.

gridr

Row and Column Positions on Microarray

# Description

Grid and spot row and column positions.

# Usage

```
gridr(layout)
gridc(layout)
spotr(layout)
spotc(layout)
```

# **Arguments**

layout

list with the components ngrid.r, ngrid.c, nspot.r and nspot.c

## Value

Vector of length prod (unlist (layout)) giving the grid rows (gridr), grid columns (gridc), spot rows (spotr) or spot columns (spotc).

## Author(s)

Gordon Smyth

heatdiagram

Stemmed Heat Diagram

## **Description**

Creates a heat diagram showing the co-regulation of genes under one condition with a range of other conditions.

# Usage

heatDiagram(results,coef,primary=1,names=NULL,treatments=colnames(coef),limit=NUheatdiagram(stat,coef,primary=1,names=NULL,treatments=colnames(stat),critical.pr

54 heatdiagram

#### **Arguments**

TestResults matrix, containing elements -1, 0 or 1, from decideTests results numeric matrix of test statistics. Rows correspond to genes and columns to stat treatments or contrasts between treatments. coef numeric matrix of the same size as stat. Holds the coefficients to be displayed in the plot. number or name of the column to be compared to the others. Genes are included primary in the diagram according to this column of stat and are sorted according to this column of coef. If primary is a name, then stat and coef must have the same column names. optional character vector of gene names names optional character vector of treatment names treatments critical.primary critical value above which the test statistics for the primary column are considered significant and included in the plot critical.other critical value above which the other test statistics are considered significant. Should usually be no larger than critical.primary although larger values are permitted. limit optional value for coef above which values will be plotted in extreme color. Defaults to max (abs (coef)). "portrait" for upright plot or "landscape" for plot orientated to be orientation wider than high. "portrait" is likely to be appropriate for inclusion in printed document while "landscape" may be appropriate for a presentation on a computer screen. color associated with repressed gene regulation low color associated with induced gene regulation high number of distinct colors used for each of up and down regulation ncolors factor to increase or decrease size of column and row text cex numeric vector of length four giving the size of the margin widths. Default is mar cex\*c(5,6,1,1) for landscape and cex\*c(1,1,4,3) for portrait. any other arguments will be passed to the image function

#### **Details**

. . .

Users are encouraged to use heatDiagram rather than heatdiagram as the later function may be removed in future versions of limma.

This function plots an image of gene expression profiles in which rows (or columns for portrait orientation) correspond to treatment conditions and columns (or rows) correspond to genes. Only genes which are significantly differentially expressed in the primary condition are included. Genes are sorted by differential expression under the primary condition.

Note: the plot produced by this function is unique to the limma package. It should not be confused with "heatmaps" often used to display results from cluster analyses.

# Value

An image is created on the current graphics device. A matrix with named rows containing the coefficients used in the plot is also invisibly returned.

helpMethods 55

## Author(s)

Gordon Smyth

#### See Also

image.

## **Examples**

```
if(require("sma")) {
  data(MouseArray)

MA <- normalizeWithinArrays(mouse.data,layout=mouse.setup)
  design <- cbind(c(1,1,1,0,0,0),c(0,0,1,1,1))
  fit <- lmFit(MA,design=design)
  contrasts.mouse <- cbind(Control=c(1,0),Mutant=c(0,1),Difference=c(-1,1))
  fit <- eBayes(contrasts.fit(fit,contrasts=contrasts.mouse))
  results <- decideTests(fit,method="global",p=0.1)
  heatDiagram(results,fit$coef,primary="Difference")
}</pre>
```

helpMethods

Prompt for Method Help Topics

# Description

For any S4 generic function, find all methods defined in currently loaded packages. Prompt the user to choose one of these to display the help document.

# Usage

```
helpMethods(genericFunction)
```

## **Arguments**

```
genericFunction
```

a generic function or a character string giving the name of a generic function

## Author(s)

Gordon Smyth

## See Also

showMethods

## **Examples**

```
## Not run: helpMethods(show)
```

imageplot3by2

|--|

# Description

Write imageplots to files in PNG format, six plots to a file in a 3 by 2 grid arrangement.

## Usage

```
imageplot3by2(RG, z="Gb", prefix=paste("image",z,sep="-"), path=NULL, zlim=NULL,
```

# Arguments

RG	an RGList or MAList object, or any list with component named by $\boldsymbol{z}$
Z	character string giving name of component of RG to plot
prefix	character string giving prefix to attach to file names
path	character string specifying directory for output files
zlim	numeric vector of length 2, giving limits of response vector to be associated with saturated colors
common.lim	logical, should all plots on a page use the same axis limits
	any other arguments are passed to imageplot

# **Details**

At the time of writing, this function writes plots in PNG format in an arrangement optimized for A4-sized paper.

## Value

No value is returned, but one or more files are written to the working directory. The number of files is determined by the number of columns of RG.

## Author(s)

Gordon Smyth

# See Also

An overview of diagnostic functions available in LIMMA is given in 09.Diagnostics.

imageplot 57

imageplot Image Plot of Microarray Statistics	imageplot	Image Plot of Microarray Statistics	
---	-----------	-------------------------------------	--

## **Description**

Creates an image of colors or shades of gray that represent the values of a statistic for each spot on a spotted microarray. This function can be used to explore any spatial effects across the microarray.

### Usage

```
imageplot(z, layout, low = NULL, high = NULL, ncolors = 123, zerocenter = NULL,
zlim = NULL, mar=c(2,1,1,1), legend=TRUE, ...)
```

## **Arguments**

Z	numeric vector or array. This vector can contain any spot statistics, such as log intensity ratios, spot sizes or shapes, or t-statistics. Missing values are allowed and will result in blank spots on the image. Infinite values are not allowed.
layout	a list specifying the dimensions of the spot matrix and the grid matrix.
low	color associated with low values of z. May be specified as a character string such as "green", "white" etc, or as a rgb vector in which $c(1,0,0)$ is red, $c(0,1,0)$ is green and $c(0,0,1)$ is blue. The default value is "green" if zerocenter=T or "white" if zerocenter=F.
high	color associated with high values of z. The default value is "red" if $zerocenter=T$ or "blue" if $zerocenter=F$ .
ncolors	number of color shades used in the image including low and high.
zerocenter	should zero values of z correspond to a shade exactly halfway between the colors low and high? The default is TRUE if z takes positive and negative values, otherwise FALSE.
zlim	numerical vector of length 2 giving the extreme values of $z$ to associate with colors low and high. By default zlim is the range of $z$ . Any values of $z$ outside the interval zlim will be truncated to the relevant limit.
mar	numeric vector of length 4 specifying the width of the margin around the plot. This argument is passed to par.
legend	logical, if TRUE the range of $z$ and $zlim$ is shown in the bottom margin
	any other arguments will be passed to the function image

## **Details**

This function may be used to plot the values of any spot-specific statistic, such as the log intensity ratio, background intensity or a quality measure such as spot size or shape. The image follows the layout of an actual microarray slide with the bottom left corner representing the spot (1,1,1,1). The color range is used to represent the range of values for the statistic. When this function is used to plot the red/green log-ratios, it is intended to be an in silico version of the classic false-colored red-yellow-green image of a scanned two-color microarray.

This function is related to the earlier plot. spatial function in the sma package and to the later maImage function in the marray package. It differs from plot. spatial most noticeably in that all the spots are plotted and the image is plotted from bottom left rather than from top left. It is

58 intraspotCorrelation

intended to display spatial patterns and artefacts rather than to highlight only the extreme values as does plot.spatial. It differs from maImage in that any statistic may be plotted and in its use of a red-yellow-green color scheme for log-ratios, similar to the classic false-colored jpeg image, rather than the red-black-green color scheme associated with heat maps.

### Value

An plot is created on the current graphics device.

### Author(s)

Gordon Smyth

#### See Also

```
maImage, image.
```

An overview of diagnostic functions available in LIMMA is given in 09. Diagnostics.

## **Examples**

```
 \begin{tabular}{ll} $M < - rnorm(8*4*16*16)$ \\ imageplot(M,layout=list(ngrid.r=8,ngrid.c=4,nspot.r=16,nspot.c=16)) \end{tabular}
```

```
intraspotCorrelation
```

Intra-Spot Correlation for Two Color Data

### **Description**

Estimate the within-block correlation associated with spots for spotted two color microarray data.

## Usage

```
intraspotCorrelation(object, design, trim=0.15)
```

## **Arguments**

object	an MAList object or a list from which M and A values may be extracted
design	a numeric matrix containing the design matrix for linear model in terms of the individual channels. The number of rows should be twice the number of arrays. The number of columns will determine the number of coefficients estimated for each gene.
trim	the fraction of observations to be trimmed from each end of the atanh-correlations when computing the consensus correlation. See mean.

## Details

This function estimates the correlation between two channels observed on each spot. The correlation is estimated by fitting a heteroscedastic regression model to the M and A-values of each gene. The function also returns a consensus correlation, which is a robust average of the individual correlations, which can be used as input for functions lmscFit.

The function may take long time to execute.

is.fullrank 59

#### Value

## A list with components

```
consensus.correlation
```

robust average of the estimated inter-duplicate correlations. The average is the trimmed mean of the correlations for individual genes on the atanh-transformed scale.

atanh.correlations

a numeric vector giving the individual genewise correlations on the atanh scale

df

numeric matrix of degrees of freedom associated with the correlations. The first column gives the degrees of freedom for estimating the within-spot or M-value mean square while the second gives the degrees of freedom for estimating the between spot or A-value mean square.

## Author(s)

Gordon Smyth

### References

Smyth, G. K. (2005). Individual channel analysis of two-colour microarray data. *Proceedings of the 55th Session of the International Statistics Institute*, 5-12 April 2005, Sydney, Australia, Paper 116. http://www.statsci.org/smyth/pubs/ISI2005-116.pdf

#### See Also

This function uses remlscore from the statmod package.

An overview of methods for single channel analysis in limma is given by 07. Single Channel.

## **Examples**

```
# See lmscFit
## Not run:
corfit <- intraspotCorrelation(MA, design)
all.correlations <- tanh(corfit$atanh.correlations)
boxplot(all.correlations)
## End(Not run)</pre>
```

is.fullrank

Check for Full Column Rank

# Description

Test whether a numeric matrix has full column rank.

# Usage

```
is.fullrank(x)
nonEstimable(x)
```

60 isNumeric

## **Arguments**

Х

a numeric matrix or vector

#### **Details**

is.fullrank is used to check the integrity of design matrices in limma, for example after subsetting operations.

nonEstimable is used by lmFit to report which coefficients in a linear model cannot be estimated.

## Value

```
is.fullrank returns TRUE or FALSE.
```

nonEstimable returns a character vector of names for the columns of x which are linearly dependent on previous columns. If x has full column rank, then the value is NULL.

### Author(s)

Gordon Smyth

## **Examples**

```
# TRUE
is.fullrank(1)
is.fullrank(cbind(1,0:1))

# FALSE
is.fullrank(0)
is.fullrank(matrix(1,2,2))
nonEstimable(matrix(1,2,2))
```

isNumeric

Test for Numeric Argument

## **Description**

Test whether argument is numeric or a data.frame with numeric columns.

# Usage

```
isNumeric(x)
```

### **Arguments**

Х

any object

## **Details**

This function is used to check the validity of arguments for numeric functions. It is an attempt to emulate the behavior of internal generic math functions.

is Numeric differs from is . numeric in that data.frames with all columns numeric are accepted as numeric.

kooperberg 61

### Value

TRUE or FALSE

### Author(s)

Gordon Smyth

### See Also

```
is.numeric, Math
```

### **Examples**

```
isNumeric(3)
isNumeric("a")
x <- data.frame(a=c(1,1),b=c(0,1))
isNumeric(x)  # TRUE
is.numeric(x)  # FALSE</pre>
```

kooperberg

Kooperberg Model-Based Background Correction

## **Description**

This function uses a Bayesian model to background correct GenePix microarray data.

### Usage

```
kooperberg(RG, a=TRUE, layout=RG$printer, verbose=TRUE)
```

## **Arguments**

RG	an RGList of GenePix data, read in using read.maimages, with other.columns=c("F635 SD", "B635 SD", "F532 SD", "B532 SD", "B532 Mean", "B635 Mean", "FPixels", "BPixels").
a	logical. If TRUE, the 'a' parameters in the model (equation 3 and 4) are estimated for each slide. If FALSE the 'a' parameters are set to unity.
layout	list containing print layout with components ngrid.r, ngrid.c, nspot.r

and nspot.c. Defaults to RG\$printer.

verbose logical. If TRUE, progress is reported to standard output.

### **Details**

This function is for use with GenePix data and is designed to cope with the problem of large numbers of negative intensities and hence missing values on the log-intensity scale. It avoids missing values in most cases and at the same time dampens down the variability of log-ratios for low intensity spots. See Kooperberg et al (2002) for more details.

kooperberg uses the foreground and background intensities, standard deviations and number of pixels to compute empirical estimates of the model parameters as described in equation 2 of Kooperberg et al (2002).

#### Value

An RGList containing the components

R matrix containing the background adjusted intensities for the red channel for

each spot for each array

G matrix containing the background adjusted intensities for the green channel for

each spot for each array

printer list containing print layout

### Author(s)

Matthew Ritchie

### References

Kooperberg, C., Fazzio, T. G., Delrow, J. J., and Tsukiyama, T. (2002) Improved background correction for spotted DNA microarrays. *Journal of Computational Biology* **9**, 55-66.

#### See Also

04.Background gives an overview of background correction functions defined in the LIMMA package.

## **Examples**

```
# This is example code for reading and background correcting GenePix data
# given GenePix Results (gpr) files in the working directory (data not
# provided).
## Not run:
genepixFiles <- dir(pattern="*\\.gpr$") # get the names of the GenePix image analysis out
RG <- read.maimages(genepixFiles, source="genepix", other.columns=c("F635 SD","B635 SD",'
RGmodel <- kooperberg(RG)
MA <- normalizeWithinArrays(RGmodel)
## End(Not run)</pre>
```

LargeDataObject-class

Large Data Object - class

# **Description**

A virtual class including the data classes RGList, MAList and MArrayLM, all of which typically contain large quantities of numerical data in vector, matrices and data.frames.

# Methods

A show method is defined for objects of class LargeDataObject which uses printHead to print only the leading elements or rows of components or slots which contain large quantities of data.

limmaUsersGuide 63

#### Author(s)

Gordon Smyth

#### See Also

02. Classes gives an overview of all the classes defined by this package.

### **Examples**

# see normalizeBetweenArrays

limmaUsersGuide

View Limma User's Guide

## **Description**

Finds the location of the Limma User's Guide and optionally opens it.

## Usage

limmaUsersGuide(view=TRUE)

## **Arguments**

view

logical, should the document be opened using the default PDF document reader?

## **Details**

The function vignette ("limma") will find the short limma Vignette which describes how to obtain the Limma User's Guide. The User's Guide is not itself a true vignette because it is not automatically generated using Sweave during the package build process. This means that it cannot be found using vignette, hence the need for this special function.

If the operating system is other than Windows, then the PDF viewer used is that given by Sys.getenv("R\_PDFVIEWE The PDF viewer can be changed using Sys.putenv(R\_PDFVIEWER=).

This function is used by drop-down Vignettes menu when the Rgui interface for Windows is used.

# Value

Character string giving the file location.

### Author(s)

Gordon Smyth

#### See Also

```
vignette, openPDF, openVignette, Sys.getenv, Sys.putenv
```

# **Examples**

limmaUsersGuide(view=FALSE)

64 ImFit

lmFit	Linear Model for Series of Arrays	
THIE I C	Linear Model for Series of Arrays	

## **Description**

Fit linear model for each gene given a series of arrays

#### Usage

lmFit(object,design=NULL,ndups=1,spacing=1,block=NULL,correlation,weights=NULL,m

#### **Arguments**

object	object of class numeric, matrix, MAList, EList, marrayNorm, ExpressionSet or PLMset containing log-ratios or log-values of expression for a series of microarrays ${\bf P}_{\bf p}$
design	the design matrix of the microarray experiment, with rows corresponding to arrays and columns to coefficients to be estimated. Defaults to the unit vector meaning that the arrays are treated as replicates.
ndups	positive integer giving the number of times each gene is printed on an array
spacing	positive integer giving the spacing between duplicate spots, spacing=1 for consecutive spots
block	vector or factor specifying a blocking variable on the arrays. Has length equal to the number of arrays. Must be NULL if ndups>2.
correlation	the inter-duplicate or inter-technical replicate correlation
weights	optional numeric matrix containing weights for each spot
method	character string, "ls" for least squares or "robust" for robust regression
	other optional arguments to be passed to lm.series, gls.series or mrlm

## **Details**

This function fits multiple linear models. It accepts data from a experiment involving a series of microarrays with the same set of probes. A linear model is fitted to the expression data for each probe. The expression data should be log-ratios for two-color array platforms or log-expression values for one-channel platforms. (To fit linear models to the individual channels of two-color array data, see lmscFit.) The coefficients of the fitted models describe the differences between the RNA sources hybridized to the arrays. The probe-wise fitted model results are stored in a compact form suitable for further processing by other functions in the limma package.

The function allows for missing values and accepts quantitative weights through the weights argument. It also supports two different correlation structures. If block is not NULL then different arrays are assumed to be correlated. If block is NULL and ndups is greater than one then replicate spots on the same array are assumed to be correlated. It is not possible at this time to fit models with both a block structure and a duplicate-spot correlation structure simultaneously.

If object is a matrix then it should contain log-ratios or log-expression data with rows corresponding to probes and columns to arrays. (A numeric vector is treated the same as a matrix with one column.) For objects of other classes, a matrix of expression values is taken from the appropriate component or slot of the object. If object is of class MAList or marrayNorm, then the matrix of log-ratios (M-values) is extracted. If object is of class ExpressionSet, then the

ImFit 65

expression matrix is extracted. (This may contain log-expression or log-ratio values, depending on the platform.) If object is of class PLMset then the matrix of chip coefficients chip.coefs is extracted.

The arguments design, ndups, spacing and weights will be extracted from the data object if available and do not normally need to set explicitly in the call. On the other hand, if any of these are set in the function call then they will over-ride the slots or components in the data object. If object is an PLMset, then weights are computed as 1/pmax (object@se.chip.coefs, 1e-05)^2. If object is an ExpressionSet object, then weights are not computed.

If the argument block is used, then it is assumed that ndups=1.

The correlation argument has a default value of 0.75, but in normal use this default value should not be relied on and the correlation value should be estimated using the function duplicateCorrelation. The default value is likely to be too high in particular if used with the block argument.

The actual linear model computations are done by passing the data to one the lower-level functions lm.series, gls.series or mrlm. The function mrlm is used if method="robust". If method="ls", then gls.series is used if a correlation structure has been specified, i.e., if ndups>1 or block is non-null and correlation is different from zero. If method="ls" and there is no correlation structure, lm.series is used.

#### Value

Object of class MArrayLM

#### Author(s)

Gordon Smyth

#### See Also

An overview of linear model functions in limma is given by 06.LinearModels.

### **Examples**

```
# Simulate gene expression data for 100 probes and 6 microarrays
# Microarray are in two groups
# First two probes are differentially expressed in second group
\# Std deviations vary between genes with prior df=4
sd \leftarrow 0.3*sqrt(4/rchisq(100,df=4))
y \leftarrow matrix(rnorm(100*6, sd=sd), 100, 6)
rownames(y) <- paste("Gene",1:100)</pre>
y[1:2,4:6] \leftarrow y[1:2,4:6] + 2
design <- cbind(Grp1=1, Grp2vs1=c(0,0,0,1,1,1))
options(digit=3)
# Ordinary fit
fit <- lmFit(y,design)</pre>
fit <- eBayes(fit)</pre>
fit
as.data.frame(fit[1:10,2])
# Various ways of summarising or plotting the results
topTable(fit,coef=2)
qqt(fit$t[,2],df=fit$df.residual+fit$df.prior)
abline(0,1)
volcanoplot(fit,coef=2,highlight=2)
```

66 lmscFit

```
# Various ways of writing results to file
## Not run: write.fit(fit,file="exampleresults.txt")
## Not run: write.table(fit,file="exampleresults2.txt")
# Robust fit
# (There may be some warning messages)
fit2 <- lmFit(y,design,method="robust")</pre>
# Fit with correlated arrays
# Suppose each pair of arrays is a block
block \leftarrow c(1,1,2,2,3,3)
dupcor <- duplicateCorrelation(y,design,block=block)</pre>
dupcor$consensus.correlation
fit3 <- lmFit(y,design,block=block,correlation=dupcor$consensus)</pre>
# Fit with duplicate probes
# Suppose two side-by-side duplicates of each gene
rownames(y) <- paste("Gene", rep(1:50, each=2))</pre>
dupcor <- duplicateCorrelation(y, design, ndups=2)</pre>
dupcor$consensus.correlation
fit4 <- lmFit(y,design,ndups=2,correlation=dupcor$consensus)</pre>
fit4 <- eBayes(fit3)
dim(fit4)
topTable(fit4,coef=2)
# Fold-change thresholding
fit <- lmFit(y,design)</pre>
fit <- treat(fit,lfc=0.1)</pre>
topTreat(fit,coef=2)
```

lmscFit

Fit Linear Model to Individual Channels of Two-Color Data

# Description

Fit a linear model to the individual log-intensities for each gene given a series of two-color arrays

### Usage

```
lmscFit(object, design, correlation)
```

# Arguments

object an MAList object or a list from which M and A values may be extracted

design a numeric matrix containing the design matrix for linear model in terms of the individual channels. The number of rows should be twice the number of arrays. The number of columns will determine the number of coefficients estimated for each gene.

correlation numeric value giving the intra-spot correlation

ImscFit 67

#### **Details**

For two color arrays, the channels measured on the same set of arrays are correlated. The M and A however are uncorrelated for each gene. This function fits a linear model to the set of M and A-values for each gene after re-scaling the M and A-values to have equal variances. The input correlation determines the scaling required. The input correlation is usually estimated using intraspotCorrelation before using lmscFit.

Missing values in M or A are not allowed.

#### Value

An object of class MArrayLM

#### Author(s)

Gordon Smyth

#### References

Smyth, G. K. (2005). Individual channel analysis of two-colour microarray data. *Proceedings of the 55th Session of the International Statistics Institute*, 5-12 April 2005, Sydney, Australia, Paper 116. http://www.statsci.org/smyth/pubs/ISI2005-116.pdf

#### See Also

```
lm.fit.
```

An overview of methods for single channel analysis in limma is given by 07. Single Channel.

## **Examples**

```
if(require("sma")) {
# Subset of data from ApoAI case study in Limma User's Guide
data(MouseArray)
# Avoid non-positive intensities
RG <- backgroundCorrect (mouse.data, method="normexp")</pre>
MA <- normalizeWithinArrays(RG, mouse.setup)
MA <- normalizeBetweenArrays (MA, method="Aq")
# Randomly choose 500 genes for this example
i <- sample(1:nrow(MA),500)</pre>
MA <- MA[i,]
targets <- data.frame(Cy3=I(rep("Pool",6)),Cy5=I(c("WT","WT","WT","KO","KO","KO")))
targets.sc <- targetsA2C(targets)</pre>
targets.sc$Target <- factor(targets.sc$Target,levels=c("Pool","WT","KO"))</pre>
design <- model.matrix(~Target,data=targets.sc)</pre>
corfit <- intraspotCorrelation(MA, design)</pre>
fit <- lmscFit (MA, design, correlation=corfit$consensus)</pre>
cont.matrix <- cbind(KOvsWT=c(0,-1,1))</pre>
fit2 <- contrasts.fit(fit,cont.matrix)</pre>
fit2 <- eBayes(fit2)</pre>
topTable(fit2,adjust="fdr")
```

68 lm.series

lm.series	Fit Linear Model to Microrray Data by Ordinary Least Squares

## **Description**

Fit a linear model genewise to expression data from a series of arrays. This function uses ordinary least squares and is a utility function for lmFit.

## Usage

```
lm.series(M,design=NULL,ndups=1,spacing=1,weights=NULL)
```

## **Arguments**

М	numeric matrix containing log-ratio or log-expression values for a series of mi- croarrays, rows correspond to genes and columns to arrays
design	numeric design matrix defining the linear model. The number of rows should agree with the number of columns of M. The number of columns will determine the number of coefficients estimated for each gene.
ndups	number of duplicate spots. Each gene is printed ndups times in adjacent spots on each array.
spacing	the spacing between the rows of M corresponding to duplicate spots, ${\tt spacing=1}$ for consecutive spots
weights	an optional numeric matrix of the same dimension as M containing weights for each spot. If it is of different dimension to M, it will be filled out to the same size.

## **Details**

This is a utility function used by the higher level function lmFit. Most users should not use this function directly but should use lmFit instead.

The linear model is fit for each gene by calling the function lm.fit or lm.wfit from the base library.

### Value

qr

# A list with components

coefficients numeric matrix containing the estimated coefficients for each linear model. Same number of rows as M, same number of columns as design.

stdev.unscaled

numeric matrix conformal with coef containing the unscaled standard deviations for the coefficient estimators. The standard errors are given by stdev.unscaled \* sigma.

sigma

numeric vector containing the residual standard deviation for each gene.

df.residual

numeric vector giving the degrees of freedom corresponding to sigma.

QR-decomposition of design

loessFit 69

#### Author(s)

Gordon Smyth

#### See Also

lm.fit.

An overview of linear model functions in limma is given by 06.LinearModels.

loessFit Fast Simple Loess

## **Description**

A fast version of locally weighted regression when there is only one x-variable and only the fitted values and residuals are required.

## Usage

```
loessFit(y, x, weights=NULL, span=0.3, bin=0.01/(2-is.null(weights)), iterations
```

### **Arguments**

У	numeric vector of response values. Missing values are allowed.
X	numeric vector of predictor values Missing values are allowed.
weights	numeric vector of non-negative weights. Missing values are allowed.
span	numeric parameter between 0 and 1 specifying proportion of data to be used in the local regression moving window. Larger numbers give smoother fits.
bin	numeric value between 0 and 1 giving the proportion of the data which can be grouped in a single bin when doing local regression fit. bin=0 forces an exact local regression fit with no interpolation.
iterations	number of iterations of loess fit

#### **Details**

This is a wrapper function to the Fortran and C code in the stats package which underlies the lowess and loess functions. Its is to give a streamlined common interface to lowess and loess for use in normalizeWithinArrays. When weights is null, this function is in effect a call to lowess in the stats package, with appropropriate choice of tuning parameters. When weights is non-null, it is in effect a call to loess. See the help pages for those functions for references and credits.

Note tha lowess is faster, needs less memory and is able to use a more accurate interpolation scheme than loess, so it is desirable to use lowess whenever loess is not needed to handle quantitative weights.

The arguments span, cell and iterations here have the same meaning as in loess. span is equivalent to the argument f of lowess and iterations is equivalent to iter+1. The parameter bin is intended to give a simple uniform interface to the delta argument of lowess and the cell argument of loess. bin translates to delta=bin\*diff(range(x)) in a call to lowess or to cell=bin/span in a call to loess.

70 ma3x3

Unlike lowess, loessFit returns values in original rather than sorted order. Also unlike lowess, loessFit allows missing values, the treatment being analogous to na.exclude. Unlike loess, loessFit returns a linear regression fit if there are insufficient observations to estimate the loess curve.

#### Value

### A list with components

```
fitted numeric vector of same length as y giving the loess fit residuals numeric vector of same length as x giving residuals from the fit
```

### Author(s)

Gordon Smyth, based on code from lowess and loess by BD Ripley

### See Also

```
See lowess and loess in the stats package.
```

See 05.Normalization for an outline of the limma package normalization functions.

### **Examples**

```
y <- rnorm(1000) 
 x <- rnorm(1000) 
 w <- rep(1,1000) 
 \# The following are equivalent apart from execution time 
 \# and interpolation inaccuracies 
 \$ system.time(fit <- loessFit(y,x)$fitted) 
 \$ system.time(fit <- loessFit(y,x,w)$fitted) 
 \$ system.time(fit <- fitted(loess(y~x,weights=w,span=0.3,family="symmetric",iterations=4))) 
 \# The same but with sorted x-values 
 \$ system.time(fit <- lowess(x,y,f=0.3)$y)
```

ma3x3

Two dimensional Moving Averages with 3x3 Window

# Description

Apply a specified function to each to each value of a matrix and its immediate neighbors.

### Usage

```
ma3x3.matrix(x,FUN=mean,na.rm=TRUE,...)
ma3x3.spottedarray(x,printer,FUN=mean,na.rm=TRUE,...)
```

# Arguments

X	numeric matrix
FUN	function to apply to each window of values
na.rm	logical value, should missing values be removed when applying FUN
	other arguments are passed to FUN
printer	list giving the printer layout, see PrintLayout-class

makeContrasts 71

#### **Details**

For ma3x3.matrix, x is an arbitrary function. for ma3x3.spotted, each column of x is assumed to contain the expression values of a spotted array in standard order. The printer layout information is used to re-arrange the values of each column as a spatial matrix before applying ma3x3.matrix.

### Value

Numeric matrix of same dimension as x containing smoothed values

### Author(s)

Gordon Smyth

#### See Also

An overview of functions for background correction are given in 04.Background.

# **Examples**

```
x \leftarrow matrix(c(2,5,3,1,6,3,10,12,4,6,4,8,2,1,9,0),4,4) ma3x3.matrix(x,FUN="mean") ma3x3.matrix(x,FUN="min")
```

makeContrasts

Construct Matrix of Custom Contrasts

# **Description**

Construct the contrast matrix corresponding to specified contrasts of a set of parameters.

# Usage

```
makeContrasts(..., contrasts=NULL, levels)
```

# **Arguments**

expressions, or character strings which can be parsed to expressions, specifying contrasts
 contrasts
 character vector specifying contrasts

levels character vector or factor giving the names of the parameters of which contrasts are desired, or a design matrix or other object with the parameter names as

are desired, or a design matrix or other object with the pa

column names.

72 makeUnique

#### **Details**

This function expresses contrasts between a set of parameters as a numeric matrix. The parameters are usually the coefficients from a linear model fit, so the matrix specifies which comparisons between the coefficients are to be extracted from the fit. The output from this function is usually used as input to contrasts.fit. The contrasts can be specified either as expressions using . . . or as a character vector through contrasts. (Trying to specify contrasts both ways will cause an error.)

The parameter names must be syntactically valid variable names in R and so, for example, must begin with a letter rather than a numeral. See make.names for a complete specification of what is a valid name.

#### Value

Matrix which columns corresponding to contrasts.

## Author(s)

Gordon Smyth

#### See Also

An overview of linear model functions in limma is given by the help page 06.LinearModels.

## **Examples**

```
\label{eq:makeContrasts} $$ $ (B-A, C-B, C-A, levels=c("A", "B", "C")) $$ $ makeContrasts(contrasts="A-(B+C)/2", levels=c("A", "B", "C")) $$ $ x <- c("B-A", "C-B", "C-A") $$ $ makeContrasts(contrasts=x, levels=c("A", "B", "C")) $$
```

makeUnique

Make Values of Character Vector Unique

# Description

Paste characters on to values of a character vector to make them unique.

#### Usage

```
makeUnique(x)
```

# **Arguments**

Х

object to be coerced to a character vector

## **Details**

Repeat values of x are labelled with suffixes "1", "2" etc.

#### Value

A character vector of the same length as x

MAList-class 73

#### Author(s)

Gordon Smyth

#### See Also

makeUnique is called by merge.RGList. Compare with make.unique in the base package.

#### **Examples**

```
x <- c("a","a","b")
makeUnique(x)</pre>
```

MAList-class

M-value, A-value Expression List - class

#### **Description**

A simple list-based class for storing M-values and A-values for a batch of spotted microarrays. MAList objects are usually created during normalization by the functions normalizeWithinArrays or MA.RG.

#### **Slots/List Components**

MAList objects can be created by new ("MAList", MA) where MA is a list. This class contains no slots (other than .Data), but objects should contain the following components:

m: numeric matrix containing the M-values (log-2 expression ratios). Rows correspond to spots and columns to arrays.

A: numeric matrix containing the A-values (average log-2 expression values).

# Optional components include:

other: list containing other matrices, all of the same dimensions as M.
genes: data.frame containing probe information. Should have one row for each spot. May have any number of colutargets: data.frame containing information on the target RNA samples. Rows correspond to arrays. May have any number of columns of the same dimensions as M.

numeric matrix of same dimensions as M containing relative spot quality weights. Elements should be non-n

printer: list containing information on the process used to print the spots on the arrays. See PrintLayout.

Valid MAList objects may contain other optional components, but all probe or array information should be contained in the above components.

#### Methods

weights:

This class inherits directly from class list so any operation appropriate for lists will work on objects of this class. In addition, MAList objects can be subsetted and combined. RGList objects will return dimensions and hence functions such as dim, nrow and ncol are defined. MALists also inherit a show method from the virtual class LargeDataObject, which means that RGLists will print in a compact way.

Other functions in LIMMA which operate on MAList objects include normalizeWithinArrays, normalizeBetweenArrays, normalizeForPrintorder, plotMA and plotPrintTipLoess.

74 MArrayLM-class

### Author(s)

Gordon Smyth

#### See Also

02. Classes gives an overview of all the classes defined by this package.

marrayNorm is the corresponding class in the marray package.

MArrayLM-class

Microarray Linear Model Fit - class

# Description

A list-based class for storing the results of fitting gene-wise linear models to a batch of microarrays. Objects are normally created by lmFit.

#### Slots/Components

 ${\tt MArrayLM}$  objects do not contain any slots (apart from .Data) but they should contain the following list components:

coefficients: matrix containing fitted coefficients or contrasts

stdev.unscaled: matrix containing unscaled standard deviations of the coefficients or contrasts

sigma: numeric vector containing residual standard deviations for each gene

df.residual: numeric vector containing residual degrees of freedom for each gene

Objects may also contain the following optional components:

Amean: numeric vector containing the average log-intensity for each probe over all the arrays in the original linear model fit. Note this vector does not change when a contrast is applied to the fit using contrasts.fit.

genes: data.frame containing gene names and annotation

design: design matrix of full column rank

contrasts: matrix defining contrasts of coefficients for which results are desired

F: numeric vector giving moderated F-statistics for testing all contrasts equal to zero

F.p. value: numeric vector giving p-value corresponding to F. stat

s2.prior: numeric value giving empirical Bayes estimated prior value for residual variances

df.prior: numeric vector giving empirical Bayes estimated degrees of freedom associated with s2.prior for each gene

s2.post: numeric vector giving posterior residual variances

t: matrix containing empirical Bayes t-statistics

var.prior: numeric vector giving empirical Bayes estimated prior variance for each true coefficient

cov.coefficients: numeric matrix giving the unscaled covariance matrix of the estimable coefficients mdplot 75

pivot: integer vector giving the order of coefficients in cov.coefficients. Is computed by the QR-decomposition of the design matrix.

If there are no weights and no missing values, then the MArrayLM objects returned by lmFit will also contain the QR-decomposition of the design matrix, and any other components returned by lm.fit.

### Methods

RGList objects will return dimensions and hence functions such as dim, nrow and ncol are defined. MArrayLM objects inherit a show method from the virtual class LargeDataObject.

The functions ebayes and classifyTestsF accept MArrayLM objects as arguments.

#### Author(s)

Gordon Smyth

#### See Also

02. Classes gives an overview of all the classes defined by this package.

mdplot

mdplot

### **Description**

Creates a mean-difference plot.

### Usage

```
mdplot(x, ...)
```

### **Arguments**

x numeric matrix with at least two columns

... any other arguments are passed to plot

#### **Details**

Plots differences vs means for a set of bivariate values. This is useful to contrast expression values for two microarrays.

Note that an MA-plot plotMA is a type of mean-difference plot.

### Value

A plot is created on the current graphics device.

### Author(s)

Gordon Smyth

76 merge

#### References

Chambers, J. M., Cleveland, W. S., Kleiner, B., and Tukey, P. A. (1983). Graphical Methods of Data Analysis. Wadsworth (pp. 48-57).

Cleveland, W. S., (1993). Visualizing Data. Hobart Press.

Bland, J. M., and Altman, D. G. (1986). Statistical methods for assessing agreement between two methods of clinical measurement. Lancet i, 307-310.

See also http://www.statsci.org/micrarra/refs/maplots.html

#### See Also

An overview of diagnostic functions available in LIMMA is given in 09. Diagnostics.

merge

Merge RGList or MAList Data Objects

### **Description**

Merge two microarray data sets represented by RGLists in possibly irregular order.

#### Usage

```
## S3 method for class 'RGList':
merge(x,y,...)
```

### **Arguments**

. . .

x RGList-class or MAList-class object

y RGList object, corresponding to the same genes as for x, possibly in a different order, but with different arrays.

other arguments are accepted but not used at present

### **Details**

RGList and MAList objects are list objects containing numeric matrices all of the same dimensions. The RGLists are merged by merging each of the components by row names or, if there are no row names, by IDs in the genes component. Unlike when using cbind, row names are not required to be in the same order or to be unique. In the case of repeated row names, the order of the rows with repeated names in preserved. This means that the first occurrence of each name in x is matched with the first occurrence of the same name in y. The second with the second, and so on. The final vector of row names is the same as in x.

Note: if the RGList objects contain the same number of genes in the same order then the appropriate function to combine them is cbind rather than merge.

### Value

An merged object of the same class as x and y with the same components as x. Component matrices have the same rows names as in x but columns from y as well as x.

mergeScans 77

#### Author(s)

Gordon Smyth

#### See Also

R base provides a merge method for merging data.frames.

An overview of limma commands for reading, subsetting and merging data is given in 03.Reading-Data.

### **Examples**

```
M <- A <- matrix(11:14,4,2)
rownames(M) <- rownames(A) <- c("a","a","b","c")
MA1 <- new("MAList",list(M=M,A=A))

M <- A <- matrix(21:24,4,2)
rownames(M) <- rownames(A) <- c("b","a","a","c")
MA2 <- new("MAList",list(M=M,A=A))

merge(MA1,MA2)
merge(MA2,MA1)</pre>
```

mergeScans

Merge two scans of two-color arrays

#### **Description**

Merge two sets of intensities of two-color arrays that are scanned twice at two different scanner settings, one at a lower gain setting with no saturated spot intensities and the other at a higher gain setting with a higher signal-to-noise ratio and some saturated spot intensities.

### Usage

mergeScansRG(RGlow, RGhigh, AboveNoiseLowG=NULL, AboveNoiseLowR=NULL, outlierp=0

### **Arguments**

RGlow object of class RGList containing red and green intensities constituting two-

color microarray data scanned at a lower gain setting.

RGhigh object of class RGList containing red and green intensities constituting two-

color microarray data scanned at a higher gain setting.

AboveNoiseLowG

matrix of 1 or 0 for low scan intensities of green color, 1 for spots above noise level or 0 otherwise. One column per array.

AboveNoiseLowR

matrix of 1 or 0 for low scan intensities of red color, 1 for spots above noise

level or 0 otherwise. One column per array.

outlierp p-value for outliers. 0 for no outlier detection or any value between 0 and 1.

Default p-value is 0.01.

78 mergeScans

#### **Details**

This function merges two separate scans of each fluorescent label on a two-color array scanned at two different scanner settings by using a nonlinear regression model consisting of two linear regression lines and a quadratic function connecting the two, which looks like a hockey stick. The changing point, i.e. the saturation point, in high scan is also estimated as part of model. Signals produced for certain spots can sometimes be very low (below noise) or too high (saturated) to be accurately read by the scanner. The proportions of spots that are below noise or above saturation are affected by the settings of the laser scanner used to read the arrays, with low scans minimizing saturation effects and high scans maximizing signal-to-noise ratios. Saturated spots can cause bias in intensity ratios that cannot be corrected for using conventional normalization methods.

Each fluorescent label on a two-color array can be scanned twice: for example, a high scan targeted at reaching saturation level for the brightest 1 percent of the spots on the array, and a low scan targeted at the lowest level of intensity which still allowed accurate grid placement on the arrays. By merging data from two separate laser scans of each fluorescent label on an array, we can avoid the potential bias in signal intensities due to below noise or above saturation and, thus provide better estimates of true differential expression as well as increase usable spots.

The merging process is designed to retain signal intensities from the high scan except when scanner saturation causes the high scan signal to be under-measured. The saturated spots are predicted from the corresponding low scans by the fitted regression model. It also checks any inconsistency between low and high scans.

#### Value

An object of class RGList-class with the following components:

G	numeric matrix containing the merged green (cy3) foreground intensities. Rows correspond to spots and columns to arrays.
R	numeric matrix containing the merged red (cy5) foreground intensities. Rows correspond to spots and columns to arrays.
Gb	numeric matrix containing the green (cy3) background intensities from high scan.
Rb	numeric matrix containing the red (cy5) background intensities from high scan.
other	list numeric matrices Gsaturated, Rsatured, Goutlier and Routlier. The first two contain saturation flags (1=saturated, 0=otherwise) for the green (cy3) and red (Cy5) channels of the high scan. The second two contain outlier flags (1=outlier, 0=otherwise) for the green (cy3) and red (Cy5) channels.

### Author(s)

Dongseok Choi <choid@ohsu.edu>.

### References

Choi D, O'Malley JP, Lasarev MR, Lapidus J, Lu X, Pattee P, Nagalla SR (2006). Extending the Dynamic Range of Signal Intensities in DNA Microarrays. *Online Journal of Bioinformatics*, **7**, 46-56.

#### **Examples**

```
## Not run:
#RG1: An RGList from low scan
```

modelMatrix 79

```
#RG2: An RGList from high scan
RGmerged <- mergeScansRG(RG1,RG2,AboveNoiseLowG=ANc3,AboveNoiseLowR=ANc5)

#merge two scans when all spots are above noise in low scan and no outlier detection.
RGmerged <- mergeScansRG(RG1,RG2,outlierp=0)

## End(Not run)</pre>
```

modelMatrix

Construct Design Matrix

#### **Description**

Construct design matrix from RNA target information for a two colour microarray experiment.

### Usage

```
modelMatrix(targets, parameters, ref, verbose=TRUE)
uniqueTargets(targets)
```

### **Arguments**

_	
targets	matrix or data.frame with columns $\texttt{Cy3}$ and $\texttt{Cy5}$ specifying which RNA was hybridized to each array
parameters	matrix specifying contrasts between RNA samples which should correspond to regression coefficients. Row names should correspond to unique RNA sample names found in targets.
ref	character string giving name of one of the RNA sources to be treated as reference. Exactly one argument of parameters or ref should be specified.
verbose	logical, if TRUE then unique names found in targets will be printed to standard output $\ensuremath{Largets}$

### **Details**

This function computes a design matrix for input to lmFit when analysing two-color microarray experiments in terms of log-ratios.

If the argument ref is used, then the experiment is treated as a one-way layout and the coefficients measure expression changes relative to the RNA source specified by ref. The RNA source ref is often a common reference which appears on every array or is a control sample to which all the others are compared. There is no restriction however. One can choose ref to be any of the RNA sources appearing the Cy3 or Cy5 columns of targets.

If the parameters argument is set, then the columns of this matrix specify the comparisons between the RNA sources which are of interest. This matrix must be of size n by (n-1), where n is the number of unique RNA sources found in Cy3 and Cy5, and must have row names which correspond to the RNA sources.

### Value

modelMatrix produces a numeric design matrix with row names as in targets and column names as in parameters.

unique Targets produces a character vector of unique target names from the columns Cy3 and Cy5 of targets.

80 modifyWeights

#### Author(s)

Gordon Smyth

#### See Also

```
model.matrix in the stats package.
```

An overview of linear model functions in limma is given by 06.LinearModels.

### **Examples**

```
targets <- cbind(Cy3=c("Ref","Control","Ref","Treatment"),Cy5=c("Control","Ref","Treatment
rownames(targets) <- paste("Array",1:4)

parameters <- cbind(C=c(-1,1,0),T=c(-1,0,1))
rownames(parameters) <- c("Ref","Control","Treatment")

modelMatrix(targets, parameters)
modelMatrix(targets, ref="Ref")</pre>
```

modifyWeights

modifyWeights

#### **Description**

Modify weights matrix for given gene status values.

# Usage

```
modifyWeights(weights=rep(1,length(status)), status, values, multipliers)
```

### **Arguments**

weights	numeric matrix of relative weights, rows corresponding to genes and columns to arrays
status	character vector giving the control status of each spot on the array, of same length as the number of rows of weights
values	character vector giving subset of the unique values of status
multipliers	numeric vector of same length as values giving factor by which weights will be modified

#### **Details**

The function is usually used to temporarily modify the weights matrix during normalization of data. The function can be used for example to give zero weight to spike-in ratio control spots during normalization.

### Value

Numeric matrix of same dimensions as weights with rows corresponding to values in status modified by the specified multipliers.

mrlm 81

### Author(s)

Gordon Smyth

#### See Also

An overview of normalization functions available in LIMMA is given in 05.Normalization.

### **Examples**

```
w <- matrix(runif(6*3),6,3)
status <- c("Gene", "Gene", "Ratio_Control", "Ratio_Control", "Gene", "Gene")
modifyWeights(w, status, values="Ratio_Control", multipliers=0)</pre>
```

mrlm

Fit Linear Model to Microrray Data by Robust Regression

### **Description**

Fit a linear model genewise to expression data from a series of arrays. The fit is by robust Mestimation allowing for a small proportion of outliers. This is a utility function for lmFit.

### Usage

```
mrlm(M, design=NULL, ndups=1, spacing=1, weights=NULL, ...)
```

### **Arguments**

М	numeric matrix containing log-ratio or log-expression values for a series of mi- croarrays, rows correspond to genes and columns to arrays.
design	numeric design matrix defining the linear model, with rows corresponding to arrays and columns to comparisons to be estimated. The number of rows must match the number of columns of M. Defaults to the unit vector meaning that the arrays are treated as replicates.
ndups	a positive integer giving the number of times each gene is printed on an array. nrow (M) must be divisible by ndups.
spacing	the spacing between the rows of M corresponding to duplicate spots, $\verb"spacing=1"$ for consecutive spots.
weights	numeric matrix of the same dimension as ${\tt M}$ containing weights. If it is of different dimension to ${\tt M}$ , it will be filled out to the same size. ${\tt NULL}$ is equivalent to equal weights.
• • •	any other arguments are passed to rlm.default.

#### **Details**

This is a utility function used by the higher level function lmFit. Most users should not use this function directly but should use lmFit instead.

This function fits a linear model for each gene by calling the function rlm from the MASS library.

Warning: don't use weights with this function unless you understand how rlm treats weights. The treatment of weights is somewhat different from that of lm.series and gls.series.

82 neqc

#### Value

### A list with components

coefficients numeric matrix containing the estimated coefficients for each linear model. Same number of rows as M, same number of columns as design.

stdev.unscaled

numeric matrix conformal with <code>coef</code> containing the unscaled standard deviations for the coefficient estimators. The standard errors are given by <code>stdev.unscaled</code>

\* sigma.

numeric vector containing the residual standard deviation for each gene.

df.residual numeric vector giving the degrees of freedom corresponding to sigma.

qr QR decomposition of design.

### Author(s)

Gordon Smyth

#### See Also

rlm.

An overview of linear model functions in limma is given by 06.LinearModels.

neqc NormExp and Quantile by Control (NEQC)

**Description** 

Perform normexp background correction and quantile normalization using control probes.

### Usage

```
neqc(x, status=NULL, negctrl="negative", regular="regular", offset=16, ...)
```

### **Arguments**

X	object of class EListRaw-class or matrix containing raw intensities for regular and control probes from a series of microarrays.
status	character vector giving probe types.
negctrl	character string giving the ${\tt status}$ identifier of negative control probes. Default is ${\tt negative}$ .
regular	character string giving the status identifier of regular probes. Default is regular.
offset	numeric value added to the intensities after background correction.
	any other parameters are passed to normalizeBetweenArrays.

#### **Details**

This function calls the function normexp.fit.control to estimate the parameters required by normal+exponential convolution model with the help of negative control probes. normexp.signal is then called to background correct the raw data. An offset is added to the data after the background correction. This function will then call the function normalizeBetweenArrays to perform quantile between-array normalization and log2 transformation.

For more descriptions to parameters x, status, negctrl and regular, please refer to functions normexp.fit.control and read.ilmn.

#### Value

An EList-class object containing normalized log2 expression values. Control probes are removed.

### Author(s)

Wei Shi and Gordon Smyth

#### References

Shi W, Oshlack A and Smyth GK. Calibrating the noise versus bias trade-off: normalization of Illumina Whole Genome Expression BeadChips. Submitted.

#### See Also

An overview of LIMMA functions for normalization is given in 05.Normalization.

An overview of background correction functions is given in 04.Background.

normexp.fit.control estimates the parameters in the normal+exponential convolution model using the negative control probes.

normexp.fit estimates parameters in the normal+exponential convolution model using a saddle-point approximation or other methods.

### **Examples**

```
## Not run:
x <- read.ilmn(files="sample probe profile.txt",ctrlfiles="control probe profile.txt")
y <- neqc(x)
## End(Not run)</pre>
```

normalizeBetweenArrays

Normalize Between Arrays

### **Description**

Normalizes expression intensities so that the intensities or log-ratios have similar distributions across a set of arrays.

#### **Usage**

```
normalizeBetweenArrays(object, method=NULL, targets=NULL, ...)
```

#### **Arguments**

object	a numeric matrix, EListRaw, RGList or MAList object.
method	character string specifying the normalization method to be used. Choices are "none", "scale", "quantile", "Aquantile", "Gquantile", "Rquantile", "Tquantile" or "vsn". A partial string sufficient to uniquely identify the choice is permitted. Default is "Aquantile" for two-color data objects or "quantile" for single-channel objects.
targets	vector, factor or matrix of length twice the number of arrays, used to indicate target groups if method="Tquantile"
• • •	other arguments are passed to normalizeQuantiles if one of the quantile methods are used or to vsn if method="vsn"

#### **Details**

normalizeWithinArrays normalizes expression values to make intensities consistent within each array. normalizeBetweenArrays normalizes expression values to achieve consistency between arrays. For two-color arrays, normalization between arrays usually occurs after normalization within arrays. An exception is method="vsn", see below. For single-channel arrays, within array normalization is not usually relevant.

The scale normalization method was proposed by Yang et al (2001, 2002) and is further explained by Smyth and Speed (2003). The idea is simply to scale the log-ratios to have the same median-abolute-deviation (MAD) across arrays. This idea has also been implemented by the manormScale function in the marray package. The implementation here is slightly different in that the MAD scale estimator is replaced with the median-absolute-value and the A-values are normalized as well as the M-values.

Quantile normalization was proposed by Bolstad et al (2003) for Affymetrix-style single-channel arrays and by Yang and Thorne (2003) for two-color cDNA arrays. method="quantile" ensures that the intensities have the same empirical distribution across arrays and across channels. method="Aquantile" ensures that the A-values (average intensities) have the same empirical distribution across arrays leaving the M-values (log-ratios) unchanged. These two methods are called "q" and "Aq" respectively in Yang and Thorne (2003).

method="Tquantile" performs quantile normalization separately for the groups indicated by targets. targets may be a target frame such as read by readTargets or can be a vector indicating green channel groups followed by red channel groups.

method="Gquantile" ensures that the green (first) channel has the same empirical distribution across arrays, leaving the M-values (log-ratios) unchanged. This method might be used when the green channel is a common reference throughout the experiment. In such a case the green channel represents the same target throughout, so it makes compelling sense to force the distribution of intensities to be same for the green channel on all the arrays, and to adjust to the red channel accordingly. method="Rquantile" ensures that the red (second) channel has the same empirical distribution across arrays, leaving the M-values (log-ratios) unchanged. Both Gquantile and Rquantile normalization have the implicit effect of changing the red and green log-intensities by equal amounts.

If object is a matrix then the scale, quantile or vsn normalization will be applied to the columns. Applying method="Aquantile" when object is a matrix will produce an error.

If object is an EListRaw object, then normalization will be applied to the matrix object\$E of expression values, which will then be log2-transformed.

method="vsn" uses the vsnMatrix function from the vsn package. For this option the input object should contain raw intensities, i.e., prior to background correction, log-transformation or any normalization.

See the limma User's Guide for more examples of use of this function.

#### Value

If object is a matrix then normalizeBetweenArrays produces a matrix of the same size. If object is an EListRaw object, then an EList object with expression values on the log2 scale is produced. For two-color data, normalizeBetweenArrays produces an MAList object with M and A-values on the log2 scale.

#### Author(s)

Gordon Smyth

#### References

Bolstad, B. M., Irizarry R. A., Astrand, M., and Speed, T. P. (2003), A comparison of normalization methods for high density oligonucleotide array data based on bias and variance. *Bioinformatics* **19**, 185-193.

Smyth, G. K., and Speed, T. P. (2003). Normalization of cDNA microarray data. *Methods* 31, 265-273.

Yang, Y. H., Dudoit, S., Luu, P., and Speed, T. P. (2001). Normalization for cDNA microarray data. In *Microarrays: Optical Technologies and Informatics*, M. L. Bittner, Y. Chen, A. N. Dorsel, and E. R. Dougherty (eds), Proceedings of SPIE, Volume 4266, pp. 141-152.

Yang, Y. H., Dudoit, S., Luu, P., Lin, D. M., Peng, V., Ngai, J., and Speed, T. P. (2002). Normalization for cDNA microarray data: a robust composite method addressing single and multiple slide systematic variation. *Nucleic Acids Research* **30**(4):e15.

Yang, Y. H., and Thorne, N. P. (2003). Normalization for two-color cDNA microarray data. In: D. R. Goldstein (ed.), *Science and Statistics: A Festschrift for Terry Speed*, IMS Lecture Notes - Monograph Series, Volume 40, pp. 403-418.

#### See Also

An overview of LIMMA functions for normalization is given in 05.Normalization.

See also maNormScale in the marray package, normalize in the affy package, and vsn and vsnMatrix in the vsn package.

# Examples

```
ngenes <- 100
narrays <- 4
x <- matrix(rnorm(ngenes*narrays),100,4)
y <- normalizeBetweenArrays(x)</pre>
```

normalizeMedianAbsValues

Normalize Columns of a Matrix to have the Median Absolute Value

### **Description**

Performs scale normalization of an M-value matrix or an A-value matrix across a series of arrays. Users do not normally need to call these functions directly - use normalizeBetweenArrays instead.

### Usage

```
normalizeMedianAbsValues(x)
```

#### Arguments

Х

numeric matrix

#### **Details**

If x is a matrix of log-ratios of expression (M-values) then normalizeMedianAbsValues is very similar to scaling to equalize the median absolute deviation (MAD) as in Yang et al (2001, 2002). Here the median-absolute value is used for preference to as to not re-center the M-values.

normalizeMedianAbsValues is also used to scale the A-values when scale-normalization is applied to an MAList object.

#### Value

A numeric matrix of the same size as that input which has been scaled so that each column as the same median-absolute value.

### Author(s)

Gordon Smyth

### See Also

An overview of LIMMA functions for normalization is given in 05.Normalization.

#### **Examples**

```
M <- cbind(Arrayl=rnorm(10),Array2=2*rnorm(10))
normalizeMedianAbsValues(M)</pre>
```

normalizeForPrintorder 87

normalizeForPrintorder

Print-Order Normalization

### **Description**

Normalize intensity values on one or more spotted microarrays to adjust for print-order effects.

### Usage

```
normalizeForPrintorder(object, layout, start="topleft", method = "loess", separate normalizeForPrintorder.rg(R, G, printorder, method = "loess", separate.channels plotPrintorder(object, layout, start="topleft", slide = 1, method = "loess", separate.channels
```

### Arguments

object	an RGList or list object containing components R and G which are matrices containing the red and green channel intensities for a series of arrays	
R	numeric vector containing red channel intensities for a single microarray	
G	numeric vector containing the green channel intensities for a single microarray	
layout	list specifying the printer layout, see PrintLayout-class	
start	character string specifying where printing starts in each pin group. Choices are "topleft" or "topright".	
printorder	numeric vector specifying order in which spots are printed. Can be computed from printorder (layout, start=start).	
slide	positive integer giving the column number of the array for which a plot is required	
method	character string, "loess" if a smooth loess curve should be fitted through the print-order trend or "plate" if plate effects are to be estimated	
separate.channels		
	logical, TRUE if normalization should be done separately for the red and green channel and FALSE if the normalization should be proportional for the two channels	
span	numerical constant between 0 and 1 giving the smoothing span for the loess the curve. Ignored if $method="plate"$ .	
plate.size	positive integer giving the number of consecutive spots corresponding to one plate or plate pack. Ignored if method="loess".	
plot	logical. If $\mathtt{TRUE}$ then a scatter plot of the print order effect is sent to the current graphics device.	

#### **Details**

Print-order is associated with the 384-well plates used in the printing of spotted microarrays. There may be variations in DNA concentration or quality between the different plates. The may be variations in ambient conditions during the time the array is printed.

This function is intended to pre-process the intensities before other normalization methods are applied to adjust for variations in DNA quality or concentration and other print-order effects.

88 normalizeForPrintorder

Printorder means the order in which spots are printed on a microarray. Spotted arrays are printed using a print head with an array of print-tips. Spots in the various tip-groups are printed in parallel. Printing is assumed to start in the top left hand corner of each tip-groups and to proceed right and down by rows, or else to start in the top right hand and to proceed left and down by rows. See printorder for more details. (WARNING: this is not always the case.) This is true for microarrays printed at the Australian Genome Research Facility but might not be true for arrays from other sources.

If object is an RGList then printorder is performed for each intensity in each array.

 $\verb|plotPrintorder| is a non-generic function which calls \verb|normalizeForPrintorder| with \\ \verb|plot=TRUE|.$ 

#### Value

normalizeForPrintorder produces an RGList containing normalized intensities.

The function plotPrintorder or normalizeForPrintorder.rg with plot=TRUE returns no value but produces a plot as a side-effect.

normalizeForPrintorder.rg with plot=FALSE returns a list with the following components:

R	numeric vector containing the normalized red channel intensities
G	numeric vector containing the normalized red channel intensites
R.trend	numeric vector containing the fitted printorder trend for the red channel
G.trend	numeric vector containing the fitted printorder trend for the green channe

#### Author(s)

Gordon Smyth

### References

```
Smyth, G. K. Print-order normalization of cDNA microarrays. March 2002. http://www.statsci.org/smyth/pubs/porder/porder.html
```

### See Also

```
printorder.
```

An overview of LIMMA functions for normalization is given in 05.Normalization.

### **Examples**

```
if(require("sma")) {
  data(MouseArray)
  plotPrintorder(mouse.data, mouse.setup, slide=1, separate=TRUE)
  RG <- normalizeForPrintorder(mouse.data, mouse.setup)
}</pre>
```

normalizeQuantiles 89

normalizeQuantiles Normalize Columns of a Matrix to have the same Quantiles

### Description

Normalize the columns of a matrix to have the same quantiles, allowing for missing values. Users do not normally need to call this function directly - use normalizeBetweenArrays instead.

### Usage

```
normalizeQuantiles(A, ties=TRUE)
```

#### **Arguments**

A numeric matrix. Missing values are allowed.

ties logical. If TRUE, ties in each column of A are treated in careful way. tied values

will be normalized to the mean of the corresponding pooled quantiles.

#### **Details**

This function is intended to normalize single channel or A-value microarray intensities between arrays. Each quantile of each column is set to the mean of that quantile across arrays. The intention is to make all the normalized columns have the same empirical distribution. This will be exactly true if there are no missing values and no ties within the columns: the normalized columns are then simply permutations of one another.

If there are ties amongst the intensities for a particular array, then with ties=FALSE the ties are broken in an unpredictable order. If ties=TRUE, all the tied values for that array will be normalized to the same value, the average of the quantiles for the tied values.

### Value

A matrix of the same dimensions as A containing the normalized values.

### Author(s)

Gordon Smyth

### References

Bolstad, B. M., Irizarry R. A., Astrand, M., and Speed, T. P. (2003), A comparison of normalization methods for high density oligonucleotide array data based on bias and variance. *Bioinformatics* **19**, 185-193.

#### See Also

An overview of LIMMA functions for normalization is given in 05.Normalization.

```
normalizeRobustSpline
```

Normalize Single Microarray Using Shrunk Robust Splines

### **Description**

Normalize the M-values for a single microarray using robustly fitted regression splines and empirical Bayes shrinkage.

### Usage

```
normalizeRobustSpline(M, A, layout, df=5, method="M")
```

### **Arguments**

M	numeric vector of M-values
A	numeric vector of A-values
layout	list specifying the dimensions of the spot matrix and the grid matrix
df	degrees of freedom for regression spline, i.e., the number of regression coefficients and the number of knots

choices are "M" for M-estimation or "MM" for high breakdown point regression

#### **Details**

method

This function implements an idea similar to print-tip loess normalization but uses regression splines in place of the loess curves and uses empirical Bayes ideas to shrink the individual prtin-tip curves towards a common value. This allows the technique to introduce less noise into good quality arrays with little spatial variation while still giving good results on arrays with strong spatial variation.

#### Value

Numeric vector containing normalized M-values.

### Author(s)

Gordon Smyth

#### References

The function is based on unpublished work by the author.

#### See Also

An overview of LIMMA functions for normalization is given in 05.Normalization.

### **Examples**

```
if(require("sma")) {
  data(MouseArray)
MA <- MA.RG(mouse.data)
normM <- normalizeRobustSpline(MA$M[,1],MA$A[,1],mouse.setup)
}</pre>
```

```
normalizeWithinArrays
```

Normalize Within Arrays

### **Description**

Normalize the expression log-ratios for one or more two-colour spotted microarray experiments so that the log-ratios average to zero within each array or sub-array.

### Usage

```
normalizeWithinArrays(object, layout, method="printtiploess", weights=object$wei
MA.RG(object, bc.method="subtract", offset=0)
RG.MA(object)
```

#### **Arguments**

object	object of class list, RGList or MAList containing red and green intensities constituting two-color microarray data.
layout	list specifying the dimensions of the spot matrix and the grid matrix. For details see PrintLayout-class.
method	character string specifying the normalization method. Choices are "none", "median", "loess", "printtiploess", "composite", "control" and "robustspline". A partial string sufficient to uniquely identify the choice is permitted.
weights	numeric matrix or vector of the same size and shape as the components of object containing spot quality weights.
span	numeric scalar giving the smoothing parameter for the loess fit
iterations	number of iterations used in loess fitting. More iterations give a more robust fit.
controlspots	numeric or logical vector specifying the subset of spots which are non-differentially-expressed control spots, for use with method="composite" or method="control".
df	degrees of freedom for spline if method="robustspline".
robust	robust regression method if method="robustspline". Choices are "M" or "MM".
bc.method	character string specifying background correct method, see backgroundCorrect for options.
offset	$numeric\ value, intensity\ offset\ used\ when\ computing\ log-ratios, see \verb backgroundCorrect .$

### Details

Normalization is intended to remove from the expression measures any systematic trends which arise from the microarray technology rather than from differences between the probes or between the target RNA samples hybridized to the arrays.

This function normalizes M-values (log-ratios) for dye-bias within each array. Apart from method="none" and method="median", all the normalization methods make use of the relationship between dye-bias and intensity. Method "none" computes M-values and A-values but does no normalization. Method "median" subtracts the weighted median from the M-values for each array.

The loess normalization methods ("loess", "printtiploess" and "composite") were proposed by Yang et al (2001, 2002). Smyth and Speed (2003) review these methods and describe how the methods are implemented in the limma package, including choices of tuning parameters. More information on the loess control parameters span and iterations can be found under loessFit. The default values used here are equivalent to those for the older function stat.ma in the sma package.

Oshlack et al (2004) consider the special issues that arise when a large proportion of probes are differentially expressed. They propose an improved version of composite loess normalization, which is implemented in the "control" method. This fits a global loess curve through a set of control spots, such as a whole-library titration series, and applies that curve to all the other spots.

The "robustspline" method calls normalizeRobustSpline. See that function for more documentation.

MA.RG converts an unlogged RGList object into an MAList object. MA.RG (object) is equivalent to normalizeWithinArrays (object, method="none").

RG.MA(object) converts back from an MAList object to a RGList object with unlogged intensities.

weights is normally a matrix giving a quality weight for every spot on every array. If weights is instead a vector or a matrix with only one column, then the weights will be assumed to be the same for every array, i.e., the weights will be probe-specific rather than spot-specific.

#### Value

An object of class MAList. Any components found in object will preserved except for R, G, Rb, Gb and other.

### Author(s)

Gordon Smyth

### References

Oshlack, A., Emslie, D., Corcoran, L., and Smyth, G. K. (2007). Normalization of boutique two-color microarrays with a high proportion of differentially expressed probes. *Genome Biology* **8**, R2.

Smyth, G. K., and Speed, T. P. (2003). Normalization of cDNA microarray data. *Methods* 31, 265-273.

Yang, Y. H., Dudoit, S., Luu, P., and Speed, T. P. (2001). Normalization for cDNA microarray data. In *Microarrays: Optical Technologies and Informatics*, M. L. Bittner, Y. Chen, A. N. Dorsel, and E. R. Dougherty (eds), Proceedings of SPIE, Vol. 4266, pp. 141-152.

Yang, Y. H., Dudoit, S., Luu, P., Lin, D. M., Peng, V., Ngai, J., and Speed, T. P. (2002). Normalization for cDNA microarray data: a robust composite method addressing single and multiple slide systematic variation. *Nucleic Acids Research* **30**(4):e15.

#### See Also

An overview of limma functions for normalization is given in 05.Normalization. In particular, see normalizeBetweenArrays for between-array normalization.

The original loess normalization function was the statma funtion in the sma package. normalizeWithinArrays is a direct generalization of that function, with more options and with support for quantitative spot quality weights.

normexp.fit.control 93

A different implementation of loess normalization methods, with potentially different behavior, is provided by the maNorm in the marray package.

### **Examples**

# See normalizeBetweenArrays

```
normexp.fit.control
```

Normexp Model Parameter Estimation Aided by Negative Controls

#### **Description**

The mean and log-standard-deviation of the background-normal part of the normexp+exponential convolution model is estimated as the mean and log-standard deviation of intensities from negative control probes. The log-mean of the signal-exponential part is estimated as the log of the difference between signal mean and background mean.

### Usage

```
normexp.fit.control(x, status=NULL, negctrl="negative")
```

### **Arguments**

X	object of class EListRaw-class or matrix containing raw intensities for
	regular and control probes for a series of microarrays
status	character vector giving probe types.

negctrl character string giving the type identifier of negative control probes.

#### **Details**

x has to contain raw expression intensities from both regular probes and negative control probes.

The probe type information for an object of <code>EListRaw-class</code> is normally saved in the <code>Status</code> column of its <code>genes</code> component. However, it will be overriden by the <code>status</code> parameter if it is explicitly provided to this function. If <code>x</code> is a <code>matrix</code> object, the probe type information has to be provided through the <code>status</code> parameter of this function. Regular probes have the status <code>regular</code>. Negative control probes have the status indicated by <code>negctrl</code>, which is <code>negative</code> by default.

This function estimates parameters of the normal+exponential convolution model with the help of negative control probes. The mean and log-standard-deviation of the background-normal part of the normexp+exponential convolution model are estimated as the mean and log-standard deviation of intensities from negative control probes respectively. The log-mean of the signal-exponential part is estimated as the log of the difference between signal mean and background mean. The signal mean is simply the mean of intensities from regular probes.

### Value

A matrix containing estimated parameters with rows being arrays and with columns being parameters. Column names are mu, logsigma and logalpha.

94 normexp.fit

#### Author(s)

Wei Shi and Gordon Smyth

#### References

Wei Shi and Gordon K Smyth. Normalizing Illumina Whole Genome Expression BeadChips. In preparation.

#### See Also

neqc calls this function to get the parameters of the normal+exponential convolution model and then calls normexp.signal to perform the background correction.

normexp.fit estimates parameters in the normal+exponential convolution model using a saddle-point approximation or other mothods.

An overview of background correction functions is given in 04.Background.

### **Examples**

```
## Not run:
x <- read.profile(files="sample probe profile", ctrlfiles="control probe profile")
par <- normexp.fit.control(x)
## End(Not run)</pre>
```

normexp.fit

Fit Normal+Exp Convolution Model to Observed Intensities

### **Description**

Fit the normal+exponential convolution model to a vector of observed intensities. The normal part represents the background and the exponential part represents the signal intensities. This function is called by backgroundCorrect and is not normally called directly by users.

#### Usage

```
normexp.fit(x, method="saddle", n.pts=NULL, trace=FALSE)
```

# Arguments

X	numeric vector of (background corrected) intensities
method	method used to estimate the three parameters. Choices for normexp.fit are "mle", "saddle", "rma" and "rma75".
n.pts	number of quantiles of $x$ to use for the fit. If $\mathtt{NULL}$ then all values of $x$ will be used.
trace	logical, if TRUE, tracing information on the progress of the optimization is given.

normexp.fit 95

#### **Details**

The Normal+Exp (normexp) convolution model is a mathematical model representing microarray intensity data for the purposes of background correction. It was proposed originally as part of the RMA algorithm for Affymetrix microarray data. For two-color microarry data, the normexp background correction method was introduced and compared with other methods by Ritchie et al (2007).

This function uses maximum likelihood estimation to fit the normexp model to background-corrected intensities. The model assumes that the observed intensities are the sum of background and signal components, the background being normal and the signal being exponential distributed.

The likelihood may be computed exactly (method="mle") or approximated using a saddle-point approximation (method="saddle"). The saddle-point approximation was proposed by Ritchie et al (2007). Silver et al (2008) added some computational refinements to the saddle-point approximation, making it more reliable in practice, and developed the exact likelihood maximization algorithm. The "mle" method uses the best performing algorithm from Silver et al (2008), which calls the optimization function nlminb with analytic first and second derivatives. Derivatives are computed with respect to the normal-mean, the log-normal-variance and the log-exponential-mean.

Two ad-hoc estimators are also available which do not require iterative estimation. "rma" results in a call to the bg.parameters function of the affy package. This provides the kernel estimation method that is part of the RMA algorithm for Affymetrix data. "rma75" uses the similar but less biased RMA-75 method from McGee and Chen (2006).

If the length x is very large, it may be worth saving computation time by setting n.pts to a value less than the total number of probes, for example n.pts= $2^14$ .

#### Value

A list containing the components

par numeric vector giving estimated values of the mean and log-standard-deviation

of the background-normal part and the log-mean of the signal-exponential part.

m2loglik numeric scalar giving minus twice the maximized log-likelihood

convergence integer code indicating successful convergence or otherwise of the optimization.

#### Author(s)

Gordon Smyth and Jeremy Silver

#### References

McGee, M., and Chen, Z. (2006). Parameter estimation for the exponential-normal convolution model for background correction of Affymetrix GeneChip data. *Stat Appl Genet Mol Biol*, 5(1), Article 24.

Ritchie, M. E., Silver, J., Oshlack, A., Silver, J., Holmes, M., Diyagama, D., Holloway, A., and Smyth, G. K. (2007). A comparison of background correction methods for two-colour microarrays. *Bioinformatics* http://bioinformatics.oxfordjournals.org/cgi/content/abstract/btm412

Silver, JD, Ritchie, ME, and Smyth, GK (2009). Microarray background correction: maximum likelihood estimation for the normal-exponential convolution. *Biostatistics* 10, 352-363. http://biostatistics.oxfordjournals.org/cgi/content/abstract/kxn042

96 normexp.signal

#### See Also

```
normexp.signal. Also bg.parameters in the affy package.
```

An overview of background correction functions is given in 04.Background.

### **Examples**

```
x \leftarrow c(2,3,1,10,3,20,5,6)
out \leftarrow normexp.fit(x)
normexp.signal(out$par, x=x)
```

normexp.signal

Expected Signal Given Observed Foreground Under Normal+Exp Model

### Description

Adjust foreground intensities for observed background using Normal+Exp Model. This function is called by backgroundCorrect and is not normally called directly by the user.

### Usage

```
normexp.signal(par, x)
```

#### **Arguments**

#### **Details**

In general the vector normmean is computed conditional on background at each spot.

#### Value

Numeric vector containing adjusted intensities.

### Author(s)

Gordon Smyth

#### References

Ritchie, M. E., Silver, J., Oshlack, A., Silver, J., Holmes, M., Diyagama, D., Holloway, A., and Smyth, G. K. (2007). A comparison of background correction methods for two-colour microarrays. Bioinformatics http://bioinformatics.oxfordjournals.org/cgi/content/abstract/btm412

Silver, JD, Ritchie, ME, and Smyth, GK (2009). Microarray background correction: maximum likelihood estimation for the normal-exponential convolution. *Biostatistics* 10, 352-363. http://biostatistics.oxfordjournals.org/cgi/content/abstract/kxn042

plotDensities 97

#### See Also

```
normexp.fit
```

An overview of background correction functions is given in 04.Background.

### **Examples**

```
# See normexp.fit
```

plotDensities

Individual-channel Densities Plot

### **Description**

Plots the densities of individual-channel intensities for two-color microarray data.

### Usage

```
plotDensities(object, log=TRUE, arrays=NULL, singlechannels=NULL, groups=NULL, c
```

an RGList or MAList object. RGList objects containing logged or unlogged

### **Arguments**

object

_	intensities can be accommodated using the log.transform argument.
log	logical, should densities be formed and plotted for the log-intensities (TRUE) or raw intensities (FALSE)?
arrays	vector of integers giving the arrays from which the individual-channels will be selected to be plotted. Corresponds to columns of $M$ and $A$ (or $R$ and $G$ ). Defaults to all arrays.
singlechanne	ls
	vector of integers indicating which individual-channels will be selected to be plotted. Values correspond to the columns of the matrix of $cbind(R,G)$ and range between 1: $ncol(R)$ for red channels and $(ncol(R)+1):(ncol(R)+ncol(G))$ ) for the green channels in object. Defaults to all channels.
groups	vector of consecutive integers beginning at 1 indicating the groups of arrays or individual-channels (depending on which of arrays or singlechannels

col

vector of colors of the same length as the number of different groups. If  $\mathtt{NULL}$  (default) the  $\mathtt{col}$  equals  $\mathtt{c("red", "green")}$ . See details for more specifi-

groups (if specified) must correspond to the arrays.

are non NULL). This is used to color any groups of the individual-channel densities. If NULL (default), groups correspond to the red and green channels. If both arrays and singlechannels are NULL all arrays are selected and

cations.

98 plotDensities

#### **Details**

This function is used as a data display technique associated with between-array normalization, especially individual-channel normalization methods such as quantile-normalization. See the section on between-array normalization in the LIMMA User's Guide.

If no col is specified, the default is to color individual channels according to red and green. If both arrays and groups are non-NULL, then the length of groups must equal the length of arrays and the maximum of groups (i.e. the number of groups) must equal the length of col otherwise the default color of black will be used for all individual-channels. If arrays is NULL and both singlechannels and groups are non-NULL, then the length of groups must equal the length of singlechannels and the maximum of groups (i.e. the number of groups) must equal the length of col otherwise the default color of black will be used for all individual-channels.

#### Value

A plot is created on the current graphics device.

#### Author(s)

Natalie Thorne

#### See Also

An overview of diagnostic plots in LIMMA is given in 09.Diagnostics. There is a section using plotDensities in conjunction with between-array normalization in the LIMMA User's Guide.

### **Examples**

```
if(require("sma")) {
data (MouseArray)
# no normalization but background correction is done
MA.n <- MA.RG(mouse.data)</pre>
# Default settings for plotDensities.
plotDensities(MA.n)
# One can reproduce the default settings.
plotDensities (MA.n, arrays=c(1:6), groups=c(rep(1,6), rep(2,6)),
col=c("red", "green"))
# Color R and G individual-channels by blue and purple.
plotDensities(MA.n,arrays=NULL,groups=NULL,col=c("blue","purple"))
# Indexing individual-channels using singlechannels (arrays=NULL).
plotDensities(MA.n, singlechannels=c(1,2,7))
# Change the default colors from c("red", "green") to c("pink", "purple")
plotDensities(MA.n, singlechannels=c(1,2,7), col=c("pink", "purple"))
# Specified too many colors since groups=NULL defaults to two groups.
plotDensities(MA.n,singlechannels=c(1,2,7),col=c("pink","purple","blue"))
# Three individual-channels, three groups, three colors.
plotDensities (MA.n, singlechannels=c(1,2,7), groups=c(1,2,3),
col=c("pink", "purple", "blue"))
```

plotFB 99

```
# Three individual-channels, one group, one color.
plotDensities(MA.n,singlechannels=c(1,2,7),groups=c(1,1,1),
col=c("purple"))

# All individual-channels, three groups (ctl,tmt,reference), three colors.
plotDensities(MA.n,singlechannels=c(1:12),
groups=c(rep(1,3),rep(2,3),rep(3,6)),col=c("darkred","red","green"))
}
```

plotFB

FB-Plot

### **Description**

Creates foreground-background plots.

### Usage

```
plotFB(RG, array=1, lim="separate", pch=16, cex=0.2, ...)
```

### **Arguments**

RG	an RGList object.
array	integer giving the array to be plotted. Corresponds to columns of ${\tt R},{\tt G},{\tt Rb}$ and ${\tt Gb}.$
lim	character string indicating whether the red and green plots should have "separate" or "common" $x$ - and $y$ - co-ordinate limits.
pch	vector or list of plotting characters. Defaults to integer code 16.
cex	numeric vector of plot symbol expansions.
• • •	any other arguments are passed to plot

### Details

A foreground-background plot is a plot of log2-foreground vs log2-background for a particular channel on a particular two-color array. This function produces a pair of plots, one for green and one for red, for a specified array.

See points for possible values for pch, col and cex.

#### Value

A plot is created on the current graphics device.

### Author(s)

Gordon Smyth

### See Also

An overview of diagnostic functions available in LIMMA is given in 09.Diagnostics.

100 plotlines

# Description

Time course style plot of expression data.

### Usage

```
plotlines(x, first.column.origin=FALSE, xlab="Column", ylab="x", col="black", lwd=1,.
```

### **Arguments**

```
numeric matrix or object containing expression data.

first.column.origin
logical, should the lines be started from zero?

xlab x-axis label
ylab y-axis label
col vector of colors for lines

lwd line width multiplier
... any other arguments are passed to plot
```

# Details

Plots a line for each probe.

### Value

A plot is created on the current graphics device.

# Author(s)

Gordon Smyth

# See Also

An overview of modeling functions and associated plots available in LIMMA is given in 06.Linear-Models.

plotMA3by2

plotMA3by2	Write MA-Plots to Files	

### **Description**

Write MA-plots to files in PNG format, six plots to a file in a 3 by 2 grid arrangement.

### Usage

```
plotMA3by2(MA, prefix="MA", path=NULL, main=colnames(MA), zero.weights=FALSE, collaboration
```

### **Arguments**

MA	an MAList or RGList object, or any list with components M containing log- ratios and A containing average intensities
prefix	character string giving prefix to attach to file names
path	character string specifying directory for output files
main	character vector giving titles for plots
zero.weights	logical, should points with non-positive weights be plotted
common.lim	logical, should all plots on a page use the same axis limits
device	$device \ driver \ for \ the \ plot. \ Choices \ are \ "png", "jpeg", "pdf", "postscript".$
	any other arguments are passed to plotMA

### **Details**

This function writes a series of graphic files to disk. Each file contains six MA-plots in three rows and two columns. The layout is optimized for A4-sized paper.

The graph format can be "png" or "jpeg", which are screen-resolution formats, or "pdf" or "postscript", which are loss-less formats. "png" is not available on every R platform. Note that "pdf" or "postscript" may produce very large files.

### Value

No value is returned, but one or more files are written to the working directory. The number of files is determined by the number of columns of MA.

### Author(s)

Gordon Smyth

### See Also

An overview of diagnostic functions available in LIMMA is given in 09. Diagnostics.

102 plotMA

|--|

# Description

Creates an MA-plot with color coding for control spots.

# Usage

```
plotMA(MA, array=1, xlab="A", ylab="M", main=colnames(MA)[array], xlim=NULL, yli
```

# Arguments

MA	an RGList, MAList or MArrayLM object, or any list with components M containing log-ratios and A containing average intensities. Alternatively a matrix or ExpressionSet object.
array	integer giving the array to be plotted. Corresponds to columns of $\ensuremath{\mathtt{M}}$ and $\ensuremath{\mathtt{A}}.$
xlab	character string giving label for x-axis
ylab	character string giving label for y-axis
main	character string giving title for plot
xlim	numeric vector of length 2 giving limits for x-axis, defaults to min and max of the data
ylim	numeric vector of length 2 giving limits for y-axis, defaults to min and max of the data
status	character vector giving the control status of each spot on the array, of same length as the number of rows of MA\$M. If omitted, all points are plotted in the default color, symbol and size.
values	character vector giving values of status to be highlighted on the plot. Defaults to unique values of status. Ignored if there is no status vector.
pch	vector or list of plotting characters. Default is integer code 16 which gives a solid circle. Ignored is there is no status vector.
col	numeric or character vector of colors, of the same length as values. Defaults to 1:length(values). Ignored if there is no status vector.
cex	numeric vector of plot symbol expansions, of the the same length as values. Defaults to 0.2 for the most common status value and 1 for the others. Ignored if there is no status vector.
legend	logical, should a legend of plotting symbols and colors be included. Ignored if there is no status vector.
zero.weights	logical, should spots with zero or negative weights be plotted?
• • •	any other arguments are passed to plot

plotMA 103

#### **Details**

An MA-plot is a plot of log-intensity ratios (M-values) versus log-intensity averages (A-values). If MA is an RGList or MAList then this function produces an ordinary within-array MA-plot. If MA is an MArrayLM object, then the plot is an fitted model MA-plot in which the estimated coefficient is on the y-axis and the average A-value is on the x-axis.

If MA is a matrix or ExpressionSet object, then this function produces a between-array MA-plot. In this case the A-values in the plot are the average log-intensities across the arrays and the M-values are the deviations of the log-intensities for the specified array from the average. If there are more than five arrays, then the average is computed robustly using medians. With five or fewer arrays, it is computed by means.

The status vector is intended to specify the control status of each spot, for example "gene", "ratio control", "house keeping gene", "buffer" and so on. The vector is usually computed using the function controlStatus and a spot-types file. However the function may be used to highlight any subset of spots.

The status can be included as the component MA\$genes\$Status instead of being passed as an argument to plotMA. The arguments values, pch, col and cex can be included as attributes to status instead of being passed as arguments to plotMA.

See points for possible values for pch, col and cex.

#### Value

A plot is created on the current graphics device.

#### Author(s)

Gordon Smyth

#### References

```
See http://www.statsci.org/micrarra/refs/maplots.html
```

#### See Also

An overview of diagnostic functions available in LIMMA is given in 09. Diagnostics.

#### **Examples**

```
MA <- new("MAList")
MA$A <- runif(300,4,16)
MA$M <- rt(300,df=3)
status <- rep("Gene",300)
status[1:3] <- "M=0"
MA$M[1:3] <- 0
status[4:6] <- "M=3"
MA$M[4:6] <- 3
status[7:9] <- "M=-3"
MA$M[7:9] <- -3
plotMA(MA,main="MA-Plot with Simulated Data",status=status,values=c("M=0","M=3","M=-3"),c

# Same as above
attr(status,"values") <- c("M=0","M=3","M=-3")
attr(status,"col") <- c("blue","red","green")
plotMA(MA,main="MA-Plot with Simulated Data",status=status)
```

104 plotMDS

```
# Same as above
MA$genes$Status <- status
plotMA(MA,main="MA-Plot with Simulated Data")</pre>
```

plotMDS

Multidimensional scaling plot of microarray data

#### **Description**

Plot the sample relations based on MDS.

### Usage

```
plotMDS(x, top=500, labels=colnames(x), col=NULL, cex=1, dim.plot=c(1,2), ndim=multiples(x), cex=1, dim.plot=c(1,2), dim.plot=c(
```

### **Arguments**

X	any data object which can be coerced to a matrix, such as ${\tt ExpressionSet}$ or ${\tt EList}.$
top	number of top genes used to calculate pairwise distances.
labels	character vector of sample names or labels. If $\mathbf{x}$ has no column names, then defaults the index of the samples.
col	numeric or character vector of colors for the plotting characters.
cex	numeric vector of plot symbol expansions.
dim.plot	which two dimensions should be plotted, numeric vector of length two.
ndim	number of dimensions in which data is to be represented
gene.selection	
	character, "pairwise" to choose the top genes separately for each pairwise comparison between the samples or "common" to select the same genes for all comparisons

### **Details**

. . .

This function is a variation on the usual multdimensional scaling (or principle coordinate) plot, in that a distance measure particularly appropriate for the microarray context is used. The distance between each pair of samples (columns) is the root-mean-square deviation (Euclidean distance) for the top top genes.

any other arguments are passed to plot.

If gene.selection is "common", then the top genes are those with the largest standard deviations between samples. If gene.selection is "pairwise", then a different set of top genes is selected for each pair of samples. The pairwise feature selection may be appropriate for microarray data when different molecular pathways are relevant for distinguishing different pairs of samples.

See text for possible values for col and cex.

#### Value

A plot is created on the current graphics device.

plotPrintTipLoess 105

#### Author(s)

Di Wu and Gordon Smyth

#### See Also

```
cmdscale
```

An overview of diagnostic functions available in LIMMA is given in 09. Diagnostics.

### **Examples**

```
# Simulate gene expression data for 1000 probes and 6 microarrays.
# Samples are in two groups
# First 50 probes are differentially expressed in second group
sd <- 0.3*sqrt(4/rchisq(1000,df=4))
x <- matrix(rnorm(1000*6,sd=sd),1000,6)
rownames(x) <- paste("Gene",1:1000)
x[1:50,4:6] <- x[1:50,4:6] + 2
# without labels, indexes of samples are plotted.
plotMDS(x, col=c(rep("black",3), rep("red",3)))
# with labels as groups, group indicators are plotted.
plotMDS(x, col=c(rep("black",3), rep("red",3)), labels= c(rep("Grp1",3), rep("Grp2",3)))</pre>
```

```
plotPrintTipLoess MA Plots by Print-Tip Group
```

#### **Description**

Creates a coplot giving MA-plots with loess curves by print-tip groups.

### Usage

```
plotPrintTipLoess(object,layout,array=1,span=0.4,...)
```

### **Arguments**

object	MAList or RGList object or list with components M containing log-ratios and A containing average intensities
layout	a list specifying the number of tip rows and columns and the number of spot rows and columns printed by each tip. Defaults to MA\$printer if that is non-null.
array	integer giving the array to be plotted. Corresponds to columns of M and A.
span	span of window for lowess curve
	other arguments passed to panel.smooth

### **Details**

Note that spot quality weights in object are not used for computing the loess curves for this plot even though such weights would be used for loess normalization using normalizeWithinArrays.

### Value

A plot is created on the current graphics device. If there are missing values in the data, then the vector of row numbers for spots with missing values is invisibly returned, as for coplot.

106 plotRLDF

### Author(s)

Gordon Smyth

#### See Also

An overview of diagnostic functions available in LIMMA is given in 09. Diagnostics.

plotRLDF

Plot of regularized linear discriminant functions for microarray data

### **Description**

Plot of regularized linear discriminant functions for microarray data.

### Usage

```
plotRLDF(y, design=NULL, z=NULL, labels.y=NULL, labels.z=NULL, col.y=1, col.z=1,
df.prior=5, show.dimensions=c(1,2), main=NULL, nprobes=500,...)
```

### **Arguments**

У	any data object which can be coerced to a matrix, such as <code>ExpressionSet</code> or <code>EList</code> . The training dataset.
Z	any data object which can be coerced to a matrix, such as <code>ExpressionSet</code> or <code>EList</code> . The dataset to be classified.
design	the design matrix of the microarray experiment for y, with rows corresponding to arrays and columns to coefficients to be estimated. Defaults to the unit vector meaning that the arrays are treated as replicates.
labels.y	character vector of sample names or labels in ${\bf y}$ . Default is integers starting from 1.
labels.z	character vector of sample names or labels in z. Default is letters.
col.y	numeric or character vector of colors for the plotting characters of y. Default is black.
col.z	numeric or character vector of colors for the plotting characters of z. Default is black.
df.prior	prior degrees of freedom for residual variances. Used in gene selection.
show.dimensi	ons
	which two dimensions should be plotted, numeric vector of length two.
main	title of the plot.
nprobes	number of probes to be used for the calculations. Selected by moderated F tests.
	any other arguments are passed to plot.

### **Details**

This function is a variation on the plot of usual linear discriminant fuction, in that the within-group covariance matrix is regularized to ensure that it is invertible, with eigenvalues bounded away from zero. A diagonal regulation using df.prior and the median within-group variance is used.

The calculations are based on a filtered list of probes. The nprobes probes with largest moderated F statistics are used to discriminate.

See text for possible values for col and cex.

plotSA 107

#### Value

A list containing metagene information is (invisibly) returned. A plot is created on the current graphics device.

#### Author(s)

Di Wu and Gordon Smyth

#### See Also

lda in package MASS

### **Examples**

```
# Simulate gene expression data for 1000 probes and 6 microarrays.
# Samples are in two groups
# First 50 probes are differentially expressed in second group
sd <- 0.3*sqrt(4/rchisq(1000,df=4))
y <- matrix(rnorm(1000*6,sd=sd),1000,6)
rownames(y) <- paste("Gene",1:1000)
y[1:50,4:6] <- y[1:50,4:6] + 2

z <- matrix(rnorm(1000*6,sd=sd),1000,6)
rownames(z) <- paste("Gene",1:1000)
z[1:50,4:6] <- z[1:50,4:6] + 1.8
z[1:50,1:3] <- z[1:50,1:3] - 0.2

design <- cbind(Grp1=1,Grp2vs1=c(0,0,0,1,1,1))
options(digit=3)

plotRLDF(y,z, design=design)</pre>
```

plotSA

Sigma vs A plot for microarray linear model

### **Description**

Plot log residual standard deviation versus average log expression for a fitted microarray linear model.

#### Usage

```
plotSA(fit, xlab="Average log-expression", ylab="log2(sigma)", zero.weights=FALS
```

# Arguments

fit	an MArrayLM object.
xlab	character string giving label for x-axis
ylab	character string giving label for y-axis
pch	vector or list of plotting characters. Default is integer code 16 which gives a solid circle.

108 poolVar

```
cex numeric expansion factor for plotting character. Defaults to 0.2.
zero.weights logical, should spots with zero or negative weights be plotted?
... any other arguments are passed to plot
```

#### **Details**

This plot is used to check the mean-variance relationship of the expression data, after fitting a linear model.

See points for possible values for pch and cex.

#### Value

A plot is created on the current graphics device.

#### Author(s)

Gordon Smyth

#### See Also

An overview of diagnostic functions available in LIMMA is given in 09. Diagnostics.

poolVar	Pool Sample Variances with Unequal Variances

#### **Description**

Compute the Satterthwaite (1946) approximation to the distribution of a weighted sum of sample variances.

# Usage

```
poolVar(var, df=n-1, multiplier=1/n, n)
```

#### **Arguments**

var numeric vector of independent sample variances

df numeric vector of degrees of freedom for the sample variances

multiplier numeric vector giving multipliers for the sample variances

n numeric vector of sample sizes

#### **Details**

The sample variances var are assumed to follow scaled chi-square distributions. A scaled chi-square approximation is found for the distribution of sum(multiplier \* var) by equating first and second moments. On output the sum to be approximated is equal to multiplier \* var which follows approximately a scaled chisquare distribution on df degrees of freedom. The approximation was proposed by Satterthwaite (1946).

If there are only two groups and the degrees of freedom are one less than the sample sizes then this gives the denominator of Welch's t-test for unequal variances.

printHead 109

## Value

A list with components

var effective pooled sample variance
df effective pooled degrees of freedom

multiplier pooled multiplier

## Author(s)

Gordon Smyth

#### References

Welch, B. L. (1938). The significance of the difference between two means when the population variances are unequal. *Biometrika* **29**, 350-362.

Satterthwaite, F. E. (1946). An approximate distribution of estimates of variance components. *Biometrics Bulletin* **2**, 110-114.

Welch, B. L. (1947). The generalization of 'Student's' problem when several different population variances are involved. *Biometrika* **34**, 28-35.

Welch, B. L. (1949). Further note on Mrs. Aspin's tables and on certain approximations to the tabled function. *Biometrika* **36**, 293-296.

## See Also

10.Other

## **Examples**

```
# Welch's t-test with unequal variances
x <- rnorm(10, mean=1, sd=2)
y <- rnorm(20, mean=2, sd=1)
s2 <- c(var(x), var(y))
n <- c(10, 20)
out <- poolVar(var=s2, n=n)
tstat <- (mean(x)-mean(y)) / sqrt(out$var*out$multiplier)
pvalue <- 2*pt(-abs(tstat), df=out$df)
# Equivalent to t.test(x,y)</pre>
```

printHead

Print Leading Rows of Large Objects

# Description

Print the leading rows of a large vector, matrix or data.frame. This function is used by show methods for data classes defined in LIMMA.

# Usage

```
printHead(x)
```

110 PrintLayout

## **Arguments**

x any object

#### **Details**

If x is a vector with more than 20 elements, then printHead(x) prints only the first 5 elements. If x is a matrix or data.frame with more than 10 rows, then printHead(x) prints only the first 5 rows. Any other type of object is printed normally.

## Author(s)

Gordon Smyth

### See Also

An overview of classes defined in LIMMA is given in 02. Classes

PrintLayout

Print Layout - class

## **Description**

A list-based class for storing information about the process used to print spots on a microarray.

PrintLayout objects can be created using getLayout. The printer component of an RGList or MAList object is of this class.

# **Slots/List Components**

Objects of this class contains no slots but should contain the following list components:

ngrid.r: number of grid rows on the arrays
ngrid.c: number of grid columns on the arrays
nspot.r: number of rows of spots in each grid
number of columns of spots in each grid

ndups: number of duplicates of each DNA clone, i.e., number of times print-head dips into each well of DNA

spacing: number of spots between duplicate spots. Only applicable if ndups>1. spacing=1 for side-by-side spots actual number of pins or tips on the print-head

start: character string giving position of the spot printed first in each grid. Choices are "topleft" or "toprig

## Author(s)

Gordon Smyth

### See Also

02. Classes gives an overview of all the classes defined by this package.

# **Examples**

# Settings for Swirl and ApoAI example data sets in User's Guide

printorder 111

```
printer <- list(ngrid.r=4, ngrid.c=4, nspot.r=22, nspot.c=24, ndups=1, spacing=1, npins=1
# Typical settings at the Australian Genome Research Facility
# Full pin set, duplicates side-by-side on same row
printer <- list(ngrid.r=12, ngrid.c=4, nspot.r=20, nspot.c=20, ndups=2, spacing=1, npins=
# Half pin set, duplicates in top and lower half of slide
printer <- list(ngrid.r=12, ngrid.c=4, nspot.r=20, nspot.c=20, ndups=2, spacing=9600, npi</pre>
```

printorder

Identify Order in which Spots were Printed

## **Description**

Identify order in which spots were printed and the 384-well plate from which they were printed.

## Usage

```
printorder(layout, ndups=1, spacing="columns", npins, start="topleft")
```

## **Arguments**

layout	list with the components ngrid.r, ngrid.c, nspot.r and nspot.c, or an RGList or MAList object from which the printer layout may be extracted.	
ndups	number of duplicate spots, i.e., number of times print-head dips into each well	
spacing	character string indicating layout of duplicate spots. Choices are "columns", "rows" or "topbottom".	
npins	actual number of pins or tips on the print-head	
start	character string giving position of the spot printed first in each grid. Choices are "topleft" or "topright" and partial matches are accepted.	

### **Details**

In most cases the printer-head contains the layout \$ngrid.r times layout \$ngrid.c pins or tips and the array is printed using layout \$nspot.r times layout \$npot.c dips of the head. The plate holding the DNA to be printed is assumed to have 384 wells in 16 rows and 24 columns.

ndups indicates the number of spots printed from each well. The replicate spots from multiple dips into the same wells are assumed to be side-by-side by columns (spacing="columns"), by rows (spacing="rows") or in the top and bottom halves of the array (spacing="topbottom").

In some cases a smaller number of physical pins is used and the total number of grids is built up by effectively printing two or more sub-arrays on the same slide. In this case the number of grids should be a multiple of the number of pins.

Printing is assumed to proceed by rows within in each grid starting either from the top-left or the top-right.

112 printtipWeights

### Value

# List with components

printorder numeric vector giving printorder of each spot, i.e., which dip of the print-head was used to print it

plate numeric vector giving plate number from which each spot was printed

plate.r numeric vector giving plate-row number of the well from which each spot was printed

plate.c numeric vector giving plate-column number of the well from which each spot was printed

plateposition

character vector summarizing plate number and plate position of the well from which each spot was printed with letters for plate rows and number for columns.

For example 02B13 is second row, 13th column, of the second plate.

## Author(s)

Gordon Smyth

### See Also

normalizeForPrintorder.

An overview of LIMMA functions for reading data is given in 03.ReadingData.

## **Examples**

```
printorder(list(ngrid.r=2,ngrid.c=2,nspot.r=12,nspot.c=8))
```

printtipWeights Sub-array Quality Weights

# **Description**

Estimates relative quality weights for each sub-array in a multi-array experiment.

# Usage

```
printtipWeights(object, design = NULL, weights = NULL, method = "genebygene", la
```

# Arguments

object	object of class numeric, matrix, MAList, marrayNorm, or ExpressionSet containing log-ratios or log-values of expression for a series of spotted microarrays.
design	the design matrix of the microarray experiment, with rows corresponding to arrays and columns to coefficients to be estimated. Defaults to the unit vector meaning that the arrays are treated as replicates.
weights	optional numeric matrix containing prior weights for each spot.

printtipWeights 113

method	character string specifying the estimating algorithm to be used. Choices are "genebygene" and "reml".
layout	list specifying the dimensions of the spot matrix and the grid matrix. For details see PrintLayout-class.
maxiter	maximum number of iterations allowed.
tol	convergence tolerance.
trace	logical variable. If true then output diagnostic information at each iteration of "reml" algorithm.

### **Details**

The relative reliability of each sub-array (print-tip group) is estimated by measuring how well the expression values for that sub-array follow the linear model.

The method described in Ritchie et al (2006) and implemented in the <code>arrayWeights</code> function is adapted for this purpose. A heteroscedastic model is fitted to the expression values for each gene by calling the function <code>lm.wfit</code>. The dispersion model is fitted to the squared residuals from the mean fit, and is set up to have sub-array specific coefficients, which are updated in either full REML scoring iterations, or using an efficient gene-by-gene update algorithm. The final estimates of the sub-array variances are converted to weights.

The data object object is interpreted as for lmFit. In particular, the arguments design, weights and layout will be extracted from the data object if available and do not normally need to be set explicitly in the call; if any of these are set in the call then they will over-ride the slots or components in the data object.

### Value

A matrix of sub-array weights which can be passed to lmFit.

### Author(s)

Matthew Ritchie and Gordon Smyth

## References

Ritchie, M. E., Diyagama, D., Neilson, van Laar, R., J., Dobrovic, A., Holloway, A., and Smyth, G. K. (2006). Empirical array quality weights in the analysis of microarray data. BMC Bioinformatics 7, 261. http://www.biomedcentral.com/1471-2105/7/261/abstract

## See Also

An overview of linear model functions in limma is given by 06.LinearModels.

# **Examples**

```
if(require("sma")) {
# Subset of data from ApoAI case study in Limma User's Guide
data(MouseArray)
# Avoid non-positive intensities
RG <- backgroundCorrect(mouse.data, method="half")
MA <- normalizeWithinArrays(RG, mouse.setup)
MA <- normalizeBetweenArrays(MA, method="Aq")
targets <- data.frame(Cy3=I(rep("Pool",6)),Cy5=I(c("WT","WT","WT","KO","KO","KO")))
design <- modelMatrix(targets, ref="Pool")</pre>
```

114 propexpr

```
subarrayw <- printtipWeights(MA, design, layout=mouse.setup)
fit <- lmFit(MA, design, weights=subarrayw)
fit2 <- contrasts.fit(fit, contrasts=c(-1,1))
fit2 <- eBayes(fit2)
# Use of sub-array weights increases the significance of the top genes
topTable(fit2)
# Create an image plot of sub-array weights from each array
zlim <- c(min(subarrayw), max(subarrayw))
par(mfrow=c(3,2), mai=c(0.1,0.1,0.3,0.1))
for(i in 1:6)
imageplot(subarrayw[,i], layout=mouse.setup, zlim=zlim, main=paste("Array", i))
}</pre>
```

propexpr

Estimate Proportion of Expressed Probes

## **Description**

Estimate the proportion of microarray probes which are expressed in each array.

## Usage

```
propexpr(x, neg.x=NULL, status=x$genes$Status, labels=c("negative", "regular"))
```

### **Arguments**

X	matrix or similar object containing raw intensities for a set of arrays.
neg.x	matrix or similar object containing raw intensities for negative control probes for the same arrays. If $\mathtt{NULL}$ , then negative controls must be provided in $\mathtt{x}$ .
status	character vector giving probe types.
labels	character vector giving probe type identifiers.

## Details

This function estimates the proportion of expressed in a microarray by utilizing the negative control probes. Illumina BeadChip arrays contain 750~1600 negative control probes. The expression profile of these control probes can be saved to a separate file by the Illumina BeadStudio software when using it to output the expression profile for regular probes. The control probe profile could be re-generated if it was not generated when the regular probe profile was created by BeadStudio. Other microarray platforms can also use this function to estimate the proportion of expressed probes in each array, provided that they have a set of negative control probes.

labels can include one or two probe type identifiers. Its first element should be the identifier for negative control probes (negative by default). If labels only contains one identifier, then it will be assumed to contain the identifier for negative control probes. By default, regular is the identifier for regular probes.

## Value

Numeric vector giving the proportions of expressed probes in each array.

protectMetachar 115

### Author(s)

Wei Shi and Gordon Smyth

### References

Wei Shi, Carolyn de Graaf, Sarah A Kinkel, Tracey Baldwin, Douglas J Hilton and Gordon K Smyth. Estimating the Proportion of Microarray Probes Expressed in an RNA sample. In preparation.

## See Also

Description to the control probes in Illumina BeadChips can be found in read.ilmn.

# **Examples**

```
## Not run:
x <- read.ilmn(files="sample probe profile.txt",ctrlfiles="control probe profile.txt")
propexpr(x, )
## End(Not run)</pre>
```

protectMetachar

Protect Metacharacters

# **Description**

Add backslashes before any metacharacters found in a string.

# Usage

```
protectMetachar(x)
```

# **Arguments**

Х

character vector

## **Details**

This function is used to protect strings containing metacharacters so that the metacharacters can be treated as ordinary characters in string matching functions operations.

## Value

A character vector of the same length as x in which two backslashes have been inserted before any metacharacter.

## Author(s)

Gordon Smyth

## See Also

An overview of LIMMA functions for reading data is given in 03.ReadingData.

116 *qqt* 

## **Examples**

```
# without protectMetachar, this would be no match
grep(protectMetachar("Ch1 (mean)"),"Ch1 (mean)")
```

qqt

Student's t Quantile-Quantile Plot

# Description

Plots the quantiles of a data sample against the theoretical quantiles of a Student's t distribution.

# Usage

# **Arguments**

У	a numeric vector or array containing the data sample
df	degrees of freedom for the t-distribution. The default ${\tt df=Inf}$ represents the normal distribution.
ylim	plotting range for y
main	main title for the plot
xlab	x-axis title for the plot
ylab	y-axis title for the plot
plot.it	whether or not to produce a plot
	other arguments to be passed to plot

# **Details**

This function is analogous to qqnorm for normal probability plots. In fact qqt (y, df=Inf) is identical to qqnorm(y) in all respects except the default title on the plot.

# Value

A list is invisibly returned containing the values plotted in the QQ-plot:

```
x theoretical quantiles of the t-distribution
y the data sample, same as input y
```

# Author(s)

Gordon Smyth

# See Also

qqnorm

QualityWeights 117

## **Examples**

```
# See also the lmFit examples y \leftarrow rt(50, df=4) qqt(y, df=4) abline(0,1)
```

QualityWeights

Spot Quality Weights

## **Description**

Functions to calculate quality weights for individual spots based on image analyis output file.

# Usage

```
wtarea(ideal=c(160,170))
wtflags(weight=0,cutoff=0)
wtIgnore.Filter
```

# Arguments

ideal numeric vector giving the ideal area or range of areas for a spot in pixels

weight to be given to flagged spots

cutoff cutoff value for Flags below which spots will be downweighted

### Details

These functions can be passed as an argument to read.maimages to construct quality weights as the microarray data is read in.

wtarea downweights unusually small or large spots and is designed for SPOT output. It gives weight 1 to spots which have areas in the ideal range, given in pixels, and linearly downweights spots which are smaller or larger than this range.

wtflags is designed for GenePix output and gives the specified weight to spots with Flags value less than the cutoff value. Choose cutoff=0 to downweight all flagged spots. Choose cutoff=-50 to downweight bad or absent spots or cutoff=-75 to downweight only spots which have been manually flagged as bad.

wtIgnore.Filter is designed for QuantArray output and sets the weights equal to the column Ignore Filter produced by QuantArray. These weights are 0 for spots to be ignored and 1 otherwise.

### Value

A function which takes a dataframe or matrix as argument and produces a numeric vector of weights between 0 and 1

## Author(s)

Gordon Smyth

118 read.columns

#### See Also

An overview of LIMMA functions for reading data is given in 03.ReadingData.

### **Examples**

```
# Read in spot output files from current directory and give full weight to 165
# pixel spots. Note: for this example to run you must set fnames to the names
# of actual spot output files (data not provided).
## Not run:
RG <- read.maimages(fnames, source="spot", wt.fun=wtarea(165))
# Spot will be downweighted according to weights found in RG
MA <- normalizeWithinArrays(RG, layout)
## End(Not run)</pre>
```

read.columns

Read specified columns from a file

# **Description**

Reads specified columns from a file in table format and creates a data frame from it, with cases corresponding to lines and variables to fields in the file.

## Usage

reserved and cannot be set by the user: header, col.names, check.names

# **Arguments**

file the name of the file which the data are to be read from. required.col character vector of names of the required columns text.to.search character string. If any column names can be found in this string, those columns will also be read. the field separator character sep character string of characters to be treated as quote marks quote the number of lines of the data file to skip before beginning to read data. skip logical: if TRUE then in case the rows have unequal length, blank fields are fill implicitly added. blank.lines.skip logical: if TRUE blank lines in the input are ignored. comment.char character: a character vector of length one containing a single character or an empty string. allowEscapes logical. Should C-style escapes such as '\n' be processed or read verbatim (the default)? other arguments are passed to read.table, excluding the following which are . . .

and colClasses.

readGAL 119

## **Details**

This function is an interface to read.table in the base package. It uses required.col and text.to.search to set up the colClasses argument of read.table.

Note the following arguments of read.table are used by read.columns and therefore cannot be set by the user: header, col.names, check.names and colClasses.

This function is used by read.maimages.

# Value

A data frame (data.frame) containing a representation of the data in the file.

# Author(s)

Gordon Smyth

## See Also

```
read.maimages, read.table.
```

An overview of LIMMA functions for reading data is given in 03.ReadingData.

readGAL Read a GAL file

# **Description**

Read a GenePix Array List (GAL) file into a dataframe.

# Usage

readGAL(galfile=NULL,path=NULL,header=TRUE,sep="\t",quote="\"",skip=NULL,as.is=T

# **Arguments**

galfile	character string giving the name of the GAL file. If $\mathtt{NULL}$ then a file with extension .gal is found in the directory specified by $\mathtt{path}$ .
path	character string giving the directory containing the files. If ${\tt NULL}$ then assumed to be the current working directory.
header	logical variable, if TRUE then the first line after $skip$ is assumed to contain column headings. If FALSE then a value should specified for $skip$ .
sep	the field separator character
quote	the set of quoting characters
skip	number of lines of the GAL file to skip before reading data. If ${\tt NULL}$ then this number is determined by searching the file for column headings.
as.is	logical variable, if ${\tt TRUE}$ then read in character columns as vectors rather than factors.
• • •	any other arguments are passed to read.table

120 readHeader

#### **Details**

A GAL file is a list of genes IDs and associated information produced by an Axon microarray scanner. Apart from header information, the file must contain data columns labeled Block, Column, Row and ID. A Name column is usually included as well. Other columns are optional. See the Axon URL below for a detaile description of the GAL file format.

This function reads in the data columns with a minimum of user information. In most cases the function can be used without specifying any of the arguments.

### Value

A data frame with columns

Block numeric vector containing the print tip indices

Column numeric vector containing the spot columns

Row numeric vector containing the spot rows

ID character vector, for factor if as.is=FALSE, containing gene library identi-

fiers

Name character vector, for factor if as .is=FALSE, containing gene names

The data frame will be sorted so that Column is the fastest moving index, then Row, then Block.

### Author(s)

Gordon Smyth

# References

```
http://www.axon.com/gn_GenePix_File_Formats.html
```

## See Also

read.Galfile in the marray package.

An overview of LIMMA functions for reading data is given in 03.ReadingData.

# **Examples**

```
# readGAL()
# will read in the first GAL file (with suffix ".gal")
# found in the current working directory
```

readHeader

Read Header Information from Image Analysis Raw Data File

# Description

Read the header information from a GenePix Results (GPR) file or from an SMD raw data file. These functions are used internally by read.maimages and are not usually called directly by users.

readHeader 121

### Usage

```
readGenericHeader(file, columns, sep="\t")
readGPRHeader(file)
readSMDHeader(file)
```

## **Arguments**

file character string giving file name. If it does not contain an absolute path, the file

name is relative to the current working directory.

columns character vector specifying data column headings expected to be in file

sep the character string separating column names

## **Details**

Raw data files exported by image analysis programs include a number of header lines which contain information about the scanning process. This function extracts that information and locates the line where the intensity data begins. readGPRHeader is for GenePix output and readSMDHeader is for files from the Stanford Microarray Database (SMD). readGenericHeader finds the line in the file on which the data begins by searching for specified column headings.

## Value

A list with components corresponds to lines of header information. A key component is NHeaderRecords which gives the number of lines in the file before the intensity data begins. All other components are character vectors.

## Author(s)

Gordon Smyth

## References

```
See http://www.axon.com/gn_GenePix_File_Formats.html for GenePix formats.

See http://www.bluegnome.co.uk for information on BlueFuse.

See http://genome-www.stanford.edu/Microarray for the SMD.
```

# See Also

```
read.maimages
```

An overview of LIMMA functions to read data is given in 03.ReadingData.

122 read.ilmn

read.ilmn Read Illumina Expression Data	
---	--

# Description

Read Illumina summary probe profile files and summary control probe profile files

## Usage

```
read.ilmn(files=NULL, ctrlfiles=NULL, path=NULL, ctrlpath=NULL,
probeid="Probe", annotation=c("TargetID", "SYMBOL"), expr="AVG_Signal", other.c
sep="\t", verbose=TRUE, ...)
```

## **Arguments**

files	character vector giving the names of the summary probe profile files.	
ctrlfiles	character vector giving the names of the summary control probe profile files.	
path	character string giving the directory containing the summary probe profile files. The default is the current working directory.	
ctrlpath	character string giving the directory containing the summary control probe profile files. The default is the current working directory.	
probeid	character string giving the name of the probe identifier column.	
annotation	character vector giving possible names of the annotation column. It could be called "TargetID" or "SYMBOL" depending on which version of BeadStudio is	
	used.	
expr	character string giving the keyword in the names of the expression intensity columns.	
other.column	as a second of the second of t	
	character vector giving the keywords in the names of extra columns required, such as "Detection", "Avg_NBEADS", "BEAD_STDEV" etc. Each keyword corresponds to one type of columns.	
sep	the field separator character.	
verbose	logical, TRUE to report names of profile files being read.	
	any other parameters are passed on to read.columns.	

## **Details**

Illumina BeadStudio ouputs probe intensities (regular probe intensities) and control probe intensities to summary probe profile files (containing regular probes) and summary control probe profile files, respectively. If both files and ctrlfiles are not NULL, this function will combine the data read from the two file types and save them to an EListRaw-class object. If one of them is NULL, then only the required data are read in.

Probe types are indicated in the Status column of genes, a component of the returned EListRaw-class object. There are totally seven types of control probes including negative, biotin, labeling, cy3\_hyb, housekeeping, high\_stringency\_hyb or low\_stringency\_hyb. Regular probes have the probe type regular. The Status column will not be created if ctrlfiles is NULL.

To read in columns other than probeid, annotation and expr, users needs to specify keywords in other.columns. One keyword corresponds to one type of columns. Examples of keywords are "Detection", "Avg\_NBEADS", "BEAD\_STDEV" etc.

read.ilmn.targets 123

### Value

An EListRaw-class object with the following components:

E numeric matrix of raw intensities.

genes data.frame of probe annotation.

targets data.frame of sample information.

other list of other column data.

### Author(s)

Wei Shi and Gordon K Smyth

### See Also

read.ilmn.targets reads in Illumina expression data using the file information extracted from a target data frame which is often created by the readTargets function.

neqc performs normexp by control background correction, log transformation and quantile between-array normalization for Illumina expression data.

normexp.fit.control estimates the parameters of the normal+exponential convolution model with the help of negative control probes.

propexpr estimates the proportion of expressed probes in a microarray.

# **Examples**

```
## Not run:
x <- read.ilmn(files="sample probe profile.txt",ctrlfiles="control probe profile.txt")
## End(Not run)</pre>
```

```
read.ilmn.targets Read Illumina Data from a Target Dataframe
```

# Description

Read Illumina data from a target dataframe

# Usage

```
read.ilmn.targets(targets, ...)
```

# Arguments

```
data frame including names of profile files.any other parameters are passed on to read.ilmn.
```

124 readImaGeneHeader

#### **Details**

targets is often created by calling the function readTargets. Rows in targets are arrays and columns contain related array or RNA sample information.

At least one of the two columns called files and/or ctrlfiles should be present in targets, which includes names of summary probe profile files and names of summary control probe profile files respectively. This function calls read.ilmn to read in the data.

## Value

An EListRaw-class object. See return value of the function read.ilmn for details.

### Author(s)

Wei Shi

### See Also

read.ilmn

readImaGeneHeader Read ImaGene Header Information

## **Description**

Read the header information from an ImaGene image analysis output file. This function is used internally by read.maimages and is not usually called directly by users.

# Usage

```
readImaGeneHeader(file)
```

## Arguments

file

character string giving file name or path

### **Details**

The raw data files exported by the image analysis software ImaGene include a number of header lines which contain information about the printing and scanning processes. This function extracts that information and locates the line where the intensity data begins.

## Value

A list containing information read from the header of the ImaGene file. Each Begin-End environment found in the file header will become a recursive list in the output object, with components corresponding to fields in the file. See the ImaGene documentation for further information. The output object will also contain a component NHeaderRecordsgiving the number of lines in the file before the intensity data begins.

## Author(s)

Gordon Smyth

read.maimages 125

### References

```
http://www.biodiscovery.com/imagene.asp
```

## See Also

```
read.imagene
```

An overview of LIMMA functions to read data is given in 03.ReadingData.

# **Examples**

```
## Not run:
h <- readImaGeneHeader("myImaGenefile.txt")
names(h)
h$NHeaderRecords
h[["Field Dimensions"]]
## End(Not run)</pre>
```

read.maimages

Read RGList or EListRaw from Image Analysis Output Files

# **Description**

Reads an RGList from a set of two-color microarray image analysis output files, or an EListRaw from a set of one-color files.

# Usage

# Arguments

files	character vector giving the names of the files containing image analysis output or, for Imagene data, a character matrix of names of files. If omitted, then all files with extension ext in the specified directory will be read in alphabetical order.
source	character string specifying the image analysis program which produced the output files. Choices are "generic", "agilent", "agilent.median", "arrayvision", arrayvision.ARM"}, arrayvision.MTM", "bluefuse", "genepix", "genepix.custom", "genepix.median", "imagene", "quantarray", "scanarrayexpress", "smd.old", "smd", "spot" or "spot.close.open".
path	character string giving the directory containing the files. The default is the current working directory.
ext	character string giving optional extension to be added to each file name

126 read.maimages

names character vector of names to be associated with each array as column name.

Defaults to removeExt (files).

columns list with fields R, G, Rb and Gb giving the column names to be used for red and

green foreground and background or, in the case of Imagene data, a list with fields f and b. This argument is optional if source is specified, otherwise it is

required.

other.columns

character vector of names of other columns to be read containing spot-specific

information

annotation character vector of names of columns containing annotation information about

the probes

channels number of channels, channels=2 or two-color and channels=1 for one-

color.

wt.fun function to calculate spot quality weights
verbose logical, TRUE to report each time a file is read

sep the field separator character

quote character string of characters to be treated as quote marks

... any other arguments are passed to read.table

### **Details**

This is the main data input function for the LIMMA package for two-color microarray data. It extracts the foreground and background intensities from a series of files, produced by an image analysis program, and assembles them into the components of one list. The image analysis programs Agilent Feature Extraction, ArrayVision, BlueFuse, GenePix, ImaGene, QuantArray (Version 3 or later), Stanford Microarray Database (SMD) and SPOT are supported explicitly. Data from some other image analysis programs can be read if the appropriate column names containing the foreground and background intensities are specified using the columns argument. (This will work if the column names are unique and if there are no rows in the file after the last line of data. Header lines are ok.)

SMD data should consist of raw data files from the database, in tab-delimited text form. There are two possible sets of column names depending on whether the data was entered into the database before or after September 2003. source="smd.old" indicates that column headings in use prior to September 2003 should be used. In the case of Agilent and GenePix, two possible foreground estimators are supported: source="genepix" uses the mean foreground estimates while source="genepix.median" uses median foreground estimates. Similarly for Agilent. GenePix 6.0 and later also supplies some custom background options, notably morphological background. If the GPR files have been written using a custom background, you may read it using source="genepix.custom". In the case of SPOT, two possible background estimators are supported: if source="spot.close.open" then background intensities are estimated from morph.close.open rather than morph.

Spot quality weights may be extracted from the image analysis files using a weight function wt.fun. wt.fun may be any user-supplied function which accepts a data.frame argument and returns a vector of non-negative weights. The columns of the data.frame are as in the image analysis output files. There is one restriction, which is that the column names should be referred to in full form in the weight function, i.e., do not rely on name expansion for partial matches when refering to the names of the columns. See QualityWeights for suggested weight functions.

For Imagene image data the argument files should be a matrix with two columns. The first column should contain the names of the files containing green channel (cy3) data and the second column should contain names of files containing red channel (cy5) data. If source="imagene"

read.maimages 127

and files is a vector of even length instead of a matrix, then each consecutive pair of file names is assumed to correspond to the same array. The function read.imagene is called by read.maimages when source="imagene". It does not need to be called directly by users.

ArrayVision reports spot intensities in a number of different ways. read.maimages caters for ArrayVision's Artifact-removed (ARM) density values as "arrayvision.ARM" or for Median-based Trimmed Mean (MTM) density values as "arrayvision.MTM". ArrayVision users may find it useful to read the top two lines of their data file to check which version of density values they have.

The argument other.columns allows arbitrary columns of the image analysis output files to be preserved in the data object. These become matrices in the component other component. For ImaGene data, the other column headings with be prefixed with "R " or "G " as appropriate.

## Value

For one-color data, an EListRaw object. For two-color data, an RGList object containing the components

R	matrix containing the red channel foreground intensities for each spot for each array.	
Rb	matrix containing the red channel background intensities for each spot for each array.	
G	matrix containing the green channel foreground intensities for each spot for each array.	
Gb	matrix containing the green channel background intensities for each spot for each array.	
weights	spot quality weights, if wt.fun is given	
other	list containing matrices corresponding to other.columns if given	
genes	data frame containing annotation information about the probes, for example gene names and IDs and spatial positions on the array, currently set only if source is "agilent", "genepix" or source="imagene" or if the annotation argument is set	
targets	data frame with column FileName giving the names of the files read	
source	character string giving the image analysis program name	
printer	list of class PrintLayout, currently set only if source="imagene"	

## Warnings

All image analysis files being read are assumed to contain data for the same genelist in the same order. No checking is done to confirm that this is true. Probe annotation information is read from the first file only.

## Author(s)

Gordon Smyth, with speed improvements by Marcus Davy

# References

Web pages for the image analysis software packages mentioned here are listed at http://www.statsci.org/micrarra/image.html

128 readSpotTypes

### See Also

read.maimages uses read.columns for efficient reading of text files. As far as possible, it is has similar behavior to read.table in the base package.

An overview of LIMMA functions for reading data is given in 03.ReadingData.

# **Examples**

```
# Read all .gpr files from current working directory
# and give weight 0.1 to spots with negative flags

## Not run: files <- dir(pattern="*\\.gpr$")
RG <- read.maimages(files, "genepix", wt.fun=wtflags(0.1))
## End(Not run)

# Read all .spot files from current working director and down-weight
# spots smaller or larger than 150 pixels

## Not run: files <- dir(pattern="*\\.spot$")
RG <- read.maimages(files, "spot", wt.fun=wtarea(150))
## End(Not run)</pre>
```

readSpotTypes

Read Spot Types File

# Description

Read a table giving regular expressions to identify different types of spots in the gene-dataframe.

# Usage

```
readSpotTypes(file="SpotTypes.txt",path=NULL,sep="\t",check.names=FALSE,...)
```

## **Arguments**

file	character string giving the name of the file specifying the spot types.
path	character string giving the directory containing the file. Can be omitted if the file is in the current working irectory.
sep	the field separator character
check.names	logical, if FALSE column names will not be converted to valid variable names, for example spaces in column names will not be left as is
	any other arguments are passed to read.table

### **Details**

The file is a text file with rows corresponding to types of spots and the following columns: SpotType gives the name for the spot type, ID is a regular expression matching the ID column, Name is a regular expression matching the Name column, and Color is the R name for the color to be associated with this type.

readTargets 129

#### Value

### A data frame with columns

SpotType character vector giving names of the spot types

ID character vector giving regular expressions

Name character vector giving regular expressions

Color character vector giving names of colors

## Author(s)

Gordon Smyth following idea of James Wettenhall

### See Also

An overview of LIMMA functions for reading data is given in 03.ReadingData.

# **Description**

Read targets file for a microarray experiment into a dataframe.

# Usage

```
readTargets(file="Targets.txt", path=NULL, sep="\t", row.names=NULL, quote="\"",
```

# Arguments

file	character:	string	giving	the name	of the	targets file.
T T T C	CHaracter 1	ou me	Z1 V 111 Z	tiic mamic	OI LIIC	turgetts me.

path character string giving the directory containing the file. Can be omitted if the

file is in the current working irectory.

sep field separator character

row.names character string giving the name of a column from which to obtain row names

quote the set of quoting characters

... other arguments are passed to read.table

### **Details**

The targets file is a text file containing information about the RNA samples used as targets in the microarray experiment. Rows correspond to arrays and columns to covariates associated with the targets. For a two-color experiment, the targets file will normally include columns labelled Cy3 and Cy5 or similar specifying which RNA samples are hybridized to each channel of each array. Other columns may contain any other covariate information associated with the arrays or targets used in the experiment.

If row.names is non-null and there is a column by that name with unique values, then those values will be used as row names for the dataframe. If row.names is null, then the column Labels will be used if such exists or, failing that, the column FileName.

See the Limma User's Guide for examples of this function.

130 removeBatchEffect

#### Value

A dataframe. Character columns are not converted into factors.

### Author(s)

Gordon Smyth

### See Also

An overview of LIMMA functions for reading data is given in 03.ReadingData.

removeBatchEffect Remove Batch Effect

# **Description**

Remove a batch effect from expression data.

## Usage

```
removeBatchEffect(x,batch,design=matrix(1,ncol(x),1))
```

# **Arguments**

x numeric matrix containing log-expression intensity for a series of microarrays.

Each column corresponds to an array.

batch a factor or vector indicating batches.

design design matrix

# Details

This function is useful for removing batch effects, associated with hybridization time or other technical variables, prior to clustering or unsupervised analysis.

This function is intended for use with clustering or PCA, not for use prior to linear modelling. If linear modelling is intended, it is better to include the batch effect as part of the linear model.

The design matrix is used to describe comparisons between the samples, for example treatment effects, which should not be removed. The batch effect is estimated after adjusting for the design factors. By default, this is just the grand mean.

## Value

A numeric matrix of log-expression values with batch effects removed.

# Author(s)

Gordon Smyth and Carolyn de Graaf

# See Also

05.Normalization

removeExt 131

removeExt

Remove Common Extension from File Names

# **Description**

Finds and removes any common extension from a vector of file names.

# Usage

```
removeExt(x)
```

# **Arguments**

Х

character vector

## Value

A character vector of the same length as x in which any common extension has been stripped off.

## Author(s)

Gordon Smyth

## See Also

An overview of LIMMA functions for reading data is given in 03.ReadingData.

# **Examples**

```
x \leftarrow c("slide1.spot", "slide2.spot", "slide3.spot") removeExt(x)
```

 $\verb"residuals.MArrayLM" \textit{Extract Residuals from MArrayLM Fit}$ 

# **Description**

This method extracts the residuals from all the probewise linear model fits and returns them in a matrix.

# Usage

```
## S3 method for class 'MArrayLM':
residuals(object, y, ...)
```

# **Arguments**

```
object a fitted model object inheriting from class MarrayLM.

y a data object containing the response data used to compute the fit. This can be of any class for which as .matrix is defined, including MAList, ExpressionSet, marrayNorm etc.

... other arguments are not used
```

132 RGList-class

#### Value

Numeric matrix of residuals.

#### See Also

residuals.

RGList-class

Red, Green Intensity List - class

## **Description**

A simple list-based class for storing red and green channel foreground and background intensities for a batch of spotted microarrays. RGList objects are normally created by read.maimages.

## **Slots/List Components**

RGList objects can be created by new ("RGList", RG) where RG is a list. Objects of this class contains no slots (other than .Data), but objects should contain the following list components:

- numeric matrix containing the red (cy5) foreground intensities. Rows correspond to spots and columns to arrays.
- numeric matrix containing the green (cy3) foreground intensities. Rows correspond to spots and columns to arrays.

# Optional components include

Rb: numeric matrix containing the red (cy5) background intensities numeric matrix containing the green (cy3) background intensities Gb:

numeric matrix of same dimension as R containing relative spot quality weights. Elements should be non-ne weights:

list containing other matrices, all of the same dimensions as R and G. other:

data.frame containing probe information. Should have one row for each spot. May have any number of colu data.frame containing information on the target RNA samples. Rows correspond to arrays. May have any nu targets: printer: list containing information on the process used to print the spots on the arrays. See PrintLayout.

Valid RGList objects may contain other optional components, but all probe or array information should be contained in the above components.

### Methods

genes:

This class inherits directly from class list so any operation appropriate for lists will work on objects of this class. In addition, RGList objects can be subsetted, combined and merged. RGList objects will return dimensions and hence functions such as dim, nrow and ncol are defined. RGLists also inherit a show method from the virtual class LargeDataObject, which means that RGLists will print in a compact way.

RGList objects can be converted to exprSet2 objects by as (RG, "exprSet2").

Other functions in LIMMA which operate on RGList objects include normalizeBetweenArrays, normalizeForPrintorder, normalizeWithinArrays.

## Author(s)

Gordon Smyth

roast 133

### See Also

02.Classes gives an overview of all the classes defined by this package.

marrayRaw is the corresponding class in the marray package.

roast Rotation Gene Set Tests

# Description

Rotation gene set testing for linear models.

# Usage

roast(iset=NULL, y, design, contrast=ncol(design), gene.weights=NULL, array.weights=NULL,

## **Arguments**

vector specifying the rows of $y$ in the test set. This can be a vector of indices, or a logical vector of the same length as statistics, or any vector such as $y[selected,]$ contains the values for the gene set to be tested.
numeric matrix giving log-expression values. If var.prior or df.prior are null, then y should contain values for all genes on the arrays. If both prior parameters are given, then only y values for the test set are required.
design matrix
contrast for which the test is required. Can be an integer specifying a column of design, or else a contrast vector of length equal to the number of columns of design.

gene.weights optional numeric vector of weights for genes in the set.

array.weights

optional numeric vector of array weights.

block optional vector of blocks.
correlation correlation between blocks.

var.prior prior value for residual variances. If not provided, this is estimated from all the

data using squeezeVar.

df.prior prior degrees of freedom for residual variances. If not provided, this is estimated

using squeezeVar.

nrot number of rotations used to estimate the p-values.

## **Details**

This function tests whether any of the genes in the set are differentially expressed. It uses rotation, which is a smoothed version of permutation suitable for linear models (Langsrud, 2005). It can be used for any linear model with replication. and negative values and otherwise will be taken to be F-like.

This is a self-contained test in the sense that genes outside the test set do not play a role (Goeman and Buhlmann, 2007). A competitive gene set test is performed by geneSetTest, for one set, or romer, for many sets.

134 roast

p-values are given for four possible alternative hypotheses. alternative=="up" means the genes in the set tend to be up-regulated, with positive t-statistics. alternative=="down" means the genes in the set tend to be down-regulated, with negative t-statistics. alternative=="either" means the set is either up or down-regulated as a whole. alternative=="mixed" test whether the genes in the set tend to be differentially expressed, without regard for direction. In this case, the test will be significant if the set contains mostly large test statistics, even if some are positive and some are negative.

The first three alternatives are appropriate if you have a prior expectation that all the genes in the set will react in the same direction. The "mixed" alternative is appropriate if you know only that the genes are involved in the relevant pathways, without knowing the direction of effect for each gene. The "mixed" alternative is the only one possible with F-like statistics.

Note that roast estimates p-values by simulation, specifically by random rotations of the orthogonalized residuals. This means that the p-values will vary slightly from run to run. To get more precise p-values, increase the number of rotations nrot. The strategy of random rotations is due to Langsrud (2005).

Following Monte Carlo hypothesis testing theory (Barnard, 1963), the p-value is computed as (x+1) / (nrot+1) where x is the number of rotations giving a more extreme statistic than that observed. This means that the smallest possible p-value is 1/(nrot+1).

### Value

data.frame with columns Z, Active and P. Value. The Z column gives average (root mean square) z-statistics for the genes in the set. The Active gives the proportion of genes in the set contributing meaningfully to significance, defined as those with squared z-values greater than 2. The P. Value gives estimated p-values. The rows correspond to the alternative hypotheses mixed, up, down or either.

# Author(s)

Gordon Smyth and Di Wu

### References

Barnard, GA (1963). Discussion of The spectral analysis of point processes (by MS Bartlett). *Journal of the Royal Statistical Society* B 25, 294.

Goeman, JJ, and Buhlmann, P 2007. Analyzing gene expression data in terms of gene sets: methodological issues. *Bioinformatics* 23, 980-987.

Langsrud, O. (2005). Rotation tests. Statistics and Computing 15, 53-60

## See Also

```
wilcoxGST, romer
```

An overview of tests in limma is given in 08. Tests.

# **Examples**

```
y <- matrix(rnorm(100*4),100,4)
design <- cbind(Intercept=1,Group=c(0,0,1,1))
iset <- 1:5
y[iset,3:4] <- y[iset,3:4]+3
roast(iset,y,design,contrast=2)</pre>
```

romer 135

```
# Alternative approach useful if multiple gene sets are tested:
fit <- lmFit(y,design)
sv <- squeezeVar(fit$sigma^2,df=fit$df.residual)
iset1 <- 1:5
iset2 <- 6:10
roast(y=y[iset1,],design=design,contrast=2,var.prior=sv$var.prior,df.prior=sv$var.prior)
roast(y=y[iset2,],design=design,contrast=2,var.prior=sv$var.prior,df.prior=sv$var.prior)</pre>
```

romer

Rotation Gene Set Enrichment Analysis

## **Description**

Rotation-mean-rank version of GSEA (gene set enrichment analysis) for linear models.

# Usage

```
romer(iset,y,design,contrast=ncol(design),array.weights=NULL,block=NULL,correlat
romer2(iset,y,design,contrast=ncol(design),array.weights=NULL,block=NULL,correlat
```

## **Arguments**

iset list of indices specifying the rows of y in the gene sets. The list can be made

by symbols2indices and the gene sets can be retrieved from the molecular sig-

natures database in Broad Institute.

y numeric matrix giving log-expression values.

design matrix

contrast contrast for which the test is required. Can be an integer specifying a column of

design, or else a contrast vector of length equal to the number of columns of

design.

array.weights

optional numeric vector of array weights.

block optional vector of blocks.

correlation correlation between blocks.

nrot number of rotations used to estimate the p-values.

floor logical, if TRUE use floor-mean statistic instead of mean-rank statistic for sum-

marizing t-statistics

## **Details**

This function implements a GSEA of a battery of gene sets similar in motivation to Subramanian et al (2005) but designed for use with linear models. It is a competitive test, in that the different gene sets are pitted against one another. Instead of permutation, it uses rotation, a parametric resampling method suitable for linear models (Langsrud, 2005). romer can be used with any linear model with some level of replication.

Number of genes is given for each gene set.

p-values are given for four possible alternative hypotheses. The alternative "up" means the genes in the set tend to be up-regulated, with positive t-statistics. The alternative "down" means the genes in the set tend to be down-regulated, with negative t-statistics. The alternative "either" means the

136 romer

set is either up or down-regulated as a whole. The alternative "mixed" test whether the genes in the set tend to be differentially expressed, without regard for direction. In this case, the test will be significant if the set contains mostly large test statistics, even if some are positive and some are negative.

The first three alternatives are appropriate if you have a prior expection that all the genes in the set will react in the same direction. The "mixed" alternative is appropriate if you know only that the genes are involved in the relevant pathways, without knowing the direction of effect for each gene. The "mixed" alternative is the only one possible with F-like statistics.

Note that romer estimates p-values by simulation, specifically by random rotations of the orthogonalized residuals. This means that the p-values will vary slightly from run to run. To get more precise p-values, increase the number of rotations nrot. The strategy of random rotations is due to Langsrud (2005).

romer and romer2 differ in the way that t-statistics are summarized to form a set-level test statistic. The genes are ranked by moderated t-statistic and romer uses the mean-ranks for each set as the statistic. If floor=TRUE then negative t-statistics are put to zero before ranking for the up test, and vice versa for the down test. This improves the power for detecting genes with a subset of responding genes. romer2 using the mean rank of the top 50 This statistic performs well in practice but is slightly slower to compute.

#### Value

Matrix with the rows corresponding to estimated p-values for each gene set and the columns corresponding to the number of genes for each gene set and the alternative hypotheses mixed, up, down or either.

#### Author(s)

Yifang Hu and Gordon Smyth

# References

Langsrud, O, 2005. Rotation tests. Statistics and Computing 15, 53-60

Subramanian, A, Tamayo, P, Mootha, VK, Mukherjee, S, Ebert, BL, Gillette, MA, Paulovich, A, Pomeroy, SL, Golub, TR, Lander, ES and Mesirov JP, 2005. Gene set enrichment analysis: a knowledge-based approach for interpreting genome-wide expression profiles. *Proc Natl Acad Sci U S A* 102, 15545-15550

# See Also

```
roast, wilcoxGST
```

An overview of tests in limma is given in 08. Tests.

# Examples

```
y <- matrix(rnorm(100*4),100,4)
design <- cbind(Intercept=1,Group=c(0,0,1,1))
iset <- 1:5
y[iset,3:4] <- y[iset,3:4]+3

iset1 <- 1:5
iset2 <- 6:10
romer(iset=list(iset1=iset1,iset2=iset2),y=y,design=design,contrast=2,nrot=99)</pre>
```

selectModel 137

selectModel Select Appropriate Linear Model
---

# Description

Select the best fitting linear model for each gene by minimizing an information criterion.

# Usage

```
selectModel(y, designlist, criterion="aic", df.prior=0, s2.prior=NULL, s2.true=N
```

# Arguments

У	a matrix-like data object, containing log-ratios or log-values of expression for a series of microarrays. Any object class which can be coerced to matrix is acceptable including numeric, matrix, MAList, marrayNorm, ExpressionSet or PLMset.
designlist	list of design matrices
criterion	information criterion to be used for model selection, "aic", "bic" or "mallowscp".
df.prior	prior degrees of freedom for residual variances. See squeezeVar
s2.prior	prior value for residual variances, to be used if df.prior>0.
s2.true	numeric vector of true variances, to be used if criterion="mallowscp".
	other optional arguments to be passed to lmFit

## **Details**

This function chooses, for each probe, the best fitting model out of a set of alternative models represented by a list of design matrices. Selection is by Akaike's Information Criterion (AIC), Bayesian Information Criterion (BIC) or by Mallow's Cp.

The criteria have been generalized slightly to accommodate an information prior on the variances represented by s2.prior and df.prior or by s2.post. Suitable values for these parameters can be estimated using squeezeVar.

## Value

List with components

IC	matrix of information criterion scores, rows for probes and columns for models
pref	factor indicating the model with best (lowest) information criterion score

# Author(s)

Alicia Oshlack and Gordon Smyth

## See Also

An overview of linear model functions in limma is given by 06.LinearModels.

138 squeezeVar

### **Examples**

```
nprobes <- 100
narrays <- 5
y <- matrix(rnorm(nprobes*narrays), nprobes, narrays)
A <- c(0,0,1,1,1)
B <- c(0,1,0,1,1)
designlist <- list(
   None=cbind(Int=c(1,1,1,1,1)),
   A=cbind(Int=1,A=A),
   B=cbind(Int=1,A=A),
   Both=cbind(Int=1,A=A,B=B),
   Add=cbind(Int=1,A=A,B=B),
   Full=cbind(Int=1,A=A,B=B),
   Full=cbind(Int=1,A=A,B=B,AB=A*B))
out <- selectModel(y,designlist)
table(out$pref)</pre>
```

squeezeVar

Smooth Sample Variances

## **Description**

Smooth a set of sample variances by computing empirical Bayes posterior means.

### Usage

```
squeezeVar(var, df)
```

## **Arguments**

var numeric vector of independent sample variances

df numeric vector of degrees of freedom for the sample variances

## **Details**

The sample variances var are assumed to follow scaled chi-squared distributions. An inverse chi-squared prior is assumed for the true variances. The scale and degrees of freedom for the prior distribution are estimated from the data.

The effect of this function is to smooth or shrink the variances towards a common value. The smoothed variances have a smaller expected mean square error to the true variances than do the sample variances themselves.

This function is called by <code>eBayes</code>, but beware a possible confusion with the output from that function. The values <code>var.prior</code> and <code>var.post</code> output by <code>squeezeVar</code> correspond to the quantities <code>s2.prior</code> and <code>s2.post</code> output by <code>eBayes</code>, whereas <code>var.prior</code> output by <code>eBayes</code> relates to a different parameter.

## Value

## A list with components

var.post numeric vector of posterior variances
var.prior location of prior distribution
df.prior degrees of freedom of prior distribution

strsplit2

## Author(s)

Gordon Smyth

### References

Smyth, G. K. (2004). Linear models and empirical Bayes methods for assessing differential expression in microarray experiments. *Statistical Applications in Genetics and Molecular Biology*, **3**, No. 1, Article 3. http://www.bepress.com/sagmb/vol3/iss1/art3

## **Examples**

```
s2 <- rchisq(20,df=5)/5
squeezeVar(s2, df=5)</pre>
```

strsplit2

Split Composite Names

# Description

Split a vector of composite names into a matrix of simple names.

# Usage

```
strsplit2(x, split, ...)
```

## **Arguments**

```
x character vector
split character to split each element of vector on, see strsplit
... other arguments are passed to strsplit
```

# **Details**

This function is the same as strsplit except that the output value is a matrix instead of a list. The first column of the matrix contains the first component from each element of x, the second column contains the second components etc. The number of columns is equal to the maximum number of components for any element of x.

The motivation for this function in the limma package is handle input columns which are composites of two or more annotation fields.

### Value

A list containing components

Name character vector of the same length as x contain first splits of each element character vector of the same length as x contain second splits of each element

## Author(s)

## Gordon Smyth

subsetting

## See Also

```
strsplit.
```

An overview of LIMMA functions for reading data is given in 03.ReadingData.

# **Examples**

```
x <- c("AA196000; actinin, alpha 3", "AA464163; acyl-Coenzyme A dehydrogenase, very long chain", "3E7; W15277; No Annotation") strsplit2(x, split=";")
```

subsetting

Subset RGList, MAList, EList or MArrayLM Objects

# Description

Extract a subset of an RGList, MAList, EList or MArrayLM object.

# Usage

```
## S3 method for class 'RGList':
object[i, j, ...]
```

# **Arguments**

```
object of class RGList, MAList, EList or MArrayLM

i, j elements to extract. i subsets the probes or spots while j subsets the arrays

not used
```

## Details

i, j may take any values acceptable for the matrix components of object. See the Extract help entry for more details on subsetting matrices.

### Value

An object of the same class as object holding data from the specified subset of genes and arrays.

# Author(s)

Gordon Smyth

## See Also

```
Extract in the base package.
```

03.ReadingData gives an overview of data input and manipulation functions in LIMMA.

summary 141

# **Examples**

```
M <- A <- matrix(11:14,4,2)
rownames(M) <- rownames(A) <- c("a","b","c","d")
colnames(M) <- colnames(A) <- c("A","B")
MA <- new("MAList",list(M=M,A=A))
MA[1:2,]
MA[1:2,2]
MA[,2]</pre>
```

summary

Summaries of Microarray Data Objects

# Description

Briefly summarize microarray data objects.

# Usage

```
## S3 method for class 'RGList':
summary(object, ...)
```

# Arguments

```
object an object of class RGList, MAList or MArrayLM ... other arguments are not used
```

## **Details**

The data objects are summarized as if they were lists, i.e., brief information about the length and type of the components is given.

## Value

A table.

# Author(s)

Gordon Smyth

## See Also

```
summary in the base package.
```

02. Classes gives an overview of data classes used in LIMMA.

142 targetsA2C

symbols2indices

Make Gene Set Indices from Symbols

## **Description**

Make a list of gene sets symbols into a list of gene sets indices used to create input for romer function.

# Usage

```
symbols2indices(gmtl.official, symbol)
```

## **Arguments**

gmtl.official

list of gene sets containing official gene symbols. The gene sets can be the ones from the molecular signatures database in Broad Institute.

symbol

vector of symbols for the probesets.

## **Details**

This function takes a list of official gene sets symbols and returns a list of indices of them according to the vector symbol.

## Value

list of indices.

## Author(s)

Gordon Smyth and Yifang Hu

# See Also

romer

targetsA2C

Convert Two-Color Targets Dataframe from One-Row-Per-Array to One-Row-Per-Channel

# Description

Convert a two-color targets dataframe with one row per array to one with one row per channel.

# Usage

```
targets \texttt{A2C} (targets, \ channel.codes = \texttt{c(1,2)}, \ channel.columns = \texttt{list} (\texttt{Target} = \texttt{c("Cy3","Columns}) = \texttt{columns} = \texttt{
```

targetsA2C 143

## **Arguments**

targets data.frame with one row per array giving information about target samples as-

sociated covariates.

channel.codes

numeric or character vector of length 2 giving codes for the channels

channel.columns

named list of character vectors of length 2. Each entry gives a pair of names of columns in targets which contain channel-specific information. This pair of

columns should be assembled into one column in the output.

grep logical, if TRUE then the channel column names are found by greping, i.e., the

actual column names need only contain the names given by channel.columns

as substrings

### **Details**

The targets dataframe holds information about the RNA samples used as targets in the microarray experiment. It is often read from a file using readTargets. This function is used to convert the dataframe from an array-orientated format with one row for each array and two columns for the two channels into a channel-orientated format with one row for each individual channel observations. In statistical terms, the first format treats the arrays as cases and treats the channels as repeated measurements. The second format treats the individual channel observations as cases. The second format may be more appropriate if the data is to be analyzed in terms of individual log-intensities.

# Value

data.frame with twice as many rows as targets. Any pair of columns named by channel.columns will now be one column.

### Author(s)

Gordon Smyth

## See Also

targetsA2C is used by the coerce method from RGList for ExpressionSet in the convert package.

An overview of methods for single channel analysis in limma is given by 07. Single Channel.

# **Examples**

```
targets <- data.frame(FileName=c("file1.gpr","file2.gpr"),Cy3=c("WT","KO"),Cy5=c("KO","WT
targetsA2C(targets)
```

144 TestResults-class

```
TestResults-class Matrix of Test Results - class
```

## **Description**

A matrix-based class for storing the results of simultaneous tests. TestResults objects are normally created by classifyTestsF, classifyTestsT or classifyTestsP.

# Usage

```
## S3 method for class 'TestResults':
summary(object, ...)
```

## **Arguments**

```
object of class TestResults
... other arguments are not used
```

# **Slots/List Components**

TestResults objects can be created by new ("TestResults", results) where results is a matrix. Objects of this class contain no slots (other than .Data), although the attributes dim and dimnames may be treated as slots.

## Methods

This class inherits directly from class matrix so any operation appropriate for matrices will work on objects of this class. show and summary methods are also implemented.

Functions in LIMMA which operate on TestResults objects include heatDiagram, vennCounts, vennDiagram, write.fit.

## Author(s)

Gordon Smyth

## See Also

02.Classes gives an overview of all the classes defined by this package. 08.Tests gives an overview of multiple testing.

# **Examples**

```
## Not run:
# Assume a data object MA and a design matrix
fit <- lmFit(MA, design)
fit <- eBayes(fit)
results <- decideTests(fit)
summary(results)
## End(Not run)</pre>
```

tmixture 145

tmixture

Estimate Scale Factor in Mixture of t-Distributions

#### **Description**

This function estimates the unscaled standard deviation of the log fold change for differentially expressed genes. It is called by the function <code>ebayes</code> and is not intended to be called by users.

#### Usage

```
tmixture.vector(tstat,stdev.unscaled,df,proportion,v0.lim=NULL)
tmixture.matrix(tstat,stdev.unscaled,df,proportion,v0.lim=NULL)
```

## **Arguments**

numeric vector or matrix of t-statistics

stdev.unscaled

numeric matrix conformal with tstatf containing the unscaled standard deviations for the coefficient estimators

df

numeric vector giving the degrees of freedom associated with tstat

proportion

assumed proportion of genes which are differentially expressed

v0.lim

numeric vector of length 2, assumed lower and upper limits for the estimated

unscaled standard deviation

## **Details**

The values in each column of tstat are assumed to follow a mixture of an ordinary t-distribution, with mixing proportion 1-proportion, and (v0+v1)/v1 times a t-distribution, with mixing proportion proportion. Here v1=stdev.unscaled^2 and v0 is the value to be estimated.

## Value

Numeric vector of length equal to the number of columns of tstat and stdev.unscaled.

# Author(s)

Gordon Smyth

## See Also

ebayes

146 toptable

toptable

Table of Top Genes from Linear Model Fit

#### **Description**

Extract a table of the top-ranked genes from a linear model fit.

## Usage

## Arguments

fit	list containing a linear model fit produced by lmFit, lm.series, gls.series or mrlm. For topTable, fit should be an object of class MArrayLM as produced by lmFit and eBayes.
coef	column number or column name specifying which coefficient or contrast of the linear model is of interest. For topTable, can also be a vector of column subscripts, in which case the gene ranking is by F-statistic for that set of contrasts.
number	maximum number of genes to list
genelist	data frame or character vector containing gene information. For topTable only, this defaults to fit\$genes.
A	matrix of A-values or vector of average A-values. For topTable only, this defaults to fit\$Amean.
eb	output list from ebayes (fit). If NULL, this will be automatically generated.
adjust.metho	od
	method used to adjust the p-values for multiple testing. Options, in increasing conservatism, include "none", "BH", "BY" and "holm". See p.adjust for the complete list of options. A NULL value will result in the default adjustment method, which is "BH".
sort.by	character string specifying statistic to rank genes by. Possibilities for topTable and toptable are "logFC", "AveExpr", "t", "P", "p", "B" or "none".  "M" is allowed as a synonym for "logFC" for backward compatibility. Other permitted synonyms are "A" or "Amean" for "AveExpr", "T" for "t" and "p" for "P". Possibilities for topTableF are "F" or "none". Possibilities for topTreat are as for topTable minus "B".
resort.by	character string specifying statistic to sort the selected genes by in the output data.frame. Possibilities are the same as for sort.by.
p.value	cutoff value for adjusted p-values. Only genes with lower p-values are listed.
lfc	cutoff value for log2-fold-change. Only genes with larger fold changes are listed.
	any other arguments are passed to ebayes if eb is NULL

toptable 147

#### **Details**

toptable is an earlier interface and is retained only for backward compatibility.

This function summarizes a linear model fit object produced by lmFit, lm.series, gls.series or mrlm by selecting the top-ranked genes for any given contrast. topTable and topTableF assume that the linear model fit has already been processed by eBayes. topTreat assumes that the fit has been processed by treat.

The p-values for the coefficient/contrast of interest are adjusted for multiple testing by a call to p.adjust. The "BH" method, which controls the expected false discovery rate (FDR) below the specified value, is the default adjustment method because it is the most likely to be appropriate for microarray studies. Note that the adjusted p-values from this method are bounds on the FDR rather than p-values in the usual sense. Because they relate to FDRs rather than rejection probabilities, they are sometimes called q-values. See help("p.adjust") for more information.

Note, if there is no good evidence for differential expression in the experiment, that it is quite possible for all the adjusted p-values to be large, even for all of them to be equal to one. It is quite possible for all the adjusted p-values to be equal to one if the smallest p-value is no smaller than 1/ngenes where ngenes is the number of genes with non-missing p-values.

The sort.by argument specifies the criterion used to select the top genes. The choices are: "logFC" to sort by the (absolute) coefficient representing the log-fold-change; "A" to sort by average expression level (over all arrays) in descending order; "T" or "t" for absolute t-statistic; "P" or "p" for p-values; or "B" for the lods or B-statistic.

Normally the genes appear in order of selection in the output table. If a different order is wanted, then the resort.by argument may be useful. For example, topTable(fit, sort.by="B", resort.by="logFC") selects the top genes according to log-odds of differential expression and then orders the selected genes by log-ratio in decreasing order. Or topTable(fit, sort.by="logFC", resort.by="logFC") would select the genes by absolute log-fold-change and then sort them from most positive to most negative.

topTableF ranks genes on the basis of moderated F-statistics for one or more coefficients. If topTable is called with coef has length greater than 1, then the specified columns will be extracted from fit and topTableF called on the result. topTable with coef=NULL is the same as topTableF, unless the fitted model fit has only one column.

Toptable output for all probes in original (unsorted) order can be obtained by topTable (fit, sort="none", n=In: However write.fit or write may be preferable if the intention is to write the results to a file. A related method is as.data.frame (fit) which coerces an MArrayLM object to a data.frame.

By default number probes are listed. Alternatively, by specifying p.value and number=Inf, all genes with adjusted p-values below a specified value can be listed.

The argument lfc gives the ability to filter genes by log-fold change. This argument is not available for topTreat because treat already handles fold-change thresholding in a more sophisticated way.

#### Value

Produces a dataframe with a row for the number top genes and the following columns:

genelist	one or more columns of probe annotation, if genelist was included as input
logFC	estimate of the log2-fold-change corresponding to the effect or contrast (for topTableF there may be several columns of log-fold-changes)
AveExpr	average log2-expression for the probe over all arrays and channels, same as ${\tt Amean}$ in the ${\tt MarrayLM}$ object
t.	moderated t-statistic (omitted for topTableF)

148 trigammaInverse

F moderated F-statistic (omitted for topTable unless more than one coef is

specified)

P. Value raw p-value

adj.P.Value adjusted p-value or q-value

B log-odds that the gene is differentially expressed (omitted for topTreat)

#### Author(s)

Gordon Smyth

#### See Also

An overview of linear model and testing functions is given in 06.LinearModels. See also p.adjust in the stats package.

# **Examples**

```
# See lmFit examples
```

trigammaInverse

Inverse Trigamma Function

## **Description**

The inverse of the trigamma function.

# Usage

```
trigammaInverse(x)
```

## **Arguments**

Х

numeric vector or array

## **Details**

The function uses Newton's method with a clever starting value to ensure monotonic convergence.

#### Value

Numeric vector or array y satisfying trigamma(y) == x.

#### Note

This function does not accept a data.frame as argument although the internal function trigamma does.

## Author(s)

Gordon Smyth

trimWhiteSpace 149

## See Also

```
trigamma
```

# **Examples**

```
y <- trigammaInverse(5)
trigamma(y)</pre>
```

trimWhiteSpace

Trim Leading and Trailing White Space

## Description

Trims leading and trailing white space from character strings.

## Usage

```
trimWhiteSpace(x)
```

# Arguments

х

character vector

## Value

A character vector of the same length as  $\mathbf{x}$  in which leading and trailing white space has been stripped off each value.

## Author(s)

Tim Beissbarth and Gordon Smyth

#### See Also

An overview of LIMMA functions for reading data is given in 03.ReadingData.

```
x <- c("a "," b ")
trimWhiteSpace(x)</pre>
```

150 unwrapdups

uniquegenelist

Eliminate Duplicate Names from the Gene List

## Description

Eliminate duplicate names from the gene list. The new list is shorter than the full list by a factor of ndups.

## Usage

```
uniquegenelist(genelist,ndups=2,spacing=1)
```

## **Arguments**

genelist vector of gene names

ndups number of duplicate spots. The number of rows of genelist must be divisible

by ndups.

spacing the spacing between duplicate names in genelist

#### Value

A vector of length length (genelist) / ndups containing each gene name once only.

#### Author(s)

Gordon Smyth

## See Also

```
unwrapdups
```

#### **Examples**

```
genelist <- c("A", "A", "B", "B", "C", "C", "D", "D")
uniquegenelist(genelist, ndups=2)
genelist <- c("A", "B", "A", "B", "C", "D", "C", "D")
uniquegenelist(genelist, ndups=2, spacing=2)</pre>
```

unwrapdups

Unwrap Duplicate Spot Values from Rows into Columns

## Description

Reshape a matrix so that a set of consecutive rows becomes a single row in the output.

## Usage

```
unwrapdups(M, ndups=2, spacing=1)
```

venn 151

#### **Arguments**

M a matrix.

ndups number of duplicate spots. The number of rows of M must be divisible by

ndups.

spacing the spacing between the rows of M corresponding to duplicate spots, spacing=1

for consecutive spots

#### **Details**

This function is used on matrices corresponding to a series of microarray experiments. Rows corresponding to duplicate spots are re-arranged to that all values corresponding to a single gene are on the same row. This facilitates fitting models or computing statistics for each gene.

#### Value

A matrix containing the same values as M but with fewer rows and more columns by a factor of ndups. Each set of ndups rows in M is strung out to a single row so that duplicate values originally in consecutive rows in the same column are in consecutive columns in the output.

## Author(s)

Gordon Smyth

## **Examples**

```
M <- matrix(1:12,6,2)
unwrapdups(M,ndups=2)
unwrapdups(M,ndups=3)
unwrapdups(M,ndups=2,spacing=3)</pre>
```

venn

Venn Diagrams

#### **Description**

Compute classification counts or plot classification counts in a Venn diagram.

## Usage

```
vennCounts(x, include="both")
vennDiagram(object, include="both", names, mar=rep(1,4), cex=1.5, lwd=1,
circle.col, counts.col, show.include, ...)
```

## **Arguments**

Х	numeric matrix of 0's and 1's indicating significance of a test. Usually created by decideTests.
object	either a TestResults matrix or a VennCounts object produced by vennCounts.
include	character string, of length one or two, specifying whether the diagram should give counts for genes up-regulated, down-regulated or both. See details. Choices are "both", "up" or "down".

venn venn

names	optional character vector giving names for the sets or contrasts
mar	numeric vector of length 4 specifying the width of the margins around the plot. This argument is passed to $par$ .
cex	numerical value giving the amount by which the contrast names should be scaled on the plot relative to the default.plotting text. See $par$ .
lwd	numerical value giving the amount by which the circles should be scaled on the plot. See ${\tt par.}$
circle.col	optional vector of color specifications defining the colors by which the circles should be drawn. See ${\tt par.}$
counts.col	optional vector of color specifications, of same length as $include$ , defining the colors by which the counts should be drawn. See $par$ .
show.include	logical value whether the value of include should be printed on the plot. Defaults to FALSE if include is a single value and TRUE otherwise
	any other arguments are passed to plot

#### **Details**

If a vennCounts object is given to vennDiagram, the include parameter is ignored. If a TestResults object is given, then it is possible to set include as a vector of 2 character strings and both will be shown.

#### Value

vennCounts produces a VennCounts object, which is a numeric matrix with last column "Counts" giving counts for each possible vector outcome. vennDiagram causes a plot to be produced on the current graphical device. For venDiagram, the number of columns of object should be three or fewer.

## Author(s)

Gordon Smyth, James Wettenhall and Francois Pepin

## See Also

An overview of linear model functions in limma is given by 06.LinearModels.

```
Y <- matrix(rnorm(100*6),100,6)
Y[1:10,3:4] <- Y[1:10,3:4]+3
Y[1:20,5:6] <- Y[1:20,5:6]+3
design <- cbind(1,c(0,0,1,1,0,0),c(0,0,0,1,1))
fit <- eBayes(lmFit(Y,design))
results <- decideTests(fit)
a <- vennCounts(results)
print(a)
vennDiagram(a)
vennDiagram(results,include=c("up","down"),counts.col=c("red","green"))</pre>
```

volcanoplot 153

## Description

Creates a volcano plot of log-fold changes versus log-odds of differential expression.

# Usage

```
volcanoplot(fit, coef=1, highlight=0, names=fit$genes$ID, xlab="Log Fold Change"
```

## Arguments

fit	an MArrayLM fitted linear model object
coef	integer giving the coefficient
highlight	number of top genes to be highlighted
names	character vector giving text labels for the probes to be used in highlighting
xlab	character string giving label for x-axis
ylab	character string giving label for y-axis
pch	vector or list of plotting characters. Default is integer code 16 which gives a solid circle.
cex	numeric vector of plot symbol expansions. Default is 0.35.
	any other arguments are passed to plot

## **Details**

A volcano plot is any plot which displays fold changes versus a measure of statistical significance of the change.

## Value

A plot is created on the current graphics device.

## Author(s)

Gordon Smyth

#### See Also

An overview of presentation plots following the fitting of a linear model in LIMMA is given in 06.LinearModels.

```
# See lmFit examples
```

154 weighted.median

weighted.median

Weighted Median

## **Description**

Compute a weighted median of a numeric vector.

## Usage

```
weighted.median(x, w, na.rm = FALSE)
```

## **Arguments**

X	a numeric vector containing th	ne values w	hose mean is to	be computed.

w a vector of weights the same length as x giving the weights to use for each

element of x.

na.rm a logical value indicating whether NA values in x should be stripped before the

computation proceeds.

## **Details**

If w is missing then all elements of x are given the same weight.

Missing values in w are not handled.

The weighted median is the median of the discrete distribution with values given by x and probabilities given by w/sum(w).

## Value

numeric value giving the weighted median

## See Also

```
median, weighted.mean
```

```
## GPA from Siegel 1994 wt <- c(5, 5, 4, 1)/15 x <- c(3.7, 3.3, 3.5, 2.8) xm <- weighted.median(x,wt)
```

write.fit 155

write.fit	Write MArrayLM Object to a File	
	, , ,	

## **Description**

Write a microarray linear model fit to a file.

## Usage

```
write.fit(fit, results=NULL, file, digits=3, adjust="none", method="separate", F
```

## Arguments

fit	object of class MArrayLM containing the results of a linear model fit
results	object of class TestResults
file	character string giving name of file
digits	integer indicating precision to be used
adjust	character string specifying multiple-testing adjustment method for the t-statistic P-values, e.g., "BH". See p.adjust for the available options. If NULL or "none" then the P-values are not adjusted.
method	character string, should the P-value adjustment be "global" or "separate" for each contrast.
F.adjust	character string specifying adjustment method for the F-statistic P-values.
sep	the field separator string. Values in the output file will be separated by this string.
	other arguments are passed to write.table

## **Details**

This function writes a tab-delimited text file containing for each gene (1) the average log-intensity, (2) the log-ratios, (3) moderated t-statistics, (4) t-statistic P-values, (5) F-statistic if available, (6) F-statistic P-values if available, (7) classification if available and (8) gene names and annotation.

## Value

No value is produced but a file is written to the current working directory.

## Author(s)

Gordon Smyth

## See Also

```
write in the base library.
```

An overview of linear model functions in limma is given by 06.LinearModels.

156 zscore

zscore

Z-score Equivalents

## **Description**

Compute z-score equivalents of gamma or t-distribution random deviates.

#### Usage

```
zscoreGamma(q, shape, rate = 1, scale = 1/rate)
zscoreT(x, df)
tZscore(x, df)
```

## **Arguments**

q, x	numeric matrix for vector giving deviates of a random variaable
shape	gamma shape parameter (>0)
rate	gamma rate parameter (>0)
scale	gamma scale parameter (>0)
df	<pre>degrees of freedom (&gt;0 for zscoreT or &gt;=1 for tZscore)</pre>

#### **Details**

These functions compute the standard normal deviates which have the same quantiles as the given values in the specified distribution. For example, if z <- zscoreT(x,df=df) then pnorm(z) equals pt(x,df=df).tZscore is the inverse of zscoreT.

Care is taken to do the computations accurately in both tails of the distributions.

## Value

Numeric vector giving equivalent deviates from a standard normal distribution (zscoreGamma and zscoreT) or deviates from a t-distribution (tzscore).

#### Author(s)

Gordon Smyth

#### See Also

gnorm, pgamma, pt in the standard distribution packages.

```
zscoreGamma(1, shape=1, scale=1)
zscoreT(2, df=3)
tZscore(2, df=3)
```

# Index

*Topic <b>IO</b>	*Topic <b>data</b>	
controlStatus, 31	as.MAList,16	
getLayout, 49	EList-class, 41	
getSpacing, 50	LargeDataObject-class, 62	
gridr, 53	MAList-class, 73	
printorder, 111	PrintLayout, 110	
read.ilmn, 122	RGList-class, 132	
read.ilmn.targets, 123	*Topic distribution	
readGAL, 119	fitFDist, 43	
readSpotTypes, 128	qqt, 116	
readTargets, 129	zscore, 156	
write.fit, 155	*Topic <b>documentation</b>	
*Topic <b>algebra</b>	02.Classes,2	
is.fullrank,59	03.ReadingData, $3$	
*Topic <b>array</b>	04.Background,4	
as.data.frame, 15	05.Normalization, $5$	
as.matrix, 16	06.Linear $Models, 5$	
avedups, 19	07.SingleChannel, 7	
avereps, 20	08.Tests, <b>8</b>	
blockDiag, 23	09.Diagnostics, $9$	
dim, 36	10.0ther, 10	
dimnames, 35	changeLog, 26	
exprs.MA,42	limmaUsersGuide, 63	
getEAWP,48	*Topic <b>file</b>	
normalizeMedianAbsValues,86	read.columns, 118	
uniquegenelist, 150	read.maimages, 125	
unwrapdups, 150	readHeader, 120	
*Topic character	readImaGeneHeader, 124	
alias2Symbol, $10$	removeExt, 131	
makeUnique, 72	*Topic <b>hplot</b>	
protectMetachar, 115	asMatrixWeights, 17	
removeExt, 131	heatdiagram, 53	
strsplit2,139	imageplot, 57	
trimWhiteSpace, 149	imageplot3by2,56	
*Topic <b>classes</b>	mdplot, 75	
as.MAList,16	modifyWeights, $80$	
EList-class,41	plotDensities,97	
LargeDataObject-class,62	plotFB, 99	
MAList-class, 73	plotlines, 100	
MArrayLM-class,74	plotMA, 102	
PrintLayout, 110	plotMA3by2,101	
RGList-class, 132	plotMDS, 104	
TestResults-class, 144	plotPrintTipLoess, 105	

plotRLDF, 106	normexp.fit,94
plotSA, 107	normexp.fit.control, 93
printHead, 109	normexp.signal, 96
volcanoplot, 153	printtipWeights, 112
*Topic <b>htest</b>	propexpr, 114
auROC, 18	residuals.MArrayLM,131
classifyTests, 28	selectModel, 137
contrasts.fit,29	*Topic <b>multivariate</b>
convest, 32	dupcor, 37
decideTests, 33	intraspotCorrelation, 58
ebayes, 39	normalizeBetweenArrays, 83
geneSetTest, 46	*Topic package
poolVar, 108	01.Introduction, 1
roast, 133	*Topic <b>programming</b>
romer, 135	isNumeric, 60
squeezeVar, 138	*Topic regression
targetsA2C, 142	arrayWeights, 13
TestResults-class, 144	arrayWeightsQuick, 12
tmixture, 145	designI2M, 27
toptable, 146	fitted.MArrayLM,44
venn, 151	genas, 44
*Topic <b>manip</b>	gls.series,51
cbind, 25	lm.series, 68
merge, 76	lmFit, 64
subsetting, 140	lmscFit, 66
*Topic <b>math</b>	makeContrasts, 71
trigammaInverse, 148	MArrayLM-class, 74
*Topic <b>methods</b>	modelMatrix, 79
helpMethods, 55	
summary, 141	mrlm, 81
*Topic <b>models</b>	printtipWeights, 112
-	QualityWeights, 117
anova.MAList-method, 11	removeBatchEffect, 130
arrayWeights, 13	residuals.MArrayLM, 131
backgroundCorrect, 21	selectModel, 137
bwss, 24	*Topic <b>smooth</b>
bwss.matrix, 23	ma3x3,70
fitted.MArrayLM,44	*Topic <b>univar</b>
genas, 44	weighted.median, 154
gls.series, 51	[.EList(subsetting), 140
kooperberg, 61	[.EListRaw(subsetting), 140
lm.series,68	[.MAList(subsetting), 140
lmFit,64	[.MArrayLM(subsetting), 140
lmscFit,66	[.RGList(subsetting), 140
loessFit, 69	01.Introduction, 1, 26
mergeScans, 77	02.Classes, 2, 15-17, 19, 20, 35, 36, 42,
mrlm, 81	49, 63, 74, 75, 110, 133, 141, 144
neqc, 82	03.ReadingData, 3, 13, 25, 32, 50, 51, 77
normalizeBetweenArrays,83	112, 115, 118–121, 125, 128–131,
normalizeForPrintorder,87	140, 149
normalizeQuantiles,89	04.Background, 4, 22, 62, 71, 83, 94, 96,
normalizeRobustSpline, 90	97
normalizeWithinArrays,91	05.Normalization, 5, 70, 81, 83, 85, 86,

00.00.02.100	
88–90, 92, 130	classifyTestsT(classifyTests), 28
06.LinearModels, 5, 8, 14, 15, 18, 30, 38,	cmdscale, 105
40, 53, 65, 69, 72, 80, 82, 100, 113,	coerce, <i>143</i>
137, 148, 152, 153, 155	coerce, RGList, exprSet2-method
07.SingleChannel, 6, 7, 8, 27, 59, 67,	(RGList-class), 132
143	combined, 2, 42, 73, 132
08.Tests, 7, 8, 11, 18, 29, 33, 34, 48, 134,	contrasts.fit, $6$ , $29$ , $45$ , $72$
136, 144	controlStatus, <i>4</i> , <i>31</i> , <i>103</i>
09. Diagnostics, 9, 12, 56, 58, 76, 98, 99,	convest, 9, 32
101, 103, 105, 106, 108	0.00.00.54.151
10.Other, 10, 23, 109	decideTests, 8, 28, 33, 54, 151
	designI2A (designI2M), 27
alias2Symbol, 9, 10	designI2M, 27
alias2SymbolTable,9	dim, 36, 36, 42, 73, 75, 132
alias2SymbolTable(alias2Symbol),	dimnames, 19, 20, 35, 35
10	dimnames <elist(dimnames), 35<="" td=""></elist(dimnames),>
anova, 9, 12, 23	dimnames <elistraw (dimnames),="" 35<="" td=""></elistraw>
anova.MAList, 24	dimnames <malist(dimnames), 35<="" td=""></malist(dimnames),>
anova.MAList	dimnames <rglist(dimnames), 35<="" td=""></rglist(dimnames),>
(anova.MAList-method), 11	dupcor, 37
anova.MAList-method, 11	duplicateCorrelation, $6,53$
array2channel(targetsA2C), 142	duplicateCorrelation(dupcor), 37
arrayWeights, 9, 13, 13	
arrayWeightsQuick, 12	eBayes, 6, 7, 45, 147
arrayWeightsSimple	eBayes (ebayes), 39
(arrayWeights), 13	ebayes, 6, 39, 43, 75, 145
as.data.frame, 2, 15, 15	EList, 2, 25
as.MAList, 3, 16	EList-class, 83
as.matrix, 16, 17	EList-class, 41
asMatrixWeights, 17	EListRaw, 2, 25, 127
AUC, 18	EListRaw-class, 82, 93, 122-124
auROC, 9, 18	EListRaw-class (EList-class), 41
avedups, 19	ExpressionSet, 42
avereps, 20	exprs, 17
	exprs.MA,42
backgroundCorrect, $4, 21, 91$	Extract, 140
barcodeplot, 8	
barcodeplot (geneSetTest), 46	fitFDist, 7, 40, 43
barcodeplot2, $8$	fitted, 44
barcodeplot2(geneSetTest),46	fitted.MArrayLM,44
bg.parameters,96	FStat, 8
blockDiag, 10, 23	FStat(classifyTests), 28
bwss, 9, 24, 24	
bwss.matrix, 9, 12, 23, 25	genas, 9, 44
	geneSetTest, $8, 46, 133$
cbind, 4, 25, 25, 76	getDupSpacing(getLayout),49
changeLog, $1, 26$	getEAWP,48
classifyTests, 28	getLayout, 3, 49, 110
classifyTestsF, 8, 75, 144	getLayout2 (getLayout), 49
classifyTestsF (classifyTests), 28	getSpacing, 3, 50
classifyTestsP, 8, 144	gls.series, $6,51$
classifyTestsP(classifyTests), 28	gridc, 3, 10
classifyTestsT, 8, 144	gridc(gridr), 53

gridr, 3, 10, 53	MArrayLM-class, 74
	marrayNorm, 16,74
heatDiagram, 8, 144	marrayRaw, 133
heatDiagram (heatdiagram), 53	Math, <i>61</i>
heatdiagram, $8,53$	mdplot, 9, 75
helpMethods, 55	mean, 58
	median, 154
image, 55, 58	merge, 4, 76, 77
imageplot, 9, 57	merge.RGList,73
imageplot3by2,9,56	merged, <i>132</i>
intraspotCorrelation, 7, 58, 67	mergeScans, 77
is.fullrank, 59	mergeScansRG (mergeScans), 77
is.numeric, 61	mixedModel2Fit, 38
isNumeric, 60	model.matrix, 6, 27, 80
	modelMatrix, $6, 79$
kooperberg, 4, 61	modifyWeights, 5, 18, 80
	mrlm, 6, 81
LargeDataObject, 2, 42, 73, 132	1111 ± 111, 0, 01
LargeDataObject-class, 62	ncol, 42, 73, 75, 132
length.EList (dim), 36	neqc, 4, 82, 94, 123
length.EListRaw (dim), 36	nlminb, 95
length.MAList (dim), 36	nonEstimable(is.fullrank),59
length.MArrayLM(dim), 36	normalize, 85
length.RGList (dim), 36	normalizeBetweenArrays, 5, 41, 73,
limma (01. Introduction), 1	83, 83, 89, 92, 132
LIMMA User's Guide, 5, 98	normalizeForPrintorder, 5, 73, 87,
limma-package (01.Introduction), 1	112, 132
limmaUsersGuide, 1, 63	normalizeMedianAbsValues, 5, 86
lm.fit, 67, 69	normalizeQuantiles, 5, 89
lm.series, 6, 68	normalizeQuantiles, 3, 89 normalizeRobustSpline, 5, 90, 92
lmFit, 6, 30, 45, 52, 60, 64, 68, 74, 81	normalizeWithinArrays, 4, 5, 22, 69,
lmscFit, 7, 64, 66	73, 91, 132
loess, 70	
loessFit, 5, 69, 92	normexp.fit, 4, 22, 83, 94, 94, 96, 97
lowess, 70	normexp.fit.control, 4, 83, 93, 123
10we33,70	normexp.signal, 4, 83, 94, 96, 96
MA.RG, 5, 73	nrow, 42, 73, 75, 132
MA.RG(normalizeWithinArrays), 91	
ma3x3,70	openPDF, 63
ma3x3.matrix,4	openVignette,63
ma3x3.spottedarray,4	p.adjust, 28, 34, 146-148, 155
maImage, 58	
<del>-</del>	par, 57
make.names, 72	pgamma, 156
make.unique, 73	plotDensities, 9, 10, 97
makeContrasts, 6, 71	plotFB, 9, 99
makeUnique, 4, 72	plotlines, 7, 100
MAList, 2, 16, 25, 37, 58, 66, 84, 85, 92	plotMA, 9, 73, 75, 102
MAList-class, 12, 76	plotMA3by2, 9, 101
MAList-class, 73	plotMDS, 9, 104
maNorm, 93	plotPrintorder,9
maNormScale, 85	plotPrintorder
MArrayLM, 2, 6, 30, 65, 67	(normalizeForPrintorder)
MArrayLM-class 40	X /

plotPrintTipLoess, 9, 73, 105	RG.MA(normalizeWithinArrays),91
plotRLDF, 106	RGList, 2, 21, 25, 84, 127
plotSA, 10, 107	RGList-class, 76, 78
points, 99, 103, 108	RGList-class, 132
poolVar, 10, 108	rlm, 82
printHead, 2, 109	roast, 8, 47, 48, 133, 136
PrintLayout, 73, 110, 127, 132	romer, 8, 48, 133, 134, 135, 142
PrintLayout-class, 21, 70, 87, 91, 113	romer2 (romer), 135
PrintLayout-class (PrintLayout),	
110	selectModel, 7, 8, 137
printorder, 3, 88, 111	show, 42, 73, 132, 144
printtipWeights, 112	show,LargeDataObject-method
printtipWeightsSimple	(LargeDataObject-class),62
(printtipWeights), 112	show, TestResults-method
propexpr, 114, 123	(TestResults-class), 144
protectMetachar, 115	showMethods, 55
pt, 156	spotc, 3
- /	spotc(gridr),53
gnorm, 156	spotr, 3
qqnorm, 116	spotr(gridr), 53
qqt, 116	squeezeVar, 40, 137, 138
QualityWeights, 3, 117, 126	strsplit, <i>139</i> , <i>140</i>
~	strsplit2, <i>3</i> , 139
rbind,4	subsetted, 2, 42, 73, 132
rbind.EList (cbind), 25	subsetting, 60, 140
rbind.EListRaw (cbind), 25	summary, 141, 141
rbind.MAList (cbind), 25	summary.TestResults
rbind.RGList (cbind), 25	(TestResults-class), 144
read.columns, 3, 118, 122, 128	Sweave, $63$
read.Galfile, 120	symbols2indices, 135, 142
read.ilmn, 3, 83, 115, 122, 123, 124	Sys.getenv, $63$
read.ilmn.targets, 123, 123	Sys.putenv, 63
read.imagene, 3, 125	Sys. Pacent, or
read.imagene (read.maimages), 125	targetsA2C, 8, 142
read.maimages, 3, 41, 119, 121, 125, 132	TestResults, 2, 8, 34
read.table, 119, 128, 129	TestResults-class, 144
readGAL, 3, 31, 119	text, 104, 106
readGenericHeader, 3	tmixture, 145
readGenericHeader (readHeader),	tmixture.matrix, 7, 40
120	tmixture.vector,7
readGPRHeader, 3	topTable, 7
readGPRHeader (readHeader), 120	topTable(toptable), 146
readHeader, 120	toptable, 7, 146
readImaGeneHeader, 3, 124	topTableF, 7
readSMDHeader (readHeader), 120	topTableF (toptable), 146
readSpotTypes, 4, 128	topTreat, 40
readTargets, 3, 123, 124, 129, 143	topTreat(toptable), 146
remlscore, 59	treat, 6, 147
removeBatchEffect, 130	treat (ebayes), 39
removeExt, 3, 131	trigamma, 149
residuals, 132	trigammaInverse, 43, 148
residuals, 132 residuals.MArrayLM, 131	trimWhiteSpace, 149
RG.MA, 10	tZscore (zscore), 156
110 • 1111, 10	$\cup \cup $

```
uniquegenelist, 3, 150
uniqueTargets (modelMatrix), 79
unwrapdups, 6, 150, 150
venn, 151
vennCounts, 8, 144
vennCounts (venn), 151
vennDiagram, 8, 144
vennDiagram (venn), 151
{\tt vignette}, 63
volcanoplot, 7, 153
vsn, 85
vsnMatrix, 85
weighted.mean, 154
weighted.median, 154
wilcox.test,48
wilcoxGST, 8, 134, 136
wilcoxGST (geneSetTest), 46
write, 147, 155
write.fit, 7, 8, 144, 147, 155
wtarea (QualityWeights), 117
wtflags(QualityWeights), 117
wtIgnore.Filter(QualityWeights),
       117
zscore, 156
zscoreGamma, 7
zscoreGamma (zscore), 156
zscoreT,7
zscoreT (zscore), 156
```